**English**

Instructions for Use: MDL/MILO Dental Prosthetics

2

**Español**

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L02070 Rev B FEB 2022

## English

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic Instructions for Use
	Manufacturer
	Intra-Lock products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Devices Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is printed on the product label.
	Reference/article number
	Lot/batch number
	Do not re-use
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
	Caution: U.S. Federal law restricts these devices for sale, distribution, and use by, or on the order of, a dentist or physician.
	EU Authorised Representative
	Do not use if package is damaged
	Medical Device
	Non-sterile
	Home Position

This document supersedes all prior revisions. Original language is English.

## **DESCRIPTION**

Intra-Lock MDL and MILO prosthetic components are intended for the restoration of MDL and MILO dental implants within the specific indications of each implant system. Intra-Lock MDL and MILO abutments are manufactured from Titanium Alloy (Ti-6Al-4V ELI). Temporary abutments are manufactured from PEEK. The label on each prosthetic package contains important product information. All prosthetic components are supplied non-sterile.

## **INDICATIONS FOR USE**

Intra-Lock MDL and MILO prosthetics are intended for use with Intra-Lock MDL and MILO dental implants to restore partially or fully edentulous patients. The abutments may be used in multiple tooth application in either the mandible or maxilla and to support removable or fixed prosthesis.

MILO abutments may also be used for single tooth rehabilitation of single maxillary lateral incisors and mandibular lateral and central incisors.

PEEK prosthetics are intended for short term temporary restorations (<30 days).

## **CONTRAINDICATIONS**

Patients with alcohol addiction or psychiatric disorders, blood dyscrasias, uncontrolled diabetes, hyperthyroidism, oral infections, malignancies, or patients who have had myocardial infarction within the last 12 months.

Patients with systemic diseases that compromise the immune system, such as AIDS, patients on medications that would compromise healing of an implant site, patients with a history of poor or noncompliance to oral hygiene procedures, or patients who cannot maintain oral hygiene procedures if implants are placed.

Tobacco usage increases the occurrence of complications and failures.

## **DIRECTIONS FOR USE**

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. Intra-Lock strongly recommends completion of dental implant courses and strict adherence to the Instructions for Use (IFU) that accompany Intra-Lock products.

## **WARNINGS AND PRECAUTIONS**

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any Intra-Lock IFU. Clinicians are responsible for understanding the appropriate technical use of Intra-Lock prosthetic components. Additional technical information is available upon request from Intra-Lock or may be viewed and/or downloaded at [www.intra-lock.com](http://www.intra-lock.com). Contact Intra-Lock Customer Care or your local representative with any questions you have regarding a specific IFU.

Dental abutments can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. The use of angled abutments in the posterior region of the mouth is not recommended due to possible failure of the implant.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. Intra-Lock assumes no responsibility for attempted re-use or re-sterilization between patients.

## COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with MDL and MILO prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) implant and/or abutment breakage; (3) infection requiring revision of the dental implant; (4) nerve damage that could cause permanent weakness, numbness, or pain; (5) histologic responses possibly involving macrophages and/or fibroblasts; (6) formation of fat emboli; (7) implant loosening requiring revision surgery; (8) maxillary sinus perforation; (9) labial or lingual plate perforation; and (10) bone loss possibly resulting in revision or removal.

## HANDLING AND STERILIZATION

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. Non-sterile product must be cleaned and sterilized. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations.
- 2) Brush the product to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse product under running utility tap water.
- 4) Place product in the sterile container filled with the prepared detergent solution and sonicate for two (2) minutes minimum.
- 5) Thoroughly rinse product under running utility tap water.
- 6) Spray or wipe product with 70% IPA.
- 7) Blot product dry with clean lint free cloth.

For sterilization of non-sterile product, place product in an FDA cleared sterilization bag or wrap and run through one of the following qualified sterilization cycles:


Sterilization Cycles				
<b>Reference:</b>	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	UK HTM 01-01 Part C:2016
<b>Type:</b>	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	30 minutes	20 minutes	20 minutes

It is recommended to include a 30-minute cool-down period before removing the product from the sterilization bag or wrap.

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