<b>English</b>	Instructions for Use: MDL/MILO Dental Implants	2
<b>Español</b>	Instrucciones de uso: MDL / MILO implantes dentales	5
<b>Français</b>	Instructions d'utilisation : MDL / MILO implants dentaires	8

BioHorizons  
2300 Riverchase Center  
Birmingham, AL 35244 USA  
TOLL FREE 888.246.8338  
TEL 205.967.7880 FAX 205.870.0304  
[www.biohorizons.com](http://www.biohorizons.com)

L01081 Rev B JULY 2022

## English

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic Instructions for Use
	Manufacturer
	Intra-Lock products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Devices Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is printed on the product label.
	Reference/article number
	Lot/batch number
	Do not re-use
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
	Caution: U.S. Federal law restricts these devices for sale, distribution, and use by, or on the order of, a dentist or physician.
	EU Authorized Representative
	Do not use if package is damaged
	Medical Device
	Non-sterile
	Home Position

This document supersedes all prior revisions. Original language is English.

This document applies to Intra-Lock Mini Drive-Lock (MDL) and MILO Dental Implants.

#### **DESCRIPTION**

Intra-Lock MDL and MILO implants are manufactured from biocompatible titanium alloy and are single use only.

#### **INDICATIONS MINI DRIVE-LOCK™ DENTAL IMPLANT SYSTEM**

Mini Drive-Lock™ Dental Implants are intended for use as a self-tapping titanium screw for transitional or intra-bony long-term applications.

Mini Drive-Lock™ Dental Implants are indicated for long-term maxillary and mandibular tissue-supported denture stabilization. Multiple implants should be used and may be restored after a period of time or placed in immediate function.

#### **INDICATIONS MILO™ DENTAL IMPLANT SYSTEM**

MILO™ Dental Implants are indicated for long-term maxillary and mandibular tissue-supported denture stabilization. Multiple implants may be restored after a period of time or placed in immediate function. They are also indicated for the rehabilitation of single maxillary lateral incisors and mandibular lateral and central incisors. The implants may be restored after a period of time or placed in immediate function.

#### **CONTRAINDICATIONS**

Patients with alcohol addiction or psychiatric disorders, blood dyscrasias, uncontrolled diabetes, hyperthyroidism, oral infections, malignancies, or patients who have had myocardial infarction within the last 12 months.

Patients with systemic diseases that compromise the immune system, such as AIDS, patients on medications that would compromise healing of an implant site, patients with a history of poor or non-compliance to oral hygiene procedures, or patients who cannot maintain oral hygiene procedures if implants are placed.

Tobacco usage increases the occurrence of complications and failures.

#### **DIRECTIONS FOR USE**

Proper surgical procedures and techniques are the responsibility of the medical professional. Each dentist must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. Intra-Lock strongly recommends completion of postgraduate dental implant education and strict adherence to the instructions and procedures for use that accompany Intra-Lock implant products. NOTE: A pre-operative 30-second rinse with a 0.12% Chlorhexidine Digluconate solution is recommended. (The influence of 0.12% Chlorhexidine Digluconate Rinses on the Incidence of Infectious Complications and Implant Success. Lambert, Paul M., et al, J Oral Maxillofacial Surgery 55:25-30, 1997, Suppl5).

IMPORTANT! (1) The implant procedure manual should be reviewed prior to attempting placement. (2) The clinician should verify that the proper instruments are available and at hand prior to attempting placement. (3) MDL implants are not recommended to be placed in sites with a cortical plate thickness greater than 3mm.

#### **WARNINGS AND PRECAUTIONS**

Appropriate training in proper implant surgery technique is strongly recommended prior to implant use. Improper technique can result in implant failure and/or loss of supporting bone. Appropriate x-ray films and/or CT scans should be utilized to determine (1) if adequate bone width and depth are available at the desired implant site and (2) the location of important anatomical landmarks, such as the mandibular canal, maxillary sinuses, and adjacent teeth.

Clinician judgment should be used in determining the minimum post-implantation time before placing the implants in occlusal function. An adequate number of implants should be used to provide support and to distribute the load to the abutments. Intra-Lock dental implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of Intra-Lock dental implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. A certain percentage of implants may fail to achieve or to maintain osseointegration, as demonstrated by mobility, and should be removed. Potential causes of abutment fracture include, but not limited to: casting beyond a 15-degree angle, casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics.

Devices should not be used in patients with a known allergy or sensitivity to the device material(s).

### **COMPLICATIONS AND ADVERSE EFFECTS**

The risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) infection requiring revision of the dental implant; (4) nerve damage that could cause permanent weakness, numbness, or pain; (5) histologic responses possibly involving macrophages and/or fibroblasts; (6) formation of fat emboli; (7) loosening of the implant requiring revision surgery; (8) perforation of the maxillary sinus; (9) perforation of the labial and lingual plates; and (10) bone loss possibly resulting in revision or removal.

### **HANDLING AND STERILIZATION**


This product is supplied sterile and should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove from the package only after the correct size has been determined and the surgical site has been prepared. Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. Implants are single use only, and re-use should not be attempted. Following this guidance eliminates the risk of cross-patient contamination from secondary use of this device. Intra-Lock assumes no responsibility for attempted re-use or re-sterilization.

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