INTRASPIN®

Leukocyte- and Platelet-Rich Fibrin

L02065 Rev E JUL 2022

English	Instructions for use: IntraSpin
Español	Instrucciones de uso: IntraSpin
русский	Инструкция по применению: IntraSpin
DEUTSCH	Gebrauchsanweisung: IntraSpin
FRANÇAIS	Instructions d'utilisation: IntraSpin
ITALIANO	Istruzioni per l'uso: IntraSpin
PORTUGUÊS	Instruções de utilização: IntraSpin
Türk	Kullanım talimatları: IntraSpin
简体中文	使用说明 : IntraSpin
日本語	使用説明書 : IntraSpin
한국어	사용 설명서: IntraSpin
عربي	إرشادات الاستخدام: IntraSpin
Polskie	Instrukcja użycia: IntraSpin
čeština	Návod k použití: IntraSpin

BioHorizons

2300 Riverchase Center
Birmingham, AL 35244 USA
+1-205-967-7880

EC REP Altatec GmbH

Maybachstraße 5
71299 Wimsheim, Germany
Tel +49 7044 94450

IntraSpin®, Xpression® & L-PRF® are registered trademarks of BioHorizons; Vacuette® is a registered trademark of Greiner Bio-One International AG.; Enzymax® is a registered trademark of Hu-Friedy Mfg. Co., LLC.

The symbol table below is for reference only. Refer to product or product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic instructions for use.
	Manufacturer.
	BioHorizons products carry the CE mark and fulfill the requirements of the Medical Devices Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745.
	Reference/ article number.
	Lot/ batch number.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.
	Date of manufacture.
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician.

Symbol	Symbol Description
	Home position.
	Do not use if package is damaged. Contact BioHorizons or your local representative for disposition.
	Medical Device.
	Non-sterile.
	Keep dry.
	Fragile; handle with care.
	Temperature limits.
	This way up.
	Humidity limits.
	Warning; Biological hazard.
	Important note(s).
	Separate collection of electric and electronic devices.
	Authorized representative in the European Community.

1 INDICATIONS FOR USE

The IntraSpin System is intended to be used for the safe and rapid preparation of autologous Leukocyte- and Platelet-Rich Fibrin (L-PRF) from a small sample of blood at the patient's point of care. The L-PRF is mixed with autograft and/or allograft bone prior to application to a bony defect for improving handling characteristics. Observing all information in the Instructions for Use is also a part of the intended use.

2 CONTRAINDICATIONS

The IntraSpin centrifuge is only meant for the purpose stated in the intended use of the device. Any other use of the device is considered non-intended. Use of the IntraSpin centrifuge is contraindicated in the presence of one or more of the following clinical situations:

- Patients with alcohol addiction or psychiatric disorders, blood dyscrasias, uncontrolled diabetes, hyperthyroidism, oral infections, malignancies or patients who have had myocardial infarction within the last 12 months.
- Patients with systemic diseases that compromise the immune system, such as AIDS, patients on medications that would compromise healing of an implant site, patients with a history of poor or noncompliance to oral hygiene procedures.
- Patients who are participating in anti-coagulant therapy.

3 WARNINGS AND SAFETY NOTES

- No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.
- This product is not authorized for sale in every market. Please consult with your local representative for additional information.
- The operating instructions are a part of the device. They must always be kept readily available. Instructions for use are available free of charge at <http://ifu.biohorizons.com> or in printed form upon request from BioHorizons or your local distributor. Additional technical information is available upon request from BioHorizons or may be viewed and/or downloaded at www.biohorizons.com. Contact BioHorizons Customer Care or your local representative with any questions you have regarding the IFU.
- Before operating the centrifuge system, user must read and understand the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is used. Meeting the country-specific requirements concerning occupational safety with regard to the use of centrifuges at the workplaces is the responsibility of the user.
- This centrifuge is a state-of-the-art piece of equipment which is safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- Maintain the centrifuge in a location where ambient temperature and humidity are within the ranges provided in these instructions for use in Section 12. If the centrifuge is repeatedly used, the centrifuging chamber may heat up. Allow time for the chamber to cool down.
- To avoid damage due to condensation, when switching from a cold to a warm room, the centrifuge must be allowed to warm up for at least 3 hours in the warm room before being connected to the

electrical power. When switching from warm to cold, the centrifuge must be allowed to run for approximately 30 minutes in the cold room.

- Before using the centrifuge, check the rotor for firm placement.
- A rotor or centrifuge accessory that is showing traces of corrosion or mechanical damage should not be used and should be replaced as soon as possible. The rotor should not be used past its expiration date.
- The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.
- The centrifuge should be installed on a good, stable base.
- The centrifuge must not be moved or knocked during operation.
- When the centrifuge is running, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- In case of fault or emergency release, never touch the rotor before it stops rotating.
- When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability. If balance is not achieved, an error message will be displayed by the centrifuge to warn users.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with inflammable or explosive materials or materials that react with one another producing energy.
- No biosafety systems are available for this centrifuge.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of the rotor or accessories.
- Repairs must only be carried out by personnel authorized by the manufacturer.
- To ensure highest level of clinical safety, the IntraSpin System devices with direct patient contact are manufactured with biocompatible materials.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

4 INTRASPIN SYSTEM COMPONENTS

COMPONENT	QUANTITY PER SYSTEM	LEGAL MANUFACTURER
IntraSpin Centrifuge including:	1	BioHorizons 2300 Riverchase Center Birmingham, AL 35244 USA +1-205-967-7880
Power Cable	1	
Fuse	2	
Hex Hand Wrench	1	
Greiner Bio-One Tube 9ml Serum Clot activator, red cap (single use)	100	Greiner Bio-One North America Inc. 4238 Capital Drive, Monroe 28110 North Carolina USA +1-704-261-7800
Greiner Bio-One balancing tubes White Cap 9ml No additive	50	
Greiner Safety Blood Collection Set + Holder, 21G (single use)	24	
Latex Free Tourniquet	1	Propper Manufacturing 30-04 Skillman Ave., Long Island NY 11101 USA +1-718-392-6650
Test Tube Rack	1	Heathrow Scientific LLC 620 Lakeview Pkwy, Vernon Hills, IL 60061 USA +1-847-816-5070
Tissue Regeneration Kit including:	1	BioHorizons 2300 Riverchase Center Birmingham, AL 35244 USA +1-205-967-7880
Surgical Curved Scissors	1	
Surgical Tissue Forceps	1	
Round Stainless-Steel Bowl	1	
Rectangular Stainless-Steel Bowl	1	
Dual Biomaterial Carrier Spatula	1	
Dual Biomaterial Packer	1	
Xpression® Box	1	

Only verified compatible components for direct use with the IntraSpin centrifuge are recommended and warranted:

COMPATIBLE PART #	DESCRIPTION
455092	Tube 9ml Serum Clot activator, red cap (50 pcs)
455001	White Cap 9ml No additive blood collection tube (50 pcs)
BHEXZ (E613)	IntraSpin Hex Key, 110v & 220v
BROTORZ (E3694)	IntraSpin Rotor, 100v & 220v
BPOWER110Z (E1673)	IntraSpin Power Cord, 110v
BPOWER220Z (E1669)	IntraSpin Power Cord, 220v
BTUBEHOLDZ (E872 x 1)	IntraSpin Tube Holder Replacement
BFUSE110Z (E997)	IntraSpin Fuse 110v
BFUSE220Z (E891)	IntraSpin Fuse 220v

Refer to the following table for material(s) of devices with direct patient contact:

Devices with Direct Patient Contact	Material (main elements)
Surgical Tissue Forceps	Stainless Steel (Iron, Chromium)
Dual Biomaterial Carrier Spatula	Stainless Steel (Iron, Chromium)
Dual Biomaterial Packer	Stainless Steel (Iron, Chromium)

5 BRIEF CENTRIFUGE SETUP

Remove and save transport bolts from bottom of centrifuge.

Attach AC cable and plug into electrical outlet.

Power centrifuge on by using the rocker switch on the back of the device.

Select speed and time: Speed = 2700 & Time = 12:00 min.

Press START.

The centrifuge cover will open automatically at the end of each cycle.

After the first procedure, the timing and speed are recorded in the centrifuge memory unless the settings are changed.

6 BLOOD COLLECTION TUBES AND BLOOD COLLECTION SET CAUTIONS AND INSTRUCTIONS

- Do not use tubes if foreign matter is present.
- Blood collection tubes must be allowed to fill completely.
- Handle all biological samples and blood collection “sharps” (e.g., needles, and blood collection sets) according to the policies and procedures of your facility.
- Do not bend the needle.
- Do not forcefully release or re-activate the needle safety mechanism after it has been activated.
- Obtain appropriate medical attention in the case of any exposure to biological samples (e.g., via puncture injury) due to the possible transmission of HIV (AIDS), viral hepatitis or other infectious diseases.
- Discard all blood collection “sharps” in approved biohazard containers.
- Transferring a sample from a syringe to a tube is not a recommended procedure.
- If blood is collected through an intravenous (IV) line, follow the policies and procedures of your institution to ensure that the line has been cleared of IV solution before beginning to fill the blood collection tubes.
- Blood clotting accelerant may appear white on the tube surface, which has no effect on the performance of the tubes. If any other discoloration or precipitates are present in the tube, it should not be used.
- Do not use the tubes after the expiration date.
- Store blood collection tubes at 4–25°C (40–77°F).
- Store blood collection set (needle and holder) at 4–36°C (40–97°F).
- Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the tube quality (i.e., vacuum loss, coloring, etc.).
- To prevent backflow, place the patient’s arm in a downward position, hold the tube with the cap up, release the tourniquet as soon as blood starts to flow into the tube, avoid tube contents coming in contact with cap or end of the needle during venipuncture.
- Be sure that the following materials are readily accessible before performing venipuncture: all necessary blood collection tubes, identified labels for positive patient identification of samples, blood collection needles and holders, alcohol swab for cleansing the puncture site, clean gauze, tourniquet, adhesive plaster or bandage, approved biohazard container. For protection against exposure to bloodborne pathogens, appropriate PPE (Personal Protective Equipment) is recommended (e.g., gloves, laboratory coat, goggles, etc.).

Venipuncture Technique and Blood Sample Collection:

The blood collection must be made as quickly as possible, since there is no anticoagulant in the collection tubes. The blood sample will begin to coagulate immediately. Wear gloves during venipuncture and when handling blood collection tubes to minimize exposure to hazards. Prior to the blood draw, wipe the top of the blood tube cap(s) with a disinfectant wipe of your choice. Remove the cover over the valve section of the needle. Prepare venipuncture site with an appropriate antiseptic. Do not palpate venipuncture area after cleansing. Place the patient's arm in a downward position. Remove the needle cap. Perform the venipuncture with the arm downward and tube cap up. Immobilize the needle with tape when necessary. Push the blood collection tube into the holder and onto the needle valve puncturing the rubber diaphragm of the blood collection tube. Center the blood collection tubes in the holder when penetrating the cap to prevent sidewall penetration and subsequent premature vacuum loss. Remove the tourniquet as soon as blood appears in the blood collection tube. During procedure, always hold the collection tube in place by pressing it with a thumb. This will ensure a complete vacuum draw. The blood collection tube will fill automatically. If no blood flows into collection tube or if blood flow ceases before an adequate specimen is collected, the following steps are suggested to complete a satisfactory collection:

- Push the blood collection tube forward to ensure the cap has been penetrated.
- Confirm the correct position of the needle in the vein.
- If blood still does not flow, remove and appropriately discard the collection tube. Obtain a new collection tube and push into holder.
- If the second collection tube does not draw, remove and appropriately discard the needle and the collection tube. Repeat the procedure.
- When the maximum volume fill line of the blood collection tube has been reached, gently remove it from the holder. Repeat with a second blood collection tube.

Gently invert each collection tube immediately upon removing from the holder. Do not shake the tubes filled with blood sample. Vigorous mixing may cause foaming or hemolysis. Insufficient mixing or delayed mixing in serum tubes may result in delayed clotting. Upon completion of blood sample collection, remove the needle from the vein. Activate the safety mechanism (safety shield) of the needle by pressing in both sides of the hub to engage the lock. Slide the safety mechanism backward until an audible click is heard. Do not recap the needle as this increases the risk of needle stick injury and blood exposure. Dispose of the used needle with holder using a suitable biohazard disposal container. Apply pressure to the puncture site with a dry sterile swab until the bleeding stops. If desired, apply a bandage once clotting has occurred. It is recommended that filled collection tubes, be kept in an upright position. Once the second blood collection tube is full, remove it and place the first and second tubes into the centrifuge on opposite locations to counterbalance the rotor. Close the cover of the IntraSpin centrifuge and press the START button to allow it to spin for 12 minutes.

If more than two tubes of blood are required, please follow this alternative procedure: After the first two tubes of blood are collected and gently inverted, immediately place them into the IntraSpin centrifuge, opposite from each other to ensure the centrifuge is properly balanced. Close the cover and press the START button and allow the centrifuge to run while you collect the remaining tubes of blood. Press the STOP button and allow the centrifuge to come to a full stop. The cover will pop open; immediately place the remaining tubes in

the centrifuge opposite from each other to ensure proper balance and press the START button to reset and complete recommended protocol.

Always place the tubes in pairs and place them in opposite positions to balance the centrifuge rotor. The tubes must always be balanced in the rotor before pressing the START button or this may cause serious damage to the centrifuge, improper coagulation, and/or separation. If the tubes are not properly balanced, there will be too much vibration during centrifugation and a poor L-PRF fibrin clot will result.


If you have an odd number of blood samples to centrifuge, then place a white cap balancing tube (e.g., 455001), filled with water to the indicated full line, opposite to the un-paired tube in the rotor. This will allow for proper balancing of the centrifuge.

Begin centrifugation immediately after collecting the blood samples. Delays affect the blood separation procedure and result in a poor L-PRF fibrin clot.

7 L-PRF PREPARATION

After centrifugation, three segments are visible:


1. Upper Segment = platelet poor plasma (PPP).
2. Middle Segment = fibrin clot: L-PRF.
3. Lower Segment = red blood cell clot.

L-PRF fibrin membranes or plugs must be prepared relatively quickly: 0-15 minutes after centrifugation or the clot will shrink in volume by releasing the trapped serum. After centrifugation, remove the rubber stopper from each tube. Using the Surgical Tissue Forceps remove the L-PRF clot from the tube. Gently scrape the red blood cell clot from the L-PRF® fibrin clot just below the union, using the Dual Biomaterial Carrier Spatula, so that only a minimal, residual amount of red blood cells are attached to L-PRF clot. Place the fibrin clot onto the Xpression Perforated Tray.

8 FIBRIN MATRIX PREPARATION

The Xpression Box enables the fabrication of fibrin membranes of constant thickness with ease. The exudate can be collected from the Xpression Collection Tray, underneath the Xpression Perforated Tray. The Xpression Box includes L-PRF plug fabrication cylinders and a piston to fabricate L-PRF plugs that easily fit post-extraction sockets.

Representative Xpression Box and components

Protocol #1 L-PRF Membrane

Place each of the fibrin clots on the Xpression Perforated Tray. Once all of the fibrin clots are placed, place the Xpression Compression Plate and Xpression Weighted Cover over the fibrin clots without exerting any pressure over the clots.

Allow the weight of the cover to slowly press down the fibrin clot while the exudate is filtered to the bottom of the tray. Do not apply pressure to the weighted cover. Gravitational force on the weighted cover will gently compress the clot and express the serum from the L-PRF clot without damaging the fibrin network.

Wait at least 5 minutes before removing and using any fibrin membranes. Do not remove any fibrin membranes until actual time of use. The fibrin membranes should be used as quickly as possible but may remain in the Xpression Box for a period of 2.5 to 3 hours as long as they are re-hydrated with exudate (MLD601, R43069r).

Protocol #2 L-PRF Plug

Place a fibrin clot inside the white plug fabrication cylinder. Use the piston to slowly press the clot inside the white L-PRF plug fabrication cylinder. Continue to press until the top edge of the piston is flush with the top edge of the white L-PRF plug fabrication cylinder. With this technique, one will be able to form a thick, round fibrin plug for the extraction socket. For a single tooth, one L-PRF plug may be sufficient. Pre-molars may need two L-PRF plugs, and three L-PRF plugs may be needed for molars, depending on the size of the extraction socket and the size of the fibrin clot created.

The working properties of L-PRF provide an excellent medium for use in combination with your biomaterial of preference. Utilizing any of the following mixing protocols, the biomaterial is captured in the fibrin matrix increasing its handling and biologic capacity.

Protocol #3- Biomaterial/L-PRF Mixture


To create a 'putty like' mixture that can be gently formed with the biomaterial instrument into the desired shape and thickness use the following protocol: Gently cut the L-PRF fibrin membrane into small pieces in a sterile dish with the Surgical Curved Scissors. Add the desired amount of bone graft material. Thoroughly mix the L-PRF and bone graft material. This mixture can be placed into defects using the Dual Biomaterial Carrier Spatula.

Protocol #4- Biomaterial/L-PRF Matrix Mixture

Place the predetermined amount of bone graft material into a sterile bowl or tray. Dip the expressed L-PRF membrane(s) or pieces of the L-PRF membrane into the graft material covering the entire surface

area of the L-PRF membrane with graft material. Alternatively, the graft material may be sprinkled onto the L-PRF membrane covering the entire surface area with graft material. Note: A wetter L-PRF membrane may retain slightly more graft material than a dryer L-PRF membrane. The graft material should cling to the surface of the L-PRF, however, if desired, gently press the graft material onto the L-PRF membrane. The Surgical Tissue Forceps can be used to place this mixture into the defect.

Protocol #5- Biomaterial Hydration

Add the desired amount of bone graft material into a sterile bowl or tray. Utilize the exudate from the bottom of the Xpression Collection Tray to hydrate the graft material. Thoroughly mix the exudate and bone graft material. This mixture can be placed into defects using the Dual Biomaterial Carrier Spatula.

9 TISSUE REGENERATION KIT CLEANING AND STERILIZATION

The Tissue Regeneration Kit (including Xpression® Box, Surgical Curved Scissors, Surgical Tissue Forceps, Round Stainless-Steel Bowl, Rectangular Stainless-Steel Bowl, Dual Biomaterial Carrier Spatula, and Dual Biomaterial Packer) is NOT supplied sterile. Remove and discard any shipping material before initial cleaning and sterilization. Clean and sterilize the devices before each use.

Disassemble the Xpression Box before each cleaning cycle. Remove the Xpression Compression Plate and the Xpression Perforated Tray from the Xpression Collection Tray. Remove the Piston from the Xpression Perforated Tray. The L-PRF plug fabrication cylinders and piston grommet are not intended to be removed from the Xpression Perforated Tray for cleaning and sterilization.

Cleaning Steps:

- 1 Remove any visible debris from the Xpression Box, Surgical Curved Scissors, Surgical Tissue Forceps, Round Stainless-Steel Bowl, Rectangular Stainless-Steel Bowl, Dual Biomaterial Carrier Spatula, and Dual Biomaterial Packer using a soft-bristled brush dampened with a broad-spectrum cleaning detergent such as Hu-Friedy's Enzymax® or equivalent. Pay special attention to cracks, crevices, seams and hard to reach areas. Refer to the labeling of the detergent used for additional instructions for use.
- 2 Thoroughly rinse the devices under cold, running utility (tap) water.
- 3 Fully immerse the devices in the detergent solution and sonicate for 10 minutes.
- 4 Thoroughly rinse the devices under cold, running utility (tap) water.
- 5 Prepare a bath of Isopropyl Alcohol (70% IPA).
- 6 Immerse the devices in the isopropyl alcohol to remove any soap residue and minerals.
- 7 Dry the devices with a lint-free cloth and allow them to air dry

Sterilization Steps:

- 1 Place the Surgical Curved Scissors, Surgical Tissue Forceps, Round Stainless-Steel Bowl, Rectangular Stainless-Steel Bowl, Dual Biomaterial Carrier Spatula, Dual Biomaterial Packer and reassembled Xpression Box in FDA cleared sterilization bags or wraps.
- 2 Run through one of the following qualified sterilization cycles:

Sterilization Method	Temperature	Exposure Time	Minimum Drying Time
Pre-vacuum Steam (ANSI/AAMI ST79)	132°C (270°F)	4min	20-30 minutes
Pre-vacuum Steam (UK DoH Health Technical Memorandum 01-01)	134°C (273°F)	3min	20-30 minutes

Attention! Improper cleaning may lead to inadequate sterilization. Failure to completely dry the Surgical Curved Scissors, Surgical Tissue Forceps, Round Stainless-Steel Bowl, Rectangular Stainless-Steel Bowl, Dual Biomaterial Carrier Spatula, Dual Biomaterial Packer and Xpression Box components during autoclaving may leave moisture and cause discoloration and oxidation. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the devices. Periodic testing, cleaning, and calibration of the autoclave equipment is recommended to ensure the unit remains in proper working order.

10 CENTRIFUGE CLEANING AND MAINTENANCE (HOUSING, LID AND CENTRIFUGING CHAMBER)

The centrifuge may become contaminated. Pull the mains plug before cleaning. Centrifuges, rotors and accessories must not be cleaned in rinsing machines. They may only be cleaned by hand and disinfected with liquids. The water temperature must be between 20 – 25°C. Only detergents/disinfectants with a pH between 5 – 8 and that do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions may be used. In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant must be considered.

Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required to prevent corrosion through adhering impurities. Ingredients of suitable detergents include soap, anionic surfactants and non-ionic surfactants. After using detergents, remove detergent residue by wiping with a damp cloth. The surfaces must be dried immediately after cleaning. In the event of water

condensation, dry the centrifugal chamber by wiping out with an absorbent cloth. Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning. The centrifuging chamber is to be checked for damage. If damage is found to be relevant to safety, the centrifuge may no longer be put into operation. In this case, notify BioHorizons Customer Service or your distributor.

For surface disinfection, if infectious materials penetrate into the centrifugal chamber, it must be disinfected immediately using a sporicidal surface disinfectant such as ECOLAB's Incidin OxyFoam S or equivalent. After using disinfectants, remove disinfectant residue by wiping with a damp cloth. The surfaces must be dried immediately after disinfecting.

11 ROTOR AND ACCESSORIES CLEANING AND MAINTENANCE

To avoid corrosion and changes to the materials, the rotor and accessories have to be cleaned regularly with soap or a mild cleaning agent and a moist cloth. Cleaning is recommended at least once a week. Contaminants should be removed immediately.

Ingredients of suitable detergents include soap, anionic surfactants and non-ionic surfactants. After using detergents, remove detergent residue with a damp cloth. The rotor and accessories have to be dried immediately after cleaning. Check the rotor and accessories weekly for wear and corrosion damage. The rotor and accessories must no longer be used if they show signs of wear or corrosion. Check the firm seating of the rotor on a weekly basis.

If infectious material should get on the rotor or accessories, they should be appropriately disinfected using a sporicidal surface disinfectant such as ECOLAB's Incidin OxyFoam S or equivalent. After using disinfectants, remove disinfectant residue by wiping off with a damp cloth. The rotor and accessories must be dried directly after disinfection.

The rotor may be autoclaved at 121°C/250°F for 20 minutes and dried appropriately. After 10 autoclaving cycles, the rotor must be exchanged for safety reasons. Autoclaving accelerates the ageing process of plastics and may cause discoloration. After autoclaving, wait until the rotor has cooled down to the ambient temperature before using it again. No statement can be made about the degree of sterility.

The period of use of the rotor is limited to 50,000 running cycles (centrifugation runs) or 5 years, whichever comes first. The maximum permissible number of run cycles can be seen on the rotor. For safety reasons, the rotor may no longer be used when the maximum allowed number of running cycles (marked on it) has been reached. The centrifuge is equipped with a cycle counter which counts the running cycles (centrifugation runs).


In the rare event of blood tube fracture, all broken parts and blood are to be completely removed from the centrifuge. The centrifuge is to be thoroughly cleaned and disinfected as indicated and rubber inserts and plastic sleeves of the rotor are to be replaced.

12 CENTRIFUGE TECHNICAL SPECIFICATIONS

Model Type	IS220Z	IS110Z
Mains voltage ($\pm 10\%$)	200 - 240 V 1~	100 - 127 V 1~
Mains frequency	50 - 60 Hz-	50 - 60 Hz
Connected load	100 VA	100 VA
Current consumption	0.5 A	1.0 A
Capacity	8 x 10 ml	
Maximum allowed density	1.2 kg/dm ³	
Maximum Speed (RPM)	6,000	
Force (RCF)	3,461	
Kinetic energy	750 Nm	
Set-up site	Indoors only	
Altitude	Up to 2000 m above sea level	
Ambient temperature for operation	5°C to 40°C	
Relative Humidity for operation	Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.	
Excess-voltage category	II	
Pollution degree	2	
Device protection class	I	
	Not suitable for use in explosion-endangered areas.	
Emitted interference, Interference immunity	EN / IEC 61326-1, Class B	FCC Class B
Noise level (dependent on rotor)	≤ 50 dB(A)	
Centrifuge width	261 mm	
Centrifuge Depth	353 mm	
Centrifuge Height	228 mm	
Centrifuge Weight	approx. 9 kg	

13 CENTRIFUGE OPERATING INSTRUCTIONS

It is imperative that the transport securing device, consisting of 2 screws (a) and 2 spacers (b), be removed. Keep the transport securing device in a safe place since it must be installed again before transporting the device. The device may only be transported with the transport securing device installed. To protect the device from damage during transport, the motor is fixed in place. This transport securing device must be removed when the device is put into operation.

13.1 INITIAL OPERATION

- Remove the transport securing device from the centrifuge bottom side.
- Position the centrifuge in a stable and level manner in a suitable place. When the centrifuge is running, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- Ventilation openings must not be blocked. A distance of 300 mm must be maintained from the ventilation slots and openings of the centrifuge.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the power cord to a standard mains socket.
- Switch on the mains switch.
- The following displays appear on the panel: the centrifuge model type, the software version, and the last used centrifugation data.
- If the lid is closed, the message "Open the lid" is displayed. In this case, open the lid to display the centrifugation data.

13.2 OPENING AND CLOSING THE CENTRIFUGE LID

The lid can only be opened if the centrifuge is switched on and the rotor is stationary. When the cycle counter is activated, after a centrifugation run, while opening the lid, the remaining number of running cycles (centrifugation runs) is briefly displayed.

Example:

To open the lid, Press the following key . The lid is unlocked by the motor.  indicates lid unlocked.

Example:

Do not reach with your fingers between the lid and housing. Do not slam the lid closed.


To close the lid, lightly press down the front edge of the lid.  indicates lid locked.

Example:

13.3 EMERGENCY UNLOCKING


In the event of a power failure, the lid cannot be unlocked with the motor. Emergency unlocking must be done by hand. To unlock in an emergency, switch off the mains switch (switch setting "0"). Look through the window in the lid to make sure that the rotor is at a standstill. Open the lid only when the rotor is at a standstill. Insert the Allen key horizontally in the bore (A) and turn carefully counterclockwise (to the left) until the lid opens. CAUTION! Turning the hexagon Allen key in clockwise direction (to the right) may damage the locking system. Pull the Allen key back out of the bore.


13.4 INSTALLATION AND REMOVAL OF THE ROTOR

To remove the rotor, loosen the rotor's clamping nut by turning counter-clockwise with the Allen wrench (included with the IntraSpin centrifuge system) and turn up to the lifting pressure point. After overcoming the lifting pressure point, the rotor is released from the cone of the motor shaft. Turn the clamping nut until the rotor can be lifted up from the motor shaft. Lift up the rotor from the motor shaft.


To install the rotor, clean the motor shaft (A) and the bore of the rotor and apply a thin coat of grease to the motor shaft. Dirt particles between the motor shaft and rotor prevent the rotor from having a perfect seat and cause it to run unsteadily. Place the rotor vertically onto the motor shaft. When putting on the rotor, the marking beam (B) on the rotor must be parallel to both surfaces (C) on the motor shaft. Tighten the clamping nut of the rotor with the Allen wrench by turning clockwise. Check the rotor to make sure it is seated firmly.


The rotors must be loaded symmetrically. The blood tubes have to be distributed evenly on all rotor positions.

Rotor evenly loaded


Rotor not evenly loaded
Not permitted!

The blood tubes may only be filled outside of the centrifuge. The maximum filling quantity for the blood tubes is specified on the tubes themselves and must not exceed the maximum weight listed on the rotor. The centrifuging vessels may only be filled so far that no fluid can be expelled from them while the centrifuge is running. When loading the rotor, no liquid may enter the rotor or the centrifuging chamber. In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the blood tubes is recommended.

13.5 DISPLAY ELEMENTS DESCRIPTION

Symbol / Panel Key	Description
	Lid unlocked.
	Lid locked.
	Rotation display. The rotation display lights up, rotating in a clockwise direction, when the rotor is turning.
	To input speed directly. If the key is kept pressed, the value changes with increasing speed.
	To input the runtime directly. Adjustable in steps of 1 second up to a minute, and in steps of 1 minute starting from 1 minute. To input the centrifugation parameters. If the key is kept pressed, the value changes with increasing speed.
	To activate individual parameters. Every time the key is pressed, the next parameter is activated. Keep the key pressed for 8 seconds to call up the "MACHINE MENU". In the "Machine Menu", select the menus "->Info", "->Settings" and "Time & Cycles". To scroll forward in the menus.
	To switch between the speed display (RPM) and relative centrifugal force display (RCF). RCF values are displayed between arrows ><.
	To start the centrifugation run. For short-term centrifugation. Centrifugation is run as long as the key is kept pressed. To select the menus "->Info", "->Settings" and "->Time & Cycles".
	To finish the centrifugation run. The rotor runs down with a pre-selected brake stage. Pressing the key twice triggers the Emergency Stop. To unlock the lid.
t/min	Runtime. Adjustable from 1 - 99 minutes, in steps of 1 minute.
t/sec	Runtime. Adjustable from 1 - 59 seconds, in steps of 1 second. Continuous run "--:--". Set the parameters t/min and t/sec to zero.
RPM	Speed. A number value from 200 RPM to the maximum speed of the rotor can be set. Adjustable in steps of 10.
>RCF<	Relative centrifugal force. A number value can be set which results in a speed between 200 RPM and the maximum rotor speed. Adjustable in steps of 1. It is only possible to input the relative centrifugal force (RCF) if the RCF display (>RCF<) is activated. The relative centrifugal force (RCF) depends on the centrifuging radius (RAD). After entering the RCF, check to make sure that the correct centrifuging radius has been set.
RAD/mm	Centrifuging radius. Adjustable from 10 mm to 250 mm, in steps of 1 mm. It is only possible to input the centrifuging radius if the RCF display >RCF< is activated.
	Brake stage. fast = short run-out time, slow = long run-out time.

13.6 DIRECT INPUT OF THE CENTRIFUGATION PARAMETERS

The speed (RPM), the relative centrifugal force (RCF), the centrifuging radius (RAD) and the runtime can be inputted directly with the keys   without previously having to press the  key. The set centrifugation parameters are only stored after starting the centrifugation run.

For Speed (RPM):

 Press the key to activate the RPM display (RPM) as needed.

RPM/RCF Set the desired value with the keys.

Example:

For Relative Centrifugal Force (RCF) and Centrifugal Radius (RAD):

 Press the key to activate the RCF display (>RCF<) as needed.



RPM/RCF Set the desired RCF value with the keys.

t Set the desired centrifuging radius with the keys as needed.

Example:

Runtime:

Up to 1 minute, the runtime can be set in steps of 1 second, and starting from 1 minute, it can only be set in steps of 1 minute. In order to set the continuous run, the parameters **t/min** and **t/sec** must be set to zero. In the time display (t/min:s), "--:--" appears.


Example:

>RCF<	RPM	t/min:s
≡	4500	--:--

t Set the desired value with the keys.

Example:



13.7 INPUT OF THE CENTRIFUGATION PARAMETERS WITH THE "SELECT" KEY

The runtime can be set in minutes and seconds (parameters **t/min** and **t/sec**). In order to set the continuous run, the parameters **t/min** and **t/sec** must be set to zero. In the time display (t/min:s), "--:--" appears.

The relative centrifugal force (RCF) depends on the centrifuging radius (RAD). During the input of the RCF, the set centrifuging radius is displayed. If no key is pressed for 8 seconds after selection or during parameter input, the previous values are shown on the display. The parameters must then be entered again. By pressing the button, the settings will be saved. If several parameters are entered, the key must be pressed after setting the last parameter. Entering parameters can be cancelled at any time by pressing the key. In this case, the settings are not stored.

13.8 CONTROL PANEL DISPLAY

After each centrifugation run, while the lid is unlocking, the remaining number of running cycles (centrifugation runs) will be briefly displayed.

Example:

If the permissible weight difference was exceeded within the rotor load, the centrifugation run is aborted during start-up and "IMBALANCE" will display.

Example:

A centrifugation run can be aborted at any time by pressing the key. During a centrifugation run, centrifugation parameters can be selected and changed. The changed values, however, only apply to the current centrifugation run and are not saved. With the key, you can toggle between the RPM display (RPM) and RCF display (>RCF<) at any time. If the RCF display (>RCF<) is modified, the centrifuging radius must be entered. Further operation of the centrifuge is then only possible after opening the lid. Operation errors and malfunctions will be shown.

Example:

13.9 CENTRIFUGATION WITH PRESET TIME

 Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed.

↓

Enter the desired centrifugation parameters.

↓

 Press the key to start the centrifugation run. During the centrifugation run, the rotor speed or the resulting RCF value and remaining time are displayed.


↓

After the time elapses or if the centrifugation run is aborted by pressing the  key, the run-out occurs with the set brake stage. The brake stage ~_f is displayed.

Example: (RPM) display

Example: (>RCF<) display

13.10 CONTINUOUS RUN

 Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed.

↓

Input the desired centrifugation parameters. Set the parameters t/min and t/sec to zero.



↓

 Press the key to start the centrifugation run. During the centrifugation run, the rotor speed or the resulting RCF value and elapsed time are displayed.




↓

 Press the key to end the centrifugation run. The run-out occurs with the set brake stage. The brake stage ~_f is displayed.

Example: (RPM) display

Example: (>RCF<) display

13.11 SHORT-TERM CENTRIFUGATION

Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed.

Enter the desired centrifugation parameters.


Press the key and keep it pressed. During the centrifugation run, the rotor speed or the resulting RCF value and elapsed time are displayed.

Release the key again to end the centrifugation run. The run-out occurs with the set brake stage. The brake stage ~_f is displayed.

Example: (RPM) display

Example: (>RCF<) display

13.12 EMERGENCY STOP

Press the key twice. During the Emergency Stop, the run-out occurs with the "fast" brake stage (short run-out time). Brake stage ~_f is displayed.

Example: (RPM) display

Example: (>RCF<) display

13.13 CYCLE COUNTER

The centrifuge is equipped with a cycle counter, which counts the running cycles (centrifugation runs). After each centrifugation run, while the lid is unlocking, the remaining number of running cycles (centrifugation runs) will be briefly displayed.


If the maximum permissible number of rotor running cycles has been exceeded, the following is displayed every time the centrifugation run is started. The centrifugation run must be restarted. The rotor must be immediately exchanged for a new rotor for safety reasons. After the rotor has been exchanged, the cycle counter must be reset to "0".

13.14 RESETTING THE CYCLE COUNTER TO ZERO

After the rotor has been exchanged, the cycle counter must be reset to zero again. The cycle counter may only be reset to zero if the rotor has been exchanged for a new rotor. With the rotor at a standstill, the cycle counter can be reset as shown below. The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.


Example:

Example:

13.15 SETTINGS AND QUERIES


13.16 QUERYING SYSTEM INFORMATION

The following system information can be queried:

- Device model
- Program version of the device
- Type number of the device
- Device date of manufacture
- Device serial number
- Frequency converter type
- Frequency converter program version

With the rotor at a standstill, the query can proceed as follows: The procedure can be aborted at any time by pressing the  key.


Press the key.
Type#1: Type number of the device.

Example:


Press the key.
Type#2: Continuation of the type number if this has more positions than the "Type#1" display field.

Example:

Press the key.
 Device date of manufacture.

Example:

Press the key.
 Device serial number.

Example:

Press the key.
 Frequency converter type.

Example:

Press the key.
 Frequency converter program version.


Example:

Press the key twice to exit the "-> Info" menu or press three times to exit the "* MACHINE MENU *".

Example:

13.17 ACOUSTIC SIGNAL

The acoustic signal sounds (1) after a malfunction occurs in a 2-second interval and (2) after ending the centrifugation run and the rotor is at a standstill in a 30-second interval. The acoustic signal is ended by pressing any key. If the rotor is at a standstill, the acoustic signal can be set as follows. The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.


Press and hold the button for eight seconds.


Press the key until the following is displayed.


13.18 OPTICAL SIGNAL AFTER ENDING THE CENTRIFUGATION RUN


The backlighting of the display flashes after the centrifugation run to visually signalize that the centrifugation run has finished. The optical signal can be switched on or off when the rotor is at a standstill. The procedure can be aborted at any time by pressing the key.

13.19 AUTOMATIC UNLOCKING OF THE LID AFTER THE CENTRIFUGATION RUN


It can be set whether the lid should automatically unlock or not after the centrifugation run. With the rotor at a standstill, this can be set as follows. The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.


13.20 BACKLIGHTING OF THE DISPLAY

To save energy, it can be set that, after a centrifugation run, the backlighting of the display switches off after 2 minutes. With the rotor at a standstill, this can be set as follows. The procedure can be aborted at any time by pressing the **STOP OPEN** key. In this case, the settings are not stored.


With the keys, set **off** or **on**.

Press the button to save the setting.

Press the key once to exit the "-> Settings" menu or press twice to exit the "* MACHINE MENU *".

Example:

13.21 QUERYING THE HOURS OF OPERATION AND THE NUMBER OF CENTRIFUGATION RUNS

The hours of operation are divided into internal and external hours of operation. Internal hours of operation: Total time the device was switched on. External hours of operation: Total time of the previous centrifugation runs. With the rotor at a standstill, the query can proceed as follows. The procedure can be aborted at any time by pressing the  key.

Press and hold the button for eight seconds.

Press the key until the following is displayed.

Press the key.
TimeExt: External hours of operation.

Example:


Press the key.
TimeInt: Internal hours of operation.

Example:

Press the key.
Starts: Number of all centrifugation runs.


Example:

Press the key twice to exit the "-> Time & Cycles" menu or press three times to exit the "* MACHINE MENU *".

Example:

14 FAULTS

If the fault cannot be eliminated with the help of the fault table below, please contact BioHorizons Customer Service or your local distributor and specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

To perform a MAINS RESET: Switch off the mains switch (switch position "0"). Wait at least 10 seconds and then switch on the mains switch again (switch position "I").

Message / fault		Cause	Remedy
No display		No voltage Mains input fuses defective.	<ul style="list-style-type: none"> – Check distribution voltage. – Check mains power input fuse and refer to instructions for changing mains fuse. – Mains switch ON.
IMBALANCE		The rotor is unevenly loaded.	<ul style="list-style-type: none"> – Open the lid after the rotor is at a standstill. – Check the loading of the rotor. – Repeat the centrifugation run.
MAINS INTER	11	Power failure during the centrifugation run. (The centrifugation run was not finished.)	<ul style="list-style-type: none"> – Open the lid after the rotor is at a standstill. – Push start/pause button. – Repeat the centrifugation run if necessary.
MAINS INTERRUPT			
TACHO ERROR	1, 2	Failure of speed impulses during operation.	<ul style="list-style-type: none"> – Perform a MAINS RESET when the rotor has been stationary.
LID ERROR	4.1 – 4.127	Error in lid locking or lid closure.	
OVER SPEED	5	Rotation too fast	
VERSION ERROR	12	Incorrect centrifuge model recognized. Error / defect electronics	<ul style="list-style-type: none"> – Perform a MAINS RESET when the rotor has been stationary.
UNDER SPEED	13	Rotation too slow	
CTRL ERROR	22.1, 25.2	Error / defect of electronics	
CRC ERROR	27.1		
COM ERROR	31 – 36		
FC ERROR	60, 61.1 – 61.21, 61.64 – 61.142		
FC ERROR	61.23	Speed measurement error	<ul style="list-style-type: none"> – The device may not be switched off as long as the rotation display  is lit up and rotating. Wait until the symbol  (lid locked) is displayed (after approx. 120 seconds). Afterwards, carry out a MAINS RESET.
TACHO ERR	61.22		
FC ERROR	61.153	Error / defect of electronics	<ul style="list-style-type: none"> – Perform a MAINS RESET. – Check the loading of the rotor per the instructions provided. – Repeat the centrifugation run.

15 RELATIVE CENTRIFUGAL FORCE (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance. The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius. It can be calculated using the formula below:

$$RCF = \left(\frac{RPM}{1000}\right)^2 \times r \times 1.118 \Rightarrow RPM = \sqrt{\frac{RCF}{r \times 1.118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the center of the turning axis to the bottom of the centrifuge container.

16 TRANSPORT AND STORAGE


Before transporting the device, the transport securing device must be installed. The device and the accessories may only be stored in closed and dry rooms. When the device and accessories are transported or stored, the following conditions must be complied with:

Temperature limits: -20°C to +60°C.

Relative Humidity limits: 10% to 80%, non-condensing.

17 CHANGE MAINS INPUT FUSES

Start by switching off the mains switch and separating the device from the mains. The fuse holder (A) with the mains input fuses is located next to the mains switch. Remove the connecting cable from the machine plug socket. Press the snap-fit (B) against the fuse holder (A) and remove. Exchange defective mains input fuses. Only use fuses with the rating defined for the centrifuge type provided by the manufacturer. Reinsert the fuse holder until the snap-fit clicks shut. Reconnect the device to the mains supply.

18 RETURNING DEVICES

Before returning the device, the transport securing device has to be installed. If the device or its accessories are returned to legal manufacturer, in order to provide protection for people, the environment and materials, it has to be cleaned and disinfected before being shipped. The legal manufacturer reserves the right to refuse contaminated devices or accessories. Costs incurred for cleaning and disinfection are to be charged to the customer.

19 DISPOSAL

Products to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant regulations.

Before disposal, the IntraSpin centrifuge must be decontaminated and cleaned to protect people, the environment and property. When disposing of the device, the respective statutory rules must be observed. Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed of as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field. The icon of the crossed-out trashcan shows that the device may not be disposed of as domestic waste. The waste disposal guidelines of individual countries might vary. If necessary, contact your supplier.

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例：The image shows a rectangular LCD display with a dark background and light text. At the top, there are three labels: '>RCF<', 'RPM', and 't/min:s'. Below these labels, the text 'RemCycles= 16703' is displayed in a large, bold font.

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例：The image shows a rectangular LCD display with a dark background and light text. At the top, there are three labels: '>RCF<', 'RPM', and 't/min:s'. Below these labels, there is a small icon of an open lid, followed by the text '4500' and '5:00'.

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例：The image shows a rectangular LCD display with a dark background and light text. At the top, there are three labels: '>RCF<', 'RPM', and 't/min:s'. Below these labels, there is a small icon of a closed lid, followed by the text '4500' and '5:00'.

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例：The image shows a control panel display with three indicators at the top: '>RCF<', 'RPM', and 't/min:s'. Below these, the text 'RemCycles= 16703' is displayed.

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例：The image shows a control panel display with three indicators at the top: '>RCF<', 'RPM', and 't/min:s'. Below these, the text 'IMBALANCE' is displayed.

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例：The image shows a control panel display with three indicators at the top: '>RCF<', 'RPM', and 't/min:s'. Below these, the text 'Open the lid' is displayed.

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예: The image shows a control panel display with three indicators at the top: '>RCF<', 'RPM', and 't/min:s'. Below these, the text 'RemCycles= 16703' is displayed.

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예: The image shows a control panel display with three indicators at the top: '>RCF<', 'RPM', and 't/min:s'. Below these, the text 'IMBALANCE' is displayed.

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예: The image shows a control panel display with three indicators at the top: '>RCF<', 'RPM', and 't/min:s'. Below these, the text 'Open the lid' is displayed.

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