English	Instructions for use: Dental Implants
Español	Instrucciones de uso: implantes dentales
русский	Инструкция по применению: Зубные имплантаты
Deutsch	Gebrauchsanweisung: Dentalimplantate
Français	Instructions d'utilisation: Implants Dentaires
Italiano	Istruzioni per l'uso: Impianti Dentali
Português	Instruções de utilização: Implantes Dentários
Türkçe	Kullanım talimatları: Dental İmplantlar
简体中文	使用说明: 牙种植体
日本語	使用説明書: 歯科インプラント
한국어	적응증: 치과 시술용 임플란트
عربي	إرشادات الاستخدام: غرسات الأسنان
Polski	Instrukcja obsługi: Implanty dentystyczne
Українська	Інструкція із застосування: зубні імплантати
Česky	Návod k použití: Dentální implantáty
Български	Инструкции за употреба: Дентални импланти
Ελληνικά	Οδηγίες χρήσης: Οδοντικά εμφυτεύματα
Magyar	Használati utasítás: Fogászati implantátumok
Română	Instrucțiuni de utilizare: implanturi dentare
Nederlands	Gebruiksaanwijzing: Tandheelkundige implantaten
Svenska	Bruksanvisning: Tandimplantat
Norsk	Bruksanvisning: Dentalimplantater

**BioHorizons Implant Systems Inc.**

2300 Riverchase Center  
 Birmingham, AL 35244 USA  
 TOLL-FREE 888.246.8338  
 TEL 205.967.7880  
 FAX 205.870.0304  
[www.biohorizons.com](http://www.biohorizons.com)

L01079 Rev D JAN 2025

This document supersedes all prior revisions. Original language is English.

This document applies to all Intra-Lock (Fusion) Dental Implants. Every product packaging label contains a description of the enclosed device.

## **DESCRIPTION**

Intra-Lock implants are manufactured from biocompatible titanium alloy (Ti-6AL-4V ELI) and are single use only.

<b>Element</b>	<b>Composition % (mass/mass)</b>
Nitrogen, max	0.05
Carbon, max	0.08
Hydrogen, max	0.012
Iron, max	0.25
Oxygen, max	0.13
Aluminum	5.5 – 6.50
Vanadium	3.5 – 4.5
Titanium	balance

## **INTENDED USE**

Intra-Lock implants are intended for use in the mandible or maxilla for use as an artificial root structure for single tooth replacement or for fixed bridgework and dental retention.

## **INDICATIONS FOR USE**

Intra-Lock implants are intended for use in the mandible or maxilla for use as an artificial root structure for single tooth replacement or for fixed bridgework and dental retention. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion or (2) when splinted together for multiple tooth replacement or when stabilized with an overdenture supported by multiple implants.

Intra-Lock 3.3mm diameter implants may be used as an artificial root structure for single tooth replacement of mandibular central and lateral incisors and maxillary lateral incisors. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion, (2) when splinted together as an artificial root structure for multiple tooth replacement of mandibular incisors, or (3) for denture stabilization using multiple implants in the anterior mandible and maxilla. The implants may be placed in immediate function when good primary stability has been achieved and with appropriate occlusal loading.

## **CONTRAINDICATIONS**

Intra-Lock implants should not be used in patients who have contraindicating systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity.

Intra-Lock 3.3mm diameter implants are contraindicated for single tooth replacement of maxillary central incisors.

## **PATIENT POPULATION**

BioHorizons dental implant systems are intended for use in skeletally mature, non-pediatric edentulous or partially edentulous patients as long as the defined contraindications are nonapplicable.

## INTENDED USERS

BioHorizons implant systems are only intended for use by licensed health care professionals, more specifically, BioHorizons implant systems are intended to be used by trained dentists, surgeons and dental technicians in a standard dental surgical setting, which may range from General Dentists' practice offices to Maxillofacial Surgical Operating Rooms, as well as laboratories for dental processes. Use of these products requires specialized knowledge and experience in implant dentistry. BioHorizons implant systems are marked and labeled as a medical device (MD) and are Rx only.

## DIRECTIONS FOR USE

Proper surgical procedures and techniques are the responsibility of the medical professional. Each dentist must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of postgraduate dental implant education and strict adherence to the instructions and procedures for use that accompany Intra-Lock implant products. NOTE: A pre-operative 30-second rinse with a 0.12% Chlorhexidine Digluconate solution is recommended. (The influence of 0.12% Chlorhexidine Digluconate Rinses on the Incidence of Infectious Complications and Implant Success. Lambert, Paul M., et al, J Oral Maxillofacial Surgery 55 :25-30, 1997, Suppl5). IMPORTANT! (1) The implant procedure manual should be reviewed prior to attempting placement. (2) In dense (D1- D2) bone, further preparation of the site may be required to avoid applying excessive torque to the implant connection. (3) To further prepare the site, the crestal bone drills may be used. (4) The clinician should verify that the proper instruments are available and at hand prior to attempting placement.

## WARNINGS AND PRECAUTIONS

Appropriate training in proper implant surgery technique is strongly recommended prior to implant use. Improper technique can result in implant failure and/or loss of supporting bone. Appropriate x-ray films and/or CT scans should be utilized to determine (1) if adequate bone width and depth are available at the desired implant site and (2) the location of important anatomical landmarks, such as the mandibular canal, maxillary sinuses, and adjacent teeth.

Clinician judgment should be used in determining the minimum post-implantation time before placing the implants in occlusal function. An adequate number of implants should be used to provide support and to distribute the load to the abutments. For short implants, clinicians should closely monitor patients for any of the following conditions: periimplant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting the short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration avoid immediate loading.

Devices should not be used in patients with a known allergy or sensitivity to the device material.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. A certain percentage of implants may fail to achieve or to maintain osseointegration, as demonstrated by mobility, and should be removed. Potential causes of abutment fracture include, but are not limited to: casting beyond a 30-degree angle, casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. Small diameter implants with either angled or straight abutments are intended for the anterior region of the mouth and not intended for the posterior region of the mouth due to possible failure of the implant.

Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at [www.intra-lock.com](http://www.intra-lock.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

For implantable devices, the summary of safety and clinical performance (SSCP) according to Article 32 of Regulation (EU) 2017/745 can be found in the European database for medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed>.

Implantable Device	Basic UDI-DI Number
Intra-Lock Fusion Implants	08472360IIBWETIMPL007PE

## MAGNETIC RESONANCE (MR) SAFETY INFORMATION

### MR Conditional

Intra-Lock implant systems have been demonstrated through non-clinical testing to be magnetic resonance (MR) conditional. A patient with this device can be scanned safely in an MR system under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Static magnetic field of 1.5-Tesla or 3-Tesla, only;
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40-T/m);
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode.

Under the scan conditions defined, the implant system is expected to produce a maximum temperature rise of 3.6°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant system extends approximately 30-mm from this system when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

Parameter	Condition
Nominal Values of Static Magnetic Field (T)	1.5-Tesla and 3.0-Tesla
Maximum Spatial Field Gradient (T/m and gauss/cm)	40-T/m (4,000-gauss/cm)
Type of RF Excitation	Circularly Polarized (CP) (i.e., Quadrature-Transmission)
Transmit RF Coil Information	There are no transmit RF coil restrictions. Accordingly, the following may be used: body transmit RF coil and all other RF coil combinations (i.e., body RF coil combined with any receive-only RF coil, transmit/receive head RF coil, transmit/receive knee RF coil, etc.)
Operating Mode of MR System	Normal Operating Mode
Maximum Whole Body Averaged SAR	2-W/kg (Normal Operating Mode)
Limits on Scan Duration	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks).
MR Image Artifact	The presence of this implant system produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant system is located in the area of interest.

To allow medical professionals to identify the specific medical devices a patient has, the MRI safety status of the medical devices, and the conditions for safe use in the MR environment for MR Conditional devices, BioHorizons recommends that clinicians provide the patient with the dental implant(s) and dental abutment(s) / bar(s) device specific peel-offs affixed to the patient card. Patient cards are available free of charge upon request from BioHorizons or for direct printing at <http://eifu.intra-lock.com>.

### COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

### HANDLING AND STERILIZATION

This product is supplied sterile and should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove from the package only after the correct size has been determined and the surgical site has been prepared. Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. Implants are single

use only, and re-use should not be attempted. Following this guidance eliminates the risk of cross-patient contamination from secondary use of this device. BioHorizons assumes no responsibility for attempted re-use or re-sterilization.

Products to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

### SYMBOLS AND DESCRIPTIONS

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic instructions for use.
	Manufacturer.
	Intra-Lock products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number.
	Lot/ batch number.
	Unique Device Identifier.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.
	Date of manufacture.
<b>Rx Only</b>	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician.

Symbol	Symbol Description
	European Union Authorized Representative.
	Do not use if package is damaged. Discard device and package.
	Medical Device.
	Non-sterile.
	Single sterile barrier system with protective packaging outside.
	Single sterile barrier system.
	Home.
	Magnetic resonance warning: Device is MR conditional.
	United Kingdom Responsible Person

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