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L01081CAen Rev A FEB 2023

This document supersedes all prior revisions. Original language is English.

This document applies to Intra-Lock Mini Drive-Lock (MDL) Dental Implants.

DESCRIPTION

Intra-Lock MDL implants are manufactured from biocompatible titanium alloy and are single use only.

INDICATIONS MINI DRIVE-LOCK™ DENTAL IMPLANT SYSTEM

Mini Drive-Lock™ Dental Implants are intended for use as a self-tapping titanium screw for transitional or intra-bony long-term applications. Mini Drive-Lock™ Dental Implants are indicated for long-term maxillary and mandibular tissue-supported denture stabilization. Multiple implants should be used and may be restored after a period of time or placed in immediate function.

CONTRAINDICATIONS

Patients with alcohol addiction or psychiatric disorders, blood dyscrasias, uncontrolled diabetes, hyperthyroidism, oral infections, malignancies, or patients who have had myocardial infarction within the last 12 months.

Patients with systemic diseases that compromise the immune system, such as AIDS, patients on medications that would compromise healing of an implant site, patients with a history of poor or non-compliance to oral hygiene procedures, or patients who cannot maintain oral hygiene procedures if implants are placed.

Tobacco usage increases the occurrence of complications and failures.

DIRECTIONS FOR USE

Proper surgical procedures and techniques are the responsibility of the medical professional. Each dentist must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of postgraduate dental implant education and strict adherence to the instructions and procedures for use that accompany Intra-Lock implant products. NOTE: A pre-operative 30-second rinse with a 0.12% Chlorhexidine Digluconate solution is recommended. (The influence of 0.12% Chlorhexidine Digluconate Rinses on the Incidence of Infectious Complications and Implant Success. Lambert, Paul M., et al, J Oral Maxillofacial Surgery 55:25-30, 1997, Suppl5).

IMPORTANT! (1) The implant procedure manual should be reviewed prior to attempting placement. (2) The clinician should verify that the proper instruments are available and at hand prior to attempting placement. (3) MDL implants are not recommended to be placed in sites with a cortical plate thickness greater than 3mm.

WARNINGS AND PRECAUTIONS

Appropriate training in proper implant surgery technique is strongly recommended prior to implant use. Improper technique can result in implant failure and/or loss of supporting bone. Appropriate x-ray films and/or CT scans should be utilized to determine (1) if adequate bone width and depth are available at the desired implant site and (2) the location of important anatomical landmarks, such as the mandibular canal, maxillary sinuses, and adjacent teeth.

Clinician judgment should be used in determining the minimum post-implantation time before placing the implants in occlusal function. An adequate number of implants should be used to provide support and to distribute the load to the abutments.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. A certain percentage of implants may fail to achieve or to

maintain osseointegration, as demonstrated by mobility, and should be removed. Potential causes of abutment fracture include, but not limited to: casting beyond a 15-degree angle, casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontally compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics.

Devices should not be used in patients with a known allergy or sensitivity to the device material(s).

Intra-Lock dental implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of Intra-Lock dental implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at www.intra-lock.com. Instructions for Use are available free of charge in printed form upon request from BioHorizons or your local representative / distributor. Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported immediately to BioHorizons Customer Care at 1-866-468-8338.

COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgment; (3) infection requiring revision of the dental implant; (4) nerve damage that could cause permanent weakness, numbness, or pain; (5) histologic responses possibly involving macrophages and/or fibroblasts; (6) formation of fat emboli; (7) loosening of the implant requiring revision surgery; (8) perforation of the maxillary sinus; (9) perforation of the labial and lingual plates; and (10) bone loss possibly resulting in revision or removal.

HANDLING AND STERILIZATION

This product is supplied sterile, using gamma irradiation, and should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove from the package only after the correct size has been determined and the surgical site has been prepared. Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. Implants are single use only, and re-use should not be attempted. Following this guidance eliminates the risk of cross-patient contamination from secondary use of this device. BioHorizons assumes no responsibility for attempted re-use or re-sterilization.


Products to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant regulations.

SYMBOLS AND DESCRIPTIONS

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic instructions for use.
	Manufacturer.

Symbol	Symbol Description
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number.
	Lot/ batch number.
	Unique Device Identifier.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.
	Date of manufacture.
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician.
	European Union Authorized Representative.
	Do not use if package is damaged. Discard device and package.
	Medical Device.
	Non-sterile.
	Single sterile barrier system with protective packaging outside.

Symbol	Symbol Description
	Single sterile barrier system.
	Home.
	United Kingdom Responsible Person.

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