English	Instructions for Use: MDL/MILO Dental Instruments
Español	Instrucciones de uso: MDL/MILO Instrumentos Dentales
Français	Instructions d'utilisation : MDL/MILO Instruments Dentaires

BioHorizons Implant Systems Inc.

2300 Riverchase Center
Birmingham, AL 35244 USA
TOLL FREE 888.246.8338
TEL 205.967.7880
FAX 205.870.0304

www.biohorizons.com

L02071 Rev C JUL 2022

English

This document supersedes all prior revisions. Original language is English.

DESCRIPTION

Intra-Lock MDL and MILO instruments are used for dental implant procedures such as site development, implant placement, and implant restorations within the specific indications of each implant system. MDL and MILO dental instruments are manufactured from Stainless Steel and Titanium Alloy (Ti-6Al-4V ELI). The label on each instrument contains important product information including whether the instrument is supplied sterile or non-sterile and if the instrument is for single patient use only. Instruments and kits supplied non-sterile must be cleaned and sterilized prior to first use and each use thereafter.

INDICATIONS FOR USE

Intra-Lock MDL and MILO instruments and surgical kit are indicated for use in site development, implant placement, and implant restorations within the specific indications of each implant system.

Intra-Lock MDL Single Patient Drills are provided sterile (by gamma radiation) and are indicated for use with the Intra-Lock MDL implants only. No attempt should be made to re-use or re-sterilize the drills.

CONTRAINDICATIONS

Intra-Lock MDL and MILO instruments should not be used with patients that have allergies to the specific materials used, including stainless steel and titanium alloy.

DIRECTIONS FOR USE

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the instructions pertaining to Intra-Lock products.

Drills should be replaced when wear is noticed, such as a decrease in cutting efficiency or when signs of discoloration appear. Cutting instruments should be replaced after approximately 12 to 20 osteotomy cycles, depending on bone density. Single Patient Drills: Insert the drill into a handpiece and twist protective tubing off prior to use.

IMPORTANT! The implant procedure manual should be reviewed prior to attempting placement for important drilling speed and technique recommendations.

WARNINGS AND PRECAUTIONS

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any Intra-Lock Instructions for Use (IFU). Refer to OSHA standard 29CFR1910.1030 prior to cleaning and sterilization. Additional technical information is available upon request from BioHorizons or may be viewed and/or downloaded at www.intra-lock.com. Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU.

Devices should not be used in patients with a known allergy or sensitivity to the device material(s).

COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with MDL and MILO instruments, prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant, instrument and/or abutment material; (2) implant, instrument and/or abutment breakage; (3) infection requiring revision of the dental implant; (4) nerve damage that could cause permanent weakness, numbness, or pain; (5) histologic responses possibly involving macrophages and/or fibroblasts; (6) formation of fat emboli; (7) implant loosening requiring revision surgery; (8) maxillary sinus perforation; (9) labial or lingual plate perforation; and (10) bone loss possibly resulting in revision or removal.

HANDLING AND STERILIZATION

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. If the product is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove product from the package only after the correct size has been determined. If the product is supplied non-sterile and requires sterilization prior to initial use, remove and discard any shipping material before sterilization. Reprocessing of reusable devices must include cleaning and sterilization. Disassemble the surgical kits before each cleaning. BioHorizons devices have not been validated for automated cleaning. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy’s Enzymax® per the manufacturer’s recommendations.
- 2) Brush the product to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse product under running utility tap water.
- 4) Place product in the sterile container filled with the prepared detergent solution and sonicate for two (2) minutes minimum.
- 5) Thoroughly rinse product under running utility tap water.
- 6) Spray or wipe product with 70% IPA.
- 7) Blot product dry with clean lint free cloth.

If applicable, return the instruments to the appropriate locations in the surgical tray. For sterilization, place the cleaned device in an approved sterilization bag or wrap and run through one of the following qualified sterilization cycles:


Sterilization Cycles				
Reference:	AAMI TIR12:2020	AAMI TIR12:2020	AAMI TIR12:2020	UK HTM 01-01 Part C:2016
Type:	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
Exposure Time and Temperature:	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
Minimum Dry Time:	30 minutes	30 minutes	20 minutes	20 minutes

It is recommended to include a 30-minute cool-down period before removing the product from the sterilization bag or wrap.

Attention! Improper cleaning may lead to inadequate sterilization. Failure to completely dry instruments during autoclaving may leave moisture and cause discoloration and oxidation. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the instruments. Periodic testing, cleaning, and calibration of the autoclave equipment is recommended to ensure the unit remains in proper working order.

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic Instructions for Use.
	Manufacturer.
	Intra-Lock products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Devices Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/article number.
	Lot/batch number.
	Unique Device Identifier.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.
	Date of manufacture.
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution, and use by, or on the order of, a dentist or physician.
	European Union Authorized Representative.

Symbol	Symbol Description
	Do not use if package is damaged. Discard device and package.
	Medical Device.
	Non-sterile.
	Home position.
 MR Conditional	Magnetic resonance warning: Device is MR conditional.
	United Kingdom Responsible Person.

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