English	Instructions for use: Dental Prosthetics
Español	Instrucciones de uso: Prótesis Dentales
русский	Инструкция по применению: Зубные протезы
Deutsch	Gebrauchsanweisung: Dentale Prothetik
Français	Instructions d'utilisation: Prothèses Dentaires
Italiano	Istruzioni per l'uso: Protesi Dentali
Português	Instruções de utilização: Próteses Dentárias
Türkçe	Kullanım talimatları: Dental Prostetik Ürünler
简体中文	使用说明: 牙科义齿
日本語	使用説明書: 歯科プロテーゼ
한국어	적응증: 치과 보철
عربي	إرشادات الاستخدام: البدليات السنية
Polski	Instrukcja obsługi: Protetyka Dentystyczna
Українська	Інструкції з використання: зубні протези
Čeština	Návod k použití: Zubní protetické komponenty
Български	Инструкции за употреба: Дентални протези
Ελληνικά	Οδηγίες χρήσης: Οδοντικά προσθετικά εξαρτήματα
Magyar	Használati utasítás: Fogászati protézisek
Română	Instrucțiuni de utilizare: Protetică Dentară
Nederlands	Gebruiksaanwijzing: tandheelkundige prothesen
Svenska	Bruksanvisning: Tandprotes
Norsk	Bruksanvisning: Tannproteser

**BioHorizons Implant Systems Inc.**

2300 Riverchase Center  
Birmingham AL, 35244 USA  
TOLL-FREE 888.246.8338  
TEL 205.967.7880  
FAX 205.870.0304  
www.biohorizons.com

L02035 Rev J OCT 2023

This document supersedes all prior revisions. Original language is English.

This document applies to BioHorizons Dental Abutments and Prosthetic components listed below. Every product packaging label contains a description of the enclosed device.

Devices in Scope	
BioHorizons healing caps and abutments	BioHorizons Multi-unit abutments for CONELOG
BioHorizons titanium temporary abutments provided sterile	BioHorizons abutment screws
BioHorizons titanium temporary abutments provided non-sterile	BioHorizons abutment attachments
BioHorizons PEEK temporary abutments	BioHorizons abutment copings
BioHorizons gold abutments	BioHorizons prosthetic screws
BioHorizons plastic abutments	BioHorizons impression components (conventional and digital)
BioHorizons titanium abutments *	BioHorizons analogs and replicas
BioHorizons detachable abutments (e.g., ODSecure)	BioHorizons Simple Solution system
BioHorizons detachable abutments for CONELOG	BioHorizons external abutments and prosthetics
BioHorizons Multi-unit abutments	BioHorizons ceramic abutments

\* For CAD/CAM titanium blanks and CAD/CAM titanium bases, refer to Custom Abutments Instructions for Use.

**DESCRIPTION**

BioHorizons prosthetic components are intended for the restoration of BioHorizons dental implants within the specific indications of each implant system. Refer to the following table for material(s) of implantable abutments:

Implantable Devices	Material (main elements)
BioHorizons healing caps and abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons titanium temporary abutments provided sterile	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons titanium temporary abutments provided non-sterile	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons PEEK temporary abutments	TecaPeek Classix (Polyetheretherketone, Carbon Black)
BioHorizons gold abutments	Gold 6019 (Gold, Palladium, Platinum, Iridium)
BioHorizons plastic abutments	Casting material as determined by clinician.
BioHorizons titanium abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons detachable abutments (e.g., ODSecure)	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons detachable abutments for CONELOG	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons Multi-unit abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons Multi-unit abutments for CONELOG	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons abutment screws	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons Simple Solution system	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons external abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons ceramic abutments	ZrO <sub>2</sub> Y-TZP (Ytria Stabilized Zirconia)

BioHorizons Multi-unit Abutments for CONELOG® and CONELOG ODSecure Abutments are intended for the restoration of CONELOG dental implants within the specific indications of each implant system.

The label on each prosthetic package contains important product information including whether the prosthetic is supplied sterile or non-sterile. Prosthetic connection color coding is defined in the following tables:

<b>Internal Hexed Prosthetic Color Coding</b>	Gray	Yellow	Green	Blue
<b>Internal Hexed Prosthetic Size</b>	3.0mm	3.5mm	4.5mm	5.7mm

<b>CONELOG Prosthetic Color Coding</b>	Gray	Yellow	Red	Blue
<b>CONELOG Prosthetic Size</b>	3.3mm	3.8mm	4.3mm	5.0mm

### INDICATIONS FOR USE

BioHorizons Prosthetics are abutments that include healing abutments for contouring tissue and final restorative abutments to support a prosthesis. The abutments may be used for a single or multiple unit restoration and are compatible for use with BioHorizons Internal and Tapered Internal implant systems and Zimmer® Dental Screw-Vent® and Tapered Screw-Vent® implants with 3.5mm, 4.5mm and 5.7mm internal hex-connection mating platform diameters.

BioHorizons Titanium Base Abutments and Laser-Lok Titanium Base Abutments are intended to be used as straight abutments.

The BioHorizons Multi-unit Abutments for CONELOG® are intended to function in the mandible or maxilla to support single and multiple-unit temporary or definitive restorations on CONELOG dental implants.

### CONTRAINDICATIONS

BioHorizons prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity. BioHorizons Multi-unit and Abutment for Screw prosthetics are contraindicated for single tooth replacements.

### DIRECTIONS FOR USE

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions For Use that accompany BioHorizons products. Healing abutments should be hand tightened (10-15 Ncm). Abutment screws must be torqued to 30 Ncm to prevent screw loosening.

<b>Platform</b>	<b>Zimmer Tapered ScrewVent Implants*</b>	<b>Zimmer ScrewVent Implants*</b>
3.5mm Internal Connection	TSVBx TSV4Bx TSVT4Bx	TSVHx TSV4Hx SVMBx SVBx SVMHx
4.5mm Internal Connection	TSVWBx TSVTWBx TSVWHx	SVWBx SVWHx
5.7mm Internal Connection	TSV6Bx TSVT6Bx TSV6Hx	N/A

\*Where variable x = implant length.

Secure Custom Cast Abutments to the supplied polishing protector to protect the implant/abutment interface during laboratory procedures.

Tissue Level Titanium Abutments for Laser Weld should be secured to the working model by torquing the accompanying screw to 30Ncm prior to laser welding. Secure Laser-Lok abutments onto the Laser-Lok Protective Sleeve to protect Laser-Lok during handling.

BioHorizons Multi-unit Abutment screws must be torqued to 30 Ncm.

The BioHorizons Multi-unit straight and angled abutment screws for CONELOG are packaged with the corresponding abutments and must be torqued to 20Ncm for intended function.

The BioHorizons Multi-unit Prosthetic Component screws must be torqued to 15Ncm.

BioHorizons Multi-unit Prosthetic Component design parameters are as follows:

- Minimum post height shall be equal to or greater than 4mm.
- Reduction of wall thickness is not permitted.
- Angle correction is not permitted. Angled Multi-unit abutments shall be used for angle correction.

BioHorizons Multi-unit Prosthetic Components used with the Multi-unit Abutments for CONELOG are intended to support multiple-unit temporary or definitive restorations except for the Titanium Coping. The Titanium Coping used with the Multi-unit Abutments for CONELOG is intended to support single and multiple-unit temporary or definitive restorations.

## WARNINGS AND PRECAUTIONS

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of BioHorizons prosthetic components.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include, but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the modifier of the implant/abutment interface a medical device company subject to FDA rules and regulations. Leaving bonding agent or particulate in the Laser-Lok surface may affect the connective tissue attachment. Use of the LIGHT BLUE abutment screws (packaged with Angled Multi-unit Abutments) with incompatible abutments may result in abutment and/or abutment screw failure.

Small diameter implants and angled abutments are not recommended for use in the posterior region of the mouth.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

Devices should not be used in patients with a known allergy or sensitivity to the device material(s).

BioHorizons Multi-unit Abutments and Multi-unit Abutments for CONELOG require the addition of a coping or cylinder for single and multiple-unit restorations. The coping or cylinder gingival collar, angulation, and wall thickness shall not be modified. The coping or cylinder post height shall not be less than 4mm for single-unit restorations.

Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at [www.biohorizons.com](http://www.biohorizons.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

For implantable devices, the summary of safety and clinical performance (SSCP) according to Article 32 of Regulation (EU) 2017/745 can be found in the European database for medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed>.

Implantable Device	Basic UDI-DI Number
BioHorizons healing caps and abutments	08472360IIBWETABUT001HS

Implantable Device	Basic UDI-DI Number
BioHorizons Titanium Temporary abutments provided sterile	08472360IIBWETABUT002HU
BioHorizons Titanium Temporary abutments provided non-sterile	08472360IIBWETABUT003HW
BioHorizons PEEK temporary abutments	08472360IIBWETABUT004HY
BioHorizons gold abutments	08472360IIBWETABUT005J2
BioHorizons plastic abutments	08472360IIBWETABUT006J4
BioHorizons Titanium abutments	08472360IIBWETABUT007J6
BioHorizons detachable abutments (e.g., ODSecure)	08472360IIBWETABUT008J8
BioHorizons detachable abutments for CONELOG	08472360IIBWETABUT012HX
BioHorizons Multi-unit abutments	08472360IIBWETABUT009JA
BioHorizons Multi-unit abutments for CONELOG	08472360IIBWETABUT011HV
BioHorizons abutment screws	08472360IIBWETABUT010HT
BioHorizons Simple Solution system	Not applicable.
BioHorizons external abutments	Not applicable.
BioHorizons ceramic abutments	Not applicable.

## MAGNETIC RESONANCE (MR) SAFETY INFORMATION

### MR Conditional

BioHorizons implant systems have been demonstrated through non-clinical testing to be magnetic resonance (MR) conditional. A patient with this device can be scanned safely in an MR system under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Static magnetic field of 1.5-Tesla or 3-Tesla, only;
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40-T/m);
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode.

Under the scan conditions defined, the implant system is expected to produce a maximum temperature rise of 3.6°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant system extends approximately 30-mm from this system when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

Parameter	Condition
Nominal Values of Static Magnetic Field (T)	1.5-Tesla and 3.0-Tesla
Maximum Spatial Field Gradient (T/m and gauss/cm)	40-T/m (4,000-gauss/cm)
Type of RF Excitation	Circularly Polarized (CP) (i.e., Quadrature-Transmission)
Transmit RF Coil Information	There are no transmit RF coil restrictions. Accordingly, the following may be used: body transmit RF coil and all other RF coil combinations (i.e., body RF coil combined with any receive-only RF coil, transmit/receive head RF coil, transmit/receive knee RF coil, etc.)
Operating Mode of MR System	Normal Operating Mode
Maximum Whole Body Averaged SAR	2-W/kg (Normal Operating Mode)
Limits on Scan Duration	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks).
MR Image Artifact	The presence of this implant system produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant system is located in the area of interest.

To allow medical professionals to identify the specific medical devices a patient has, the MRI safety status of the medical devices, and the conditions for safe use in the MR environment for MR Conditional devices, BioHorizons recommends that clinicians provide the patient with the dental implant(s) and dental abutment(s) / bar(s) device specific peel-offs affixed to the patient card. Patient cards are available free of charge upon request from BioHorizons or for direct printing at <https://ifu.biohorizons.com>.

### COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) abutment screw and/ or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/ or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

### HANDLING AND STERILIZATION

Always handle the device with powder-free gloves and avoid contact with hard objects that may damage the surface. If the device is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove product from the package only after the correct size has been determined and the surgical site has been prepared. For non-sterile devices, remove and discard any shipping material before initial processing. Non-sterile single use devices and devices intended for reuse must be cleaned and sterilized prior to use. BioHorizons devices have not been validated for automated cleaning. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations. Refer to the legal manufacturer's instructions for use of detergent solution preparation.
- 2) Brush the device to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse device under running utility tap water.
- 4) Place device in the sterile container filled with the prepared detergent solution and sonicate for two (2) minutes.
- 5) Thoroughly rinse device under running utility tap water.
- 6) Spray device with 70% isopropyl alcohol (IPA).
- 7) Blot device dry with clean lint free cloth.

For sterilization of non-sterile devices (not including ceramic prosthetics), place devices in an FDA cleared sterilization bag or wrap and run through one of the following qualified sterilization cycles:

Sterilization Cycles				
<b>Reference:</b>	AAMI TIR12:2020	AAMI TIR12:2020	AAMI TIR12:2020	UK HTM 01-01 Part C:2016
<b>Type:</b>	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	30 minutes	20 minutes	20 minutes

For non-sterile ceramic prosthetics, dry heat sterilize at 170°C (338°F) for sixty (60) minutes.


It is recommended to include a 30-minute cool-down period before removing the devices from the sterilization bag or wrap.

Devices to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

### SYMBOLS AND DESCRIPTIONS

The symbols table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic instructions for use.
	Manufacturer.
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number.
	Lot/ batch number.
	Unique Device Identifier.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.
	Date of manufacture.
<b>Rx Only</b>	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician.
	European Union Authorized Representative.
	Do not use if package is damaged. Discard device and package.
	Medical Device.

Symbol	Symbol Description
	Non-sterile.
	Single sterile barrier system with protective packaging outside.
	Single sterile barrier system.
	Home.
	Magnetic resonance warning: Device is MR conditional.
	United Kingdom Responsible Person.

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc. Z-CAD® is a registered trademark of Metoxit AG. CONELOG® is registered trademark of Camlog Biotechnologies GmbH.

**UK  
CA  
0086**

**BioHorizons UK Limited**  
 Reflex  
 Cain Road  
 Bracknell  
 Berkshire  
 RG12 1HL

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©Zimmer، وDental ScrewVent®، وTapered ScrewVent® عبارة عن علامات تجارية مسجلة لشركة Z-CAD®. Zimmer Dental, Inc. هي علامة تجارية مسجلة لشركة Camlog Biotechnologies GmbH. Metoxit AG. هي علامة تجارية مسجلة لشركة CONELOG®.