English	Instructions for Use: Dental Prosthetics
Español	Instrucciones de uso: Prótesis Dentales
русский	Инструкция по применению: Зубные протезы
Deutsch	Gebrauchsanweisung: Dentale Prothetik
Français	Instructions d'utilisation: Prothèses Dentaires
Italiano	Istruzioni per l'uso: Protesi Dentali
Português	Instruções de utilização: Próteses Dentárias
Türkçe	Kullanım talimatları: Dental Prostetik Ürünler
简体中文	使用说明：牙科义齿
日本語	使用説明書：歯科プロテーゼ
한국어	적응증: 치과 보철
عربي	إرشادات الاستخدام: البدليات السنية
Polski	Instrukcja obsługi: Protetyka Dentystyczna
Українська	Інструкція із застосування: зубні протези
Čeština	Návod k použití: Dentální protetické komponenty
Български	Инструкции за употреба: Дентални протези
Ελληνικά	Οδηγίες χρήσης: Οδοντικά προσθετικά εξαρτήματα
Magyar	Használati utasítás: Fogászati protézisek
Română	Instrucțiuni de utilizare: Protetică dentară
Nederlands	Gebruiksaanwijzing: tandheelkundige prothesen
Svenska	Bruksanvisning: Tandprotes
Norsk	Bruksanvisning: Tannproteser

BioHorizons Implant Systems Inc.
 2300 Riverchase Center
 Birmingham, AL 35244 USA
 TOLL-FREE 888.246.8338
 TEL 205.967.7880
 FAX 205.870.0304
www.biohorizons.com

This document supersedes all prior revisions. Original language is English.

This document applies to Intra-Lock (Fusion) Dental Abutments and Prosthetic Components. Every product packaging label contains a description of the enclosed device.

Devices in Scope
Intra-Lock Fusion healing caps and abutments
Intra-Lock Fusion PEEK temporary abutments provided sterile
Intra-Lock Fusion titanium temporary abutments
Intra-Lock Fusion gold abutments
Intra-Lock Fusion plastic abutments
Intra-Lock Fusion Multi-unit abutments
Intra-Lock Fusion titanium abutments
Intra-Lock Fusion Ti base abutments
Intra-Lock Fusion Custom milled abutments
Intra-Lock Fusion detachable abutments
Intra-Lock Fusion abutment screws
Intra-Lock Fusion prosthetic screws
Intra-Lock Fusion impression copings (conventional and digital)
Intra-Lock Fusion analogs and replicas
Intra-Lock Fusion abutment copings
Intra-Lock Fusion abutment attachments

DESCRIPTION

Intra-Lock prosthetic components are intended for the restoration of dental implants within the specific indications of each implant system. The label on each prosthetic package contains important product information including whether the prosthetic is supplied sterile or non-sterile. Refer to the following table for material(s) of implantable abutments.

Implantable Devices	Material (main elements)
Intra-Lock Fusion healing caps and abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Intra-Lock Fusion PEEK temporary abutments provided sterile	TecaPeek Classix (Polyetheretherketone)
Intra-Lock Fusion titanium temporary abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Intra-Lock Fusion gold abutments	Gold 6019 (Gold, Palladium, Platinum, Iridium)
Intra-Lock Fusion plastic abutments	Casting material as determined by clinician
Intra-Lock Fusion Multi-unit abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Intra-Lock Fusion titanium abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Intra-Lock Fusion Ti base abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Intra-Lock Fusion Custom milled abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Intra-Lock Fusion detachable abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Intra-Lock Fusion abutment screws	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)

Prosthetic connection color coding is defined in the following table:

Internal Hexed Prosthetic Color Coding	Grey	Yellow
Internal Hexed Prosthetic Size	3.0mm	3.5mm

INTENDED USE

Intra-Lock abutments are intended for the restoration of Intra-Lock dental implants in the mandible or maxilla, for single tooth replacement or for fixed bridgework and dental retention, within the specific indications of each implant system.

Intra-Lock abutment accessories are intended for use in conjunction with Intra-Lock dental abutments connected to an implanted artificial root structure in the mandible or maxilla for single tooth replacement or for fixed bridgework and dental retention, within the specific indications of each implant system.

INDICATIONS FOR USE

Intra-Lock Prosthetics are abutments that include healing abutments for contouring tissue and final restorative abutments to support a prosthesis. The abutments may be used for a single or multiple unit restoration and are compatible for use with Intra-Lock implants.

Intra-Lock Titanium Base Abutments and Laser-Lok Titanium Base Abutments are intended to be used as straight abutments.

Intra-Lock CAD/CAM Abutments are dental abutments placed onto a dental implant to provide support for dental prosthetic restorations. The abutments include: 1) Titanium abutment blanks with a pre-machined implant connection where the upper portion may be custom-milled in accordance with a patient-specific design using CAD/CAM techniques; and 2) Titanium bases with a pre-machined implant connection upon which a CAD/CAM designed superstructure may be fitted to complete a two-piece dental abutment. The abutments include an abutment screw for fixation to the underlying implant. The abutments may be used for single-unit (single-tooth) or multiple-unit (bridges and bars) restorations and are compatible for use with Intra-lock implants. All digitally designed abutments and/or copings for use with Intra-Lock CAD/CAM Abutments are intended to be sent to a BioHorizons-validated milling center for manufacture. Intra-Lock abutments designed using CAD/CAM techniques must fulfill the BioHorizons allowable range of design parameters.

CONTRAINDICATIONS

Intra-Lock prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity. Intra-Lock Multi-unit prosthetics are contraindicated for single tooth replacements.

PATIENT POPULATION

BioHorizons dental implant systems are intended for use in skeletally mature, non-pediatric edentulous or partially edentulous patients as long as the defined contraindications are nonapplicable.

INTENDED USERS

BioHorizons implant systems are only intended for use by licensed health care professionals, more specifically, BioHorizons implant systems are intended to be used by trained dentists, surgeons and dental technicians in a standard dental surgical setting, which may range from General Dentists' practice offices to Maxillofacial Surgical Operating Rooms, as well as laboratories for dental processes. Use of these products requires specialized knowledge and experience in implant dentistry. BioHorizons implant systems are marked and labeled as a medical device (MD) and are Rx only.

DIRECTIONS FOR USE

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions For

Use that accompany Intra-Lock products.

Healing abutments should be hand tightened (10-15 Ncm).

Abutment screws must be torqued to 30 Ncm to prevent screw loosening. Intra-Lock Multi-unit Abutments, Custom Milled Abutments, and Titanium Base Abutments contain a LIGHT BLUE abutment screw specifically designed for use in those abutments. Use of the LIGHT BLUE abutment screws with incompatible abutments may result in abutment and/orabutment screw failure.

Intra-Lock Multi-Unit Prosthetic Screws must be torqued to 15 Ncm to prevent screw loosening.

Intra-Lock Multi-Unit Prosthetic Component design parameters are as follows:

- Minimum post height should be equal to or greater than 4mm.
- Reduction of wall thickness is not permitted.
- Angle correction is not permitted. Angled Multi-unit Abutments shall be used for angle correction.

3.5mm Internal prosthetic connection compatibility with Zimmer® Dental Screw-Vent® and Tapered Screw-Vent® implants is provided in the following table.

Platform	Zimmer Tapered ScrewVent Implants*	Zimmer ScrewVent Implants*
3.5mm Internal Connection	TSVBx TSVHx TSV4Bx TSV4Hx TSVT4Bx	SVMBx SVHx SVBx SVHx SVMHx

*Where variable x = implant length.

Secure Custom Cast Abutments to the supplied polishing protector to protect the implant/abutment interface during laboratory procedures.

For CAD/CAM abutments, a custom abutment design file is created by the customer using a 3Shape abutment library. The abutment design file is converted to a stereolithography (.stl) file using 3Shape software. The .stl file is converted to a numerical control (.nc) file using an appropriate software application. After the custom abutment design is uploaded, the milling machine software executes the necessary commands to mill the final abutment. The digitally designed file is to be sent to a registered and listed BioHorizons contract manufacturer for milling. Intra-Lock CAD/CAM Abutments are compatible with commercially available dental CAD/CAM systems, such as 3Shape.

Design limitation parameters for the Custom Milled Abutments are included in the following table.

Feature	Custom Milled Platforms	
	Grey, Yellow	
	Min	Max
Wall thickness	0.3mm	3.7mm
Gingival Margin Diameter	3.3mm, 3.8mm respectively	9.5mm
Gingival Margin Height	0.3mm, 0.2mm respectively	12.3mm, 12.5mm respectively
Length of the Abutment Post	*	12.3mm, 12.5mm respectively
Taper of the Abutment	0°	90°
Angulation	0°	30°

* Minimum allowable abutment post dimensions are 3.0mm for screw retained abutment and 4.0mm for cement retained abutment.

Titanium base abutments – BioHorizons recommends the use of Metoxit Z-CAD® dental ceramic (or equivalent) to fabricate the custom zirconia superstructure and RelyX Unicem (or equivalent) to affix the zirconia superstructure to the titanium base. Design limitation parameters for the custom zirconia superstructure portion of the Titanium Base Abutments are included in the following table.

Feature	Titanium Base Platforms	
	Grey, Yellow	
	Min	Max
Wall Thickness	0.4mm	**
Gingival Margin Diameter	3.8mm, 3.8mm respectively	**
Length of the Abutment Post	4.0mm	**
Taper	Variable based on anatomy	**
Angulation of the Final Titanium Base	0°	20°

**Feature maximum sizes are based on patient anatomy and clinician discretion.

WARNINGS AND PRECAUTIONS

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any Intra-Lock Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of Intra-Lock prosthetic components.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include, but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics.

If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the modifier of the implant/abutment interface a medical device company subject to FDA rules and regulations.

Use of the LIGHT BLUE abutment screws (packaged with Angled Multi-unit Abutments, Custom Milled Abutments, and Titanium Bases) with incompatible abutments may result in abutment and/or abutment screw failure.

Small diameter implants and angled abutments are not recommended for use in the posterior region of the mouth.

Intra-Lock Multi-unit Abutments require the addition of a coping or cylinder for single and multiple-unit restorations. The coping or cylinder gingival collar, angulation, or wall thickness shall not be modified. The coping or cylinder post height shall not be less than 4mm for single-unit restorations.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

Devices should not be used in patients with a known allergy or sensitivity to the device material(s).

For implantable devices, the summary of safety and clinical performance (SSCP) according to Article 32 of Regulation (EU) 2017/745 can be found in the European database for medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed>.

Implantable Device	Basic UDI-DI Number
Intra-Lock Fusion healing caps and abutments	08472360IIBWETABUT001HS
Intra-Lock Fusion PEEK temporary abutments provided sterile	08472360IIBWETABUT004HY
Intra-Lock Fusion titanium temporary abutments	08472360IIBWETABUT003HW
Intra-Lock Fusion gold abutments	08472360IIBWETABUT005J2
Intra-Lock Fusion plastic abutments	08472360IIBWETABUT006J4
Intra-Lock Fusion Multi-unit abutments	08472360IIBWETABUT009JA
Intra-Lock Fusion titanium abutments	08472360IIBWETABUT007J6

Implantable Device	Basic UDI-DI Number
Intra-Lock Fusion Ti base abutments	08472360IIBWETABUT007J6
Intra-Lock Fusion Custom milled abutments	08472360IIBWETABUT007J6
Intra-Lock Fusion detachable abutments (e.g., ODSecure)	08472360IIBWETABUT008J8
Intra-Lock Fusion abutment screws	08472360IIBWETABUT010HT

MAGNETIC RESONANCE (MR) SAFETY INFORMATION

MR Conditional

Intra-Lock implant systems have been demonstrated through non-clinical testing to be magnetic resonance (MR) conditional. A patient with this device can be scanned safely in an MR system under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Static magnetic field of 1.5-Tesla or 3-Tesla, only;
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40-T/m);
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode.

Under the scan conditions defined, the implant system is expected to produce a maximum temperature rise of 3.6°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant system extends approximately 30-mm from this system when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

Parameter	Condition
Nominal Values of Static Magnetic Field (T)	1.5-Tesla and 3.0-Tesla
Maximum Spatial Field Gradient (T/m and gauss/cm)	40-T/m (4,000-gauss/cm)
Type of RF Excitation	Circularly Polarized (CP) (i.e., Quadrature-Transmission)
Transmit RF Coil Information	There are no transmit RF coil restrictions. Accordingly, the following may be used: body transmit RF coil and all other RF coil combinations (i.e., body RF coil combined with any receive-only RF coil, transmit/receive head RF coil, transmit/receive knee RF coil, etc.)
Operating Mode of MR System	Normal Operating Mode
Maximum Whole Body Averaged SAR	2-W/kg (Normal Operating Mode)
Limits on Scan Duration	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks).
MR Image Artifact	The presence of this implant system produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant system is located in the area of interest.

To allow medical professionals to identify the specific medical devices a patient has, the MRI safety status of the medical devices, and the conditions for safe use in the MR environment for MR Conditional devices, BioHorizons recommends that clinicians provide the patient with the dental implant(s) and dental abutment(s) / bar(s) device specific peel-offs affixed to the patient card. Patient cards are available free of charge upon request from BioHorizons or for direct printing at <http://eifu.intra-lock.com>.

COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2)) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) abutment screw and/ or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses

possibly involving macrophages and/ or fibroblasts; formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

HANDLING AND STERILIZATION

Always handle the device with powder-free gloves and avoid contact with hard objects that may damage the surface. If the device is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove device from the package only after the correct size has been determined and the surgical site has been prepared. Non-sterile device must be cleaned and sterilized (including CAD/CAM abutments after abutment preparation and before final placement in the restorative site). BioHorizons devices have not been validated for automated cleaning. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations. Refer to the manufacturer's instructions for use for detergent solution preparation.
- 2) Brush the device to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse device under running utility tap water.
- 4) Place product in the sterile tube filled with the prepared detergent solution and sonicate for two (2) minutes.
- 5) Thoroughly rinse device under running utility tap water.
- 6) Spray device with 70% isopropyl alcohol (IPA).
- 7) Blot device dry with clean lint free cloth.

For sterilization of non-sterile devices, place device in an FDA cleared sterilization bag or wrap, and run through one of the following qualified sterilization cycles:

Sterilization Cycles				
Reference:	AAMI TIR12:2020	AAMI TIR12:2020	AAMI TIR12:2020	UK HTM 01-01 Part C:2016
Type:	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
Exposure Time and Temperature:	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
Minimum Dry Time:	30 minutes	30 minutes	20 minutes	20 minutes

It is recommended to include a 30-minute cool-down period before removing the devices from the sterilization bag or wrap.

Devices to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

SYMBOLS AND DESCRIPTIONS

The symbols table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic Instructions for Use.
	Manufacturer.

Symbol	Symbol Description
	Intra-Lock products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/article number.
	Lot/batch number.
	Unique Device Identifier.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.
	Date of manufacture.
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution, and use by, or on the order of, a dentist or physician.
	European Union Authorized Representative.
	Do not use if package is damaged. Discard device and package.
	Medical Device.
	Non-sterile.
	Single sterile barrier system with protective packaging outside.
	Single sterile barrier system.
	Home position.
 MR Conditional	Magnetic resonance warning: Device is MR conditional.

Symbol	Symbol Description		
<table border="1"><tr><td data-bbox="274 195 358 254">UK</td><td data-bbox="358 195 443 254">RP</td></tr></table>	UK	RP	United Kingdom Responsible Person.
UK	RP		

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc. Z-CAD® is a registered trademark of Metoxit AG.

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