ENGLISH	Instructions for use: AutoTac
ESPAÑOL	Instrucciones de uso: AutoTac
РУССКИЙ	Инструкции по эксплуатации: AutoTac
DEUTSCH	Gebrauchsanweisung: AutoTac
FRANÇAIS	Instructions d'utilisation : AutoTac
ITALIANO	Istruzioni per l'uso: AutoTac
PORTUGUÊS	Instruções de utilização: AutoTac
TÜRKÇE	Kullanım talimatları: AutoTac
العربية	AutoTac: تعليمات استخدام
POLSKI	Instrukcja obsługi: AutoTac
УКРАЇНСЬКА	Інструкції із застосування: АутоТак

BioHorizons Implant Systems Inc
2300 Riverchase Center
Birmingham AL, 35244 USA
TOLL FREE 888.246.8338
TEL 205.967.7880 FAX 205.870.0304
www.biohorizons.com

This document supersedes all prior revisions. Original language is English. Instructions for Use are available free of charge in printed form upon request from BioHorizons or your local representative / distributor.

DESCRIPTION

The AutoTac System kit is a system containing fixation tacks and a tack delivery handle.

This product is NOT supplied sterile. After removing packaging and supporting material (e.g., cable ties), the AutoTac kit must be DISASSEMBLED, cleaned, reassembled, and sterilized before initial use and each subsequent re-use.

INDICATIONS

The AutoTac System is intended for use in the fixation of barrier membranes and/or soft tissue used for regeneration of bone and/or soft tissue in the oral cavity. It may be used in conjunction with commercially available membranes.

CONTRAINDICATIONS

The AutoTac System is not intended for use in the fixation of barrier membranes and/or soft tissue other than within the oral cavity. Known contraindications associated with Guided Bone Regeneration/Guided Tissue Regeneration (GBR/GTR) procedures and/or dental implant treatment include, but are not limited to: infections, insufficient available bone or bone of poor quality, vascular impairment at surgical site, poor oral hygiene, uncontrolled diabetes, heavy smoking or tobacco abuse, drug or alcohol abuse, chronic high dose steroid therapy, patients with a heart valve or other prosthetic device, heart valve defect (e.g., heart murmur, prolapsed mitral valve, history of rheumatic heart disease, etc.), medical conditions such as blood/clotting disorders, current or ongoing anticoagulant therapy, metabolic bone disease or other metabolic or systemic disorders which may adversely affect bone or wound healing or cases in which the available bone is too diminished to provide adequate width or height to adequately hold guided bone regeneration appliances or implants.

DIRECTIONS FOR USE

Place the user-sterilized delivery handle and tack cassette containing the tacks, on a sterile tray. Load the delivery handle by pressing the tip vertically into the loading recess of the tack cassette until the trigger button pops up and you hear an audible click. If you are having difficulty loading the delivery handle, try angling the tip at a steeper angle. Once the delivery handle is loaded, insert the tip into one of the recesses in the tack cassette that contains a tack and pick up a tack. You will be able to see the pointed end of the tack protruding about 1mm from the nozzle tip. The protruding shank should be coaxial to the nozzle tip for proper loading. If it is not, align the shaft axis with the tip of a sterile instrument.

Place the membrane on the defect/grafted site following the instructions for use specific to the membrane. Place the tip of the delivery handle in the desired location, press against the bone, and press the trigger button to deliver the tack. If the tack head is not flush, cock the delivery handle again. Carefully place the nozzle tip over the tack head and press the trigger to seat the tack further. Stabilize the membrane with as many tacks as required.

WARNINGS AND PRECAUTIONS

It is strongly recommended that the instructions for use provided by the manufacturers of membranes, fixation devices and, where applicable, dental implant and restorative devices be read, understood and followed prior to undertaking any treatment involving the use of such devices. Surgical techniques required to properly utilize these devices and procedures are highly specialized and complex. Therefore, it is recommended that specialized training be undertaken prior to performing any such procedure. Improper patient selection and technique can cause or contribute to case failure with possible loss of supporting bone. Improper technique can result in tack failure and/or loss of membrane stabilization. Clinicians must carefully consider the therapeutic benefit of such treatment before proceeding. Membranes should not be used and affixed with the tack in situations where there are: defects with severe horizontal bone loss and with little remaining periodontal ligament and cementum; defects that do not allow for creation and maintenance of a space; and, defects that exhibit characteristics of inadequate regeneration.

Titanium tacks can be damaged in function for a number of reasons including improper loading or improper surgical placement. An adequate number of tacks should be used to stabilize the membrane/ tissue. A certain percentage of tacks may fail to fixate rigidly to adjacent bone and provide proper retention to the membrane. Tacks demonstrating mobility should be removed.

The AutoTac System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the AutoTac System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

COMPLICATIONS AND ADVERSE EFFECTS

Possible complications with any periodontal surgery include thermal sensitivity, gingival recession, flap sloughing, resorption or ankylosis of a treated root, some loss of crestal bone height, perforation or abscess formation, infection, pain, gingival irregularities, and complications associated with the use of anesthesia. Depending on the type and severity of the complication, as judged by the clinician, membranes and support or fixation devices may require removal and/or antibiotic therapy. Visual inspection as well as panoramic and periapical radiographs are essential to determine anatomical landmarks, occlusal conditions, periodontal status, and adequacy of bone.

Preoperative treatment may include antibiotic therapy at clinician's discretion. As with any surgical procedure, careful post-operative management is important for optimal healing, including oral hygiene maintenance, plaque control, adequate flossing and brushing instructions and close patient monitoring and professional prophylaxis at least every other week for the first eight weeks. Post-surgical exposure of secured (or unsecured) membranes may occur and should be anticipated. The situation should be monitored but should not interfere with the GBR/GTR process.

BioHorizons Titanium tacks are not intended for permanent implantation and should be removed after proper healing. Tack removal is performed by exposing the surgical site and using a scalpel blade, periosteal elevator, or other similar thin flat surface, prying the head of the tack away from the underlying bone. Removed tacks should be accounted for and discarded, and the surgical site then closed and re-sutured. In the event of tissue inflammation or evidence of infection, and at the clinician's discretion, both the membrane and Titanium Tacks may be removed. Some membranes are designed to resorb and will lose structural integrity over time. Removal of the material may require curettage and should be accompanied by thorough flushing with a sterile saline solution.

SURGICAL COMPLICATIONS

Potential risk associated with the AutoTac System includes respiratory, and/or intestinal track damage as a result of component aspiration or ingestion. Periodontal surgery that involves GBR/GTR procedures and implants has risks, including localized swelling, dehiscence, tenderness of short duration, edema, hematoma, or bleeding. Numbness of the lower lip and chin region following lower jaw surgery, and of the tissue beside the nose following upper jaw surgery, is a possible side effect of the surgery; though most likely of a temporary nature, in very rare cases such numbness has been permanent. Gingival/Mucosal (gum disease) ulceration, tissue reaction, or infection may occur, but generally responds to local care.

HANDLING AND STERILIZATION

This product is NOT supplied sterile. The AutoTac kit must be cleaned and sterilized before initial use and each subsequent re-use. For non-sterile product, remove and discard any shipping material before initial processing. Non-sterile single use devices and devices intended for reuse must be cleaned and sterilized prior to use. BioHorizons devices have not been validated for automated cleaning. To clean:

- 1) DISASSEMBLE. See below for disassembly/assembly instructions.
- 2) Brush external and internal surfaces of the handle thoroughly using a broad spectrum cleaning agent such as Hu-Friedy's (800-483-7433) Enzymax® or equivalent.
- 3) Rinse handle with water.
- 4) Blow pressurized air through the nozzle to purge.
- 5) REASSEMBLE. See below for disassembly/assembly instructions.
- 6) Use a medical instrument lubricant such as STERIS' Hinge Free® after each use. Refer to the legal manufacturer's instructions for use for lubricant preparation.

Titanium tacks must be placed in the tack holder for sterilization; always handle tacks with powder-free gloves and avoid contact with hard objects that may damage the surface. The AutoTac kit can be sterilized using one of four qualified steam sterilization cycles:


Sterilization Cycles				
Type:	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
Exposure Time and Temperature:	50 minutes at 121°C (250°F)	30 minutes at 132°C (270°F)	5 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
Dry Time Range:	20 to 50 minutes*	20 to 50 minutes*	20 to 50 minutes*	20 minutes
*Recommended drying time range				

DO NOT “flash” (278°F+) sterilize, which may damage the tack cassette. Tacks are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

Product to be disposed of must be treated and decontaminated as dental surgery waste in compliance with relevant local regulations.

AUTOTAC ASSEMBLY COMPONENTS

Cleaning the AutoTac Delivery Handle requires disassembly and reassembly using a 5/64” Hex Key. Replacement 5/64” Hex Keys and Handle Screws/Button Spring may be ordered from BioHorizons, part numbers 400-000 and 400-001, respectively.

	Part Name
1	Handle
2	Button
3	Button Spring
4	Handle Screw
5	Firing Spring
6	Drive Shaft
7	Nozzle
8	Spacer


AutoTac System Kit Re-order Items:

- 400-200.....AutoTac Delivery Handle
- 400-248.....Autoclavable Tack Cassette (tacks not included)
- 400-260..... Titanium Tacks with Cassette (pack of 21)
- 400-240..... Titanium Tack Vial (pack of 5)
- HF-DPU17..... Dressing Pliers, Utility Pick-up
- 400-000.....AutoTac System IFU and Hex Key
- 400-001.....AutoTac Spare Parts Kit: (Hex Key, Handle Screws, Button Spring)

DISASSEMBLY INSTRUCTIONS

Attention: It is strongly recommended that gloves and eye protection be worn during disassembly and assembly of the AutoTac

Delivery Handle. The assistance of a second person may be required during disassembly and assembly.


1. To release the nozzle from the handle, hold on to the flange of the nozzle as shown, maintain pressure on the button and slowly loosen the two handle screws approximately two full turns using a 5/64" hex key.


Caution: Removing the screws without holding on to the flange of the nozzle can allow the nozzle to release from the handle under pressure. Never point the nozzle toward your eyes or toward others.


2. Keeping pressure on the button, slide the nozzle completely out of the handle.

3. Release pressure from the button to remove the button and the button spring.

4. Remove the handle screws completely from the handle and set aside.


5. Remove the drive shaft, firing spring and spacer from the nozzle using gentle tension.

5. Brush external and internal surfaces of the handle thoroughly using a broad spectrum cleaning agent such as Hu-Friedy's (800-483-


7433) Enzymax® or equivalent; rinse with water and blow pressurized air through the nozzle to purge. Reassemble.

ASSEMBLY INSTRUCTIONS


Attention: It is strongly recommended that gloves and eye protection be worn during assembly of the AutoTac Delivery Handle. The assistance of a second person may be required during assembly.


1. To begin, insert the drive shaft into the nozzle, taking care to align the curvature of the drive shaft with the curvature of the nozzle. Set the nozzle assembly aside for the moment.


2. Pre-assemble the handle screws into the handle 3-4 full rotations.


3. Holding the handle vertically, drop the spacer and the firing spring into the handle.

4. Carefully place the button spring into the small spring hole in the handle. Place the button on top of the button spring with the cut-out notch facing the firing end of the handle. Press the button down until fully seated and hold in place for final assembly.


5. Insert the nozzle assembly through the hole in the button and into the handle making sure the screw holes in the handle and nozzle are in alignment and the nozzle tip is angled away from the firing button (**See step 7 for final assembled**


6. Hold the nozzle firmly against the handle with one hand and use your free hand to tighten the handle screws. Hand tighten using a 5/64" hex key.

Tip: One person can use two hands to hold the nozzle flush with the handle while the other person tightens the handle screws.

7. A correctly assembled handle should have the nozzle bending away from the button (top) of the handle, as shown below. Verify correct assembly by completing a single loading/firing cycle with the handle.


8. Use a medical instrument lubricant such as STERIS' (800-548-4873) Hinge Free® after each use. The AutoTac Kit can be sterilized using one of four qualified steam sterilization cycles found on the reverse side in the section titled Handling and Sterilization.

SYMBOLS

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic instructions for use

Symbol	Symbol Description
	Manufacturer
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number
	Lot/ batch number
	Unique Device Identifier
	Do not re-use
	Do not re-sterilize
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician
	European Union Authorized Representative
	Do not use if package is damaged. Discard device and package.
	Medical Device
	Non-sterile
	Single sterile barrier system with protective packaging outside
	Single sterile barrier system

Symbol	Symbol Description
	Home
	Magnetic resonance warning: Device is MR conditional
	United Kingdom Responsible Person

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