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L02035CAen Rev B JUL 2024

ENGLISH

This document supersedes all prior revisions. Original language is English.

This document applies to BioHorizons Dental Abutments and Prosthetic components for sale in Canada. Every product packaging label contains a description of the enclosed device.

DESCRIPTION

BioHorizons prosthetic components are intended for the restoration of BioHorizons dental implants within the specific indications of each implant system.

BioHorizons Multi-unit Abutments for CONELOG® and CONELOG ODSecure Abutments are intended for the restoration of CONELOG dental implants within the specific indications of each implant system.

The label on each prosthetic package contains important product information including whether the prosthetic is supplied sterile or non-sterile.

INDICATIONS FOR USE

Tapered Internal and Internal Prosthetics are intended for the restoration of Tapered Internal and Internal dental implants and Zimmer® Dental ScrewVent® and Tapered ScrewVent® dental implants. Tissue Level Prosthetics are intended for the restoration of BioHorizons Tissue Level dental implants. Ceramic abutments are only intended for cement-retained cases in the anterior esthetic zone. BioHorizons Angled Multi-unit Abutments contain a unique LIGHT BLUE abutment screw specifically designed for use in those abutments. Multi-unit and Abutment for Screw prosthetics are intended for use as a support for multiple-implant supported prostheses. Ball Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible (3.0mm Ball Abutments being intended for use with 3.0mm dental implants). ODSecure Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible and maxilla. Easy Ti abutments are intended for long term temporary restorations (>30 days); hexed for single-unit and non-hexed for multiple-unit restorations. PEEK prosthetics are intended for short term temporary restorations (<30 days).

The BioHorizons Multi-unit Abutments for CONELOG are intended to function in the mandible or maxilla to support single and multiple-unit temporary or definitive restorations on CONELOG dental implants. CONELOG ODSecure Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible and maxilla.

CONTRAINDICATIONS

BioHorizons prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity. BioHorizons Multi-unit and Abutment for Screw prosthetics are contraindicated for single tooth replacements.

DIRECTIONS FOR USE

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions For Use that accompany BioHorizons products. Healing abutments should be hand tightened (10-15 Ncm). Abutment screws must be torqued to 30 Ncm to prevent screw loosening.

Platform	Zimmer Tapered ScrewVent Implants*	Zimmer ScrewVent Implants*
3.5mm Internal Connection	TSVBx TSV4Bx TSVT4Bx	TSVHx TSV4Hx SVMBx SVBx SVMHx SVHx
4.5mm Internal Connection	TSVWBx TSVTWBx TSVWHx	SVWBx SVWHx
5.7mm Internal Connection	TSV6Bx TSVT6Bx TSV6Hx	N/A

*Where variable x = implant length.

Secure Custom Cast Abutments to the supplied polishing protector to protect the implant/abutment interface during laboratory procedures.

Secure Laser-Lok abutments onto the Laser-Lok Protective Sleeve to protect Laser-Lok during handling.

BioHorizons Multi-unit Abutment screws must be torqued to 30 Ncm.

The BioHorizons Multi-unit straight and angled abutment screws for CONELOG are packaged with the corresponding abutments and must be torqued to 20Ncm for intended function.

The BioHorizons Multi-unit Prosthetic Component screws must be torqued to 15Ncm. BioHorizons Multi-unit Prosthetic Component design parameters are as follows:

- Minimum post height shall be equal to or greater than 4mm.
- Reduction of wall thickness is not permitted.
- Angle correction is not permitted. Angled Multi-unit abutments shall be used for angle correction.

BioHorizons Multi-unit Prosthetic Components used with the Multi-unit Abutments for CONELOG are intended to support multiple-unit temporary or definitive restorations except for the Titanium Coping. The Titanium Coping used with the Multi-unit Abutments for CONELOG is intended to support single and multiple-unit temporary or definitive restorations.

WARNINGS AND PRECAUTIONS

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of BioHorizons prosthetic components. Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at www.biohorizons.com. Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported immediately to BioHorizons Customer Care at 1-866-468-8338.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the modifier of the implant/abutment interface a medical device company subject to FDA rules and regulations. Leaving bonding agent or particulate in the Laser-Lok surface may affect the connective tissue attachment. Use of the LIGHT BLUE abutment screws (packaged with Angled Multi-unit Abutments) with incompatible abutments may result in abutment and/or abutment screw failure.

Small diameter implants and angled abutments are not recommended for use in the posterior region of the mouth.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

Devices should not be used in patients with a known allergy or sensitivity to the device material(s).

BioHorizons Multi-unit Abutments and Multi-unit Abutments for CONELOG require the addition of a coping or cylinder for single and multiple-unit restorations. The coping or cylinder gingival collar, angulation, and wall thickness shall not be modified. The coping or cylinder post height shall not be less than 4mm for single-unit restorations.

MAGNETIC RESONANCE (MR) SAFETY INFORMATION

MR-Conditional

Non-clinical testing has demonstrated BioHorizons Implant Systems are MR conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 4,000 G/cm (40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined above, the implant system is expected to produce a maximum temperature rise of 3.6 °C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant system extends approximately 30 mm from this implant when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) abutment screw and/ or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/ or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

HANDLING AND STERILIZATION

Always handle the device with powder-free gloves and avoid contact with hard objects that may damage the surface. If the device is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove device from the package only after the correct size has been determined and the surgical site has been prepared. Non-sterile devices must be cleaned and sterilized. BioHorizons devices have not been validated for automated cleaning. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a container using a broad-spectrum cleaning agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations. Refer to the legal manufacturer's instructions for use/preparation of detergent solution.
- 2) Brush the device to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.

- 3) Thoroughly rinse device under running utility tap water.
- 4) Place device in the sterile tube filled with the prepared detergent solution and sonicate for two (2) minutes.
- 5) Thoroughly rinse device under running utility tap water.
- 6) Spray device with 70% isopropyl alcohol (IPA).
- 7) Blot device dry with clean lint free cloth.

For sterilization of non-sterile devices (not including ceramic prosthetics), place devices in an FDA cleared sterilization bag or wrap and run through one of the following qualified sterilization cycles:

Sterilization Cycles				
Reference:	AAMI TIR12:2020	AAMI TIR12:2020	AAMI TIR12:2020	UK HTM 01-01 Part C:2016
Type:	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
Exposure Time and Temperature:	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
Minimum Dry Time:	30 minutes	30 minutes	20 minutes	20 minutes

For non-sterile ceramic prosthetics, dry heat sterilize at 170°C (338°F) for sixty (60) minutes.

It is recommended to include a 30-minute cool-down period before removing the devices from the sterilization bag or wrap.

Devices to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant regulations.

SYMBOLS AND DESCRIPTIONS

The symbols table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic instructions for use.
	Manufacturer.
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number.
	Lot/ batch number.

Symbol	Symbol Description
	Unique Device Identifier.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.
	Date of manufacture.
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician.
	European Union Authorized Representative.
	Do not use if package is damaged. Discard device and package.
	Medical Device.
	Non-sterile.
	Single sterile barrier system with protective packaging outside.
	Single sterile barrier system.
	Home Position.
 MR Conditional	Magnetic resonance warning: Device is MR conditional.
	United Kingdom Responsible Person.

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc. Z-CAD® is a registered trademark of Metoxit AG. Conelog® is registered trademark of Camlog Biotechnologies GmbH.

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