English	Instructions for use: Dental Implants
Español	Instrucciones de uso: Implantes Dentales
Русский	Инструкция по применению: Зубные имплантаты
Deutsch	Gebrauchsanweisung: Dentalimplantate
Français	Instructions d'utilisation: Implants Dentaires
Italiano	Istruzioni per l'uso: Impianti Dentali
Português	Instruções de utilização: Implantes Dentários
Türkçe	Kullanım talimatları: Dental İmplantlar
简体中文	使用说明: 牙种植体
日本語	使用説明書: 歯科インプラント
한국어	적응증: 치과 시술용 임플란트
عربي	إرشادات الاستخدام: غرسات الأسنان
Polski	Instrukcja obsługi: Implanty dentystyczne
Українська	Інструкції з використання: зубні імплантати.
Čeština	Návod k použití: Zubní protetické komponenty
Български	Инструкции за употреба: Дентални импланти
Ελληνικά	Οδηγίες χρήσης: Οδοντικά εμφυτεύματα
Magyar	Használati utasítás: Fogászati implantátumok
Română	Instrucțiuni de utilizare: implanturi dentare
Nederlands	Gebruiksaanwijzing: tandheelkundige implantaten
Svenska	Bruksanvisning: Tandimplantat
Norsk	Bruksanvisning: Dentalimplantater

BioHorizons Implant Systems Inc.

2300 Riverchase Center
Birmingham AL, 35244 USA
TOLLFREE 888.246.8338
TEL 205.967.7880
FAX 205.870.0304
www.biohorizons.com

This document supersedes all prior revisions. Original language is English.

This document applies to the following BioHorizons Dental Implants. Every product packaging label contains a description of the enclosed device:

Devices in Scope
BioHorizons Tapered Internal implants
BioHorizons Tapered Tissue Level implants
BioHorizons Tapered Plus implants
BioHorizons Tapered Pro implants
BioHorizons Tapered Short implants
BioHorizons Tapered PTG implants
BioHorizons Tapered IM implants
BioHorizons Internal Hex implants
BioHorizons External D3 implants

DESCRIPTION

BioHorizons implants are manufactured from biocompatible titanium alloy (Ti-6AL-4V ELI) and are single use only.

Element	Composition % (mass/mass)
Nitrogen, max	0.05
Carbon, max	0.08
Hydrogen, max	0.012
Iron, max	0.25
Oxygen, max	0.13
Aluminum	5.5 – 6.50
Vanadium	3.5 – 4.5
Titanium	balance

Some BioHorizons implants come packaged with a 3inOne abutment pre-mounted. 3inOne abutments are for single patient use only. This abutment may be used as: (1) an implant mount; (2) a closed tray transfer when used with a Ball-top Screw and (3) a temporary or final abutment. “Peel-and-stick” labels supplied on the implant package contain important product information and should be applied to the patient’s record in the event future reference is necessary.

INTENDED USE

BioHorizons implants are intended for use in the mandible or maxilla for use as an artificial root structure for single tooth replacement or for fixed bridgework and dental retention.

INDICATIONS FOR USE

BioHorizons implants are intended for use in the mandible or maxilla for use as an artificial root structure for single tooth replacement or for fixed bridgework and dental retention. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion or (2) when splinted together for multiple tooth replacement or when stabilized with an overdenture supported by multiple implants.

BioHorizons Tapered Short Implants are intended for use in the mandible or maxilla as an artificial root structure for single

tooth replacement or for fixed bridgework and dental retention. The implants may be restored using delayed loading, or with a terminal or intermediate abutment for fixed or removable bridgework, and for overdentures.

BioHorizons Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants may be used as an artificial root structure for single tooth replacement of mandibular central and lateral incisors and maxillary lateral incisors. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion, (2) when splinted together as an artificial root structure for multiple tooth replacement of mandibular incisors, or (3) for denture stabilization using multiple implants in the anterior mandible and maxilla. The implants may be placed in immediate function when good primary stability has been achieved and with appropriate occlusal loading.

CONTRAINDICATIONS

BioHorizons dental implants should not be used in patients who have contraindicating systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity. Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants are contraindicated for single tooth replacement of maxillary central incisors.

PATIENT POPULATION

BioHorizons dental implant systems are intended for use in skeletally mature, non-pediatric edentulous or partially edentulous patients as long as the defined contraindications are nonapplicable.

INTENDED USERS

BioHorizons implant systems are only intended for use by licensed health care professionals, more specifically, BioHorizons implant systems are intended to be used by trained dentists, surgeons and dental technicians in a standard dental surgical setting, which may range from General Dentists' practice offices to Maxillofacial Surgical Operating Rooms, as well as laboratories for dental processes. Use of these products requires specialized knowledge and experience in implant dentistry. BioHorizons implant systems are marked and labeled as a medical device (MD) and are Rx only.

DIRECTIONS FOR USE

Proper surgical procedures and techniques are the responsibility of the medical professional. Each dentist must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of postgraduate dental implant education and strict adherence to the instructions and procedures for use that accompany BioHorizons implant products. NOTE: A pre-operative 30-second rinse with a 0.12% Chlorhexidine Digluconate solution is recommended. (The influence of 0.12% Chlorhexidine Digluconate Rinses on the Incidence of Infectious Complications and Implant Success. Lambert, Paul M., et al, J Oral Maxillofacial Surgery 55:25-30, 1997, Suppl5). Tapered Internal Plus Implants - IMPORTANT! When using the TSK2021 (Tapered Internal Surgical Kit) to place the 3.8mm diameter implants, 3.0mm implant-level drivers are also required (TP3IDH, TP3IDR). Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants - IMPORTANT! (1) The implant procedure manual should be reviewed prior to attempting placement. (2) The osteotomy may be prepared with any BioHorizons 2.0mm and 2.5mm drill used in sequence. (3) Torque of 70Ncm or greater may compromise the mechanical integrity of the implant. In dense (D1-D2) bone further preparation of the site may be required to avoid applying excessive torque to the implant connection. (4) To further prepare the site, the crestal bone drill (TP3CBD) and the bone tap (122-900) may be used. (5) The clinician should verify that the proper instruments are available and at hand prior to attempting placement.

WARNINGS AND PRECAUTIONS

Appropriate training in proper implant surgery technique is strongly recommended prior to implant use. Improper technique can result in implant failure and/or loss of supporting bone. Appropriate x-ray films and/or CT scans should be utilized to determine (1) if adequate bone width and depth are available at the desired implant site and (2) the location of important anatomical landmarks, such as the mandibular canal, maxillary sinuses, and adjacent teeth. Clinician judgment should be used in determining the minimum post-implantation time before placing the implants in occlusal function. An adequate number of implants should be used to provide support and to distribute the load to the abutments.

For short implants, clinicians should closely monitor patients for any of the following conditions: periimplant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting the short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration to avoid immediate loading.

Devices should not be used in patients with a known allergy or sensitivity to the device material.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. A certain percentage of implants may fail to achieve or to maintain osseointegration, as demonstrated by mobility, and should be removed. Potential causes of abutment fracture include, but not limited to: casting beyond a 30-degree angle, casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. Small diameter implants with either angled or straight abutments are intended for the anterior region of the mouth and not intended for the posterior region of the mouth due to possible failure of the implant. Tapered IM Implants (1) are limited to the molar region; (2) restricted to only one implant per molar extraction site; (3) are to be placed 3-5mm from adjacent implants or teeth, with the distance from the apical extent of the implant to adjacent implant or implant to tooth contact point should be 3mm to coronal to the alveolar crest (**Figure 1**); (4) should be at least 1.5mm from buccal or lingual alveolar surfaces to permit adequate integration; and, (5) are used when the surgeon has determined that the placement of a narrower diameter implant would increase the probability of failure due to poor primary stability, or increased surgical procedures leading to complications.


Figure 1. Warnings and Precautions (3)

Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at www.biohorizons.com. Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

For implantable devices, the summary of safety and clinical performance (SSCP) according to Article 32 of Regulation (EU) 2017/745 can be found in the European database for medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed>.

Implantable Device	Basic UDI-DI Number
BioHorizons Tapered Internal implants	08472360IIBWETIMPL001P2
BioHorizons Tapered Tissue Level implants	08472360IIBWETIMPL002P4
BioHorizons Tapered Plus implants	08472360IIBWETIMPL003P6
BioHorizons Tapered Pro implants	08472360IIBWETIMPL004P8
BioHorizons Tapered Short implants	08472360IIBWETIMPL005PA
BioHorizons Tapered PTG implants	08472360IIBWETIMPL006PC
BioHorizons Tapered IM implants	Not applicable.
BioHorizons Internal Hex implants	Not applicable.
BioHorizons External D3 implants	Not applicable.

MAGNETIC RESONANCE (MR) SAFETY INFORMATION

MR Conditional

BioHorizons implant systems have been demonstrated through non-clinical testing to be magnetic resonance (MR) conditional. A patient with this device can be scanned safely in an MR system under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Static magnetic field of 1.5-Tesla or 3-Tesla, only;
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40-T/m);
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode.

Under the scan conditions defined, the implant system is expected to produce a maximum temperature rise of 3.6°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant system extends approximately 30-mm from this system when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

Parameter	Condition
Nominal Values of Static Magnetic Field (T)	1.5-Tesla and 3.0-Tesla
Maximum Spatial Field Gradient (T/m and gauss/cm)	40-T/m (4,000-guass/cm)
Type of RF Excitation	Circularly Polarized (CP) (i.e., Quadrature-Transmission)
Transmit RF Coil Information	There are no transmit RF coil restrictions. Accordingly, the following may be used: body transmit RF coil and all other RF coil combinations (i.e., body RF coil combined with any receive-only RF coil, transmit/receive head RF coil, transmit/receive knee RF coil, etc.)
Operating Mode of MR System	Normal Operating Mode
Maximum Whole Body Averaged SAR	2-W/kg (Normal Operating Mode)
Limits on Scan Duration	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks).
MR Image Artifact	The presence of this implant system produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant system is located in the area of interest.

To allow medical professionals to identify the specific medical devices a patient has, the MRI safety status of the medical

devices, and the conditions for safe use in the MR environment for MR Conditional devices, BioHorizons recommends that clinicians provide the patient with the dental implant(s) and dental abutment(s) / bar(s) device specific peel-offs affixed to the patient card. Patient cards are available free of charge upon request from BioHorizons or for direct printing at <https://ifu.biohorizons.com>.

COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

HANDLING AND STERILIZATION

This product is supplied sterile and should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove from the package only after the correct size has been determined and the surgical site has been prepared. Always handle the product with powder-free gloves, and avoid contact with hard objects that may damage the surface. Implants are single use only, and re-use should not be attempted. Following this guidance eliminates the risk of cross-patient contamination from secondary use of this device. BioHorizons assumes no responsibility for attempted re-use or re-sterilization.

Products to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

SYMBOLS AND DESCRIPTIONS

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic instructions for use
	Manufacturer
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number
	Lot/ batch number
	Unique Device Identifier

Symbol	Symbol Description
	Do not re-use
	Do not re-sterilize
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician
	European Union Authorized Representative
	Do not use if package is damaged. Discard device and package.
	Medical Device
	Non-sterile
	Single sterile barrier system with protective packaging outside
	Single sterile barrier system
	Home
 MR Conditional	Magnetic resonance warning: Device is MR conditional
	United Kingdom Responsible Person

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