<b>English</b>	Instructions for Use: MDL/MILO Dental Prosthetics
<b>Español</b>	Instrucciones de uso: Prótesis dentales MDL/MILO
<b>Français</b>	Mode d'emploi: Prothèses dentaires MDL/MILO

**BioHorizons Implant Systems Inc.**

2300 Riverchase Center  
Birmingham, AL 35244 USA  
TOLL FREE 888.246.8338  
TEL 205.967.7880  
FAX 205.870.0304  
[www.biohorizons.com](http://www.biohorizons.com)

L02070 Rev E FEB 2025

This document supersedes all prior revisions. Original language is English.

This document applied to Intra-Lock Mini Drive-Lock (MDL) and MILO Dental Prosthetics.

## **DESCRIPTION**

Intra-Lock MDL and MILO prosthetic components are intended for the restoration of MDL and MILO dental implants within the specific indications of each implant system. Intra-Lock MDL and MILO abutments are manufactured from Titanium Alloy (Ti-6Al-4V ELI). Temporary abutments are manufactured from PEEK. The label on each prosthetic package contains important product information. All prosthetic components are supplied non-sterile.

## **INTENDED USE**

Intra-Lock MDL and MILO prosthetics are intended for the restoration of removable or fixed prostheses supported by Intra-Lock MDL and MILO dental implants in partially or fully edentulous mandibles or maxillas, within the specific indications of each implant system.

## **INDICATIONS FOR USE**

The MILO and Mini Drive Lock Implant System Prosthetics have been designed to restore partially or fully edentulous patients. The abutments have been designed to be used in either the mandible or maxilla and to support removable or fixed prosthesis.

## **CONTRAINDICATIONS**

Patients with alcohol addiction or psychiatric disorders, blood dyscrasias, uncontrolled diabetes, hyperthyroidism, oral infections, malignancies, or patients who have had myocardial infarction within the last 12 months.

Patients with systemic diseases that compromise the immune system, such as AIDS, patients on medications that would compromise healing of an implant site, patients with a history of poor or noncompliance to oral hygiene procedures, or patients who cannot maintain oral hygiene procedures if implants are placed.

Tobacco usage increases the occurrence of complications and failures.

## **PATIENT POPULATION**

BioHorizons dental implant systems are intended for use in skeletally mature, non-pediatric edentulous or partially edentulous patients as long as the defined contraindications are nonapplicable.

## **INTENDED USERS**

BioHorizons implant systems are only intended for use by licensed health care professionals, more specifically, BioHorizons implant systems are intended to be used by trained dentists, surgeons and dental technicians in a standard dental surgical setting, which may range from General Dentists' practice offices to Maxillofacial Surgical Operating Rooms, as well as laboratories for dental processes. Use of these products requires specialized knowledge and experience in implant dentistry. BioHorizons implant systems are marked and labeled as a medical device (MD) and are Rx only.

## **DIRECTIONS FOR USE**

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions for Use (IFU) that accompany Intra-Lock products.

## **WARNINGS AND PRECAUTIONS**

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any Intra-Lock IFU. Clinicians are responsible for understanding the appropriate technical use of Intra-Lock prosthetic components. Additional technical information is available upon request from BioHorizons or may be viewed and/or downloaded at [www.intra-lock.com](http://www.intra-lock.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU.

Dental abutments can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. The use of angled abutments in the posterior region of the mouth is not recommended due to possible failure of the implant.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. Intra-Lock assumes no responsibility for attempted re-use or re-sterilization between patients.

Devices should not be used in patients with a known allergy or sensitivity to the device material(s).

## MAGNETIC RESONANCE (MR) SAFETY INFORMATION

### MR Conditional

BioHorizons implant systems have been demonstrated through non-clinical testing to be magnetic resonance (MR) conditional. A patient with this device can be scanned safely in an MR system under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Static magnetic field of 1.5-Tesla or 3-Tesla, only.
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40-T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode.

Under the scan conditions defined, the implant system is expected to produce a maximum temperature rise of 3.6°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant system extends approximately 30-mm from this system when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

Parameter	Condition
Nominal Values of Static Magnetic Field (T)	1.5-Tesla and 3.0-Tesla
Maximum Spatial Field Gradient (T/m and gauss/cm)	40-T/m (4,000-gauss/cm)
Type of RF Excitation	Circularly Polarized (CP) (i.e., Quadrature-Transmission)
Transmit RF Coil Information	There are no transmit RF coil restrictions. Accordingly, the following may be used: body transmit RF coil and all other RF coil combinations (i.e., body RF coil combined with any receive-only RF coil, transmit/receive head RF coil, transmit/receive knee RF coil, etc.)
Operating Mode of MR System	Normal Operating Mode
Maximum Whole Body Averaged SAR	2-W/kg (Normal Operating Mode)
Limits on Scan Duration	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks).
MR Image Artifact	The presence of this implant system produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant system is located in the area of interest.

To allow medical professionals to identify the specific medical devices a patient has, the MRI safety status of the medical devices, and the conditions for safe use in the MR environment for MR Conditional devices, BioHorizons recommends that clinicians provide the patient with the dental implant(s) and dental abutment(s) / bar(s) device specific peel-offs affixed to the patient card. Patient cards are available free of charge upon request from BioHorizons or for direct printing at <http://eifu.intra-lock.com>.

## COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with MDL and MILO prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) implant and/or abutment breakage; (3) infection requiring revision of the dental implant; (4) nerve damage that could cause permanent weakness, numbness, or pain; (5) histologic responses possibly involving macrophages and/or fibroblasts; (6) formation of fat emboli; (7) implant loosening requiring revision surgery; (8) maxillary sinus perforation; (9) labial or lingual plate perforation; and (10) bone loss possibly resulting in revision or removal.

**HANDLING AND STERILIZATION**

Always handle the device with powder-free gloves and avoid contact with hard objects that may damage the surface. Non-sterile device must be cleaned and sterilized. BioHorizons devices have not been validated for automated cleaning. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a container using a broad-spectrum cleaning agent such as Hu-Friedy’s Enzymax® per the manufacturer’s recommendations. Refer to the legal manufacturer’s instructions for use for detergent solution preparation.
- 2) Brush the device to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse device under running utility tap water.
- 4) Place device in the sterile tube filled with the prepared detergent solution and sonicate for two (2) minutes.
- 5) Thoroughly rinse device under running utility tap water.
- 6) Spray device with 70% isopropyl alcohol (IPA).
- 7) Blot device dry with clean lint free cloth.

For sterilization of non-sterile devices, place devices in an FDA cleared sterilization bag or wrap, and run through one of the following qualified sterilization cycles:

Sterilization Cycles				
<b>Reference:</b>	AAMI TIR12:2020	AAMI TIR12:2020	AAMI TIR12:2020	UK HTM 01-01 Part C:2016
<b>Type:</b>	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	30 minutes	20 minutes	20 minutes

It is recommended to include a 30-minute cool-down period before removing the device from the sterilization bag or wrap.

Devices to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

**SYMBOLS AND DESCRIPTIONS**

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic instructions for use.
	Manufacturer.
	Intra-Lock products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label.

Symbol	Symbol Description
	The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number.
	Lot/ batch number.
	Unique Device Identifier.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.
	Date of manufacture.
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician.
	European Union Authorized Representative.
	Do not use if package is damaged. Discard device and package.
	Medical Device.
	Non-sterile.
	Single sterile barrier system with protective packaging outside.
	Single sterile barrier system.
	Home.

Symbol	Symbol Description
 MR Conditional	Magnetic resonance warning: Device is MR conditional.
	United Kingdom Responsible Person.

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