English	Instructions for use: Dental Instruments
Español	Instrucciones de uso: Instrumentos Dentales
русский	Инструкция по применению: Зубоврачебные инструменты
Deutsch	Gebrauchsanweisung: Zahnmedizinische Instrumente
Français	Instructions d'utilisation: Instruments Dentaires
Italiano	Istruzioni per l'uso: Strumenti dentali
Português	Instruções de utilização: Instrumentos Dentários
Türkçe	Kullanım talimatları: Dental Aletler
简体中文	使用说明: 牙科器械
日本語	使用説明書: 歯科用機器
한국어	적용증: 치과 기구
العربية	إرشادات الاستخدام: أدوات طب الأسنان
Polski	Instrukcja obsługi: Narzędzia Dentystyczne
ČEŠTINA	Návod k použití: Zubařské nástroje
УКРАЇНСЬКА	Інструкція по застосуванню: Стоматологічні інструменти
Български	Инструкции за употреба: Дентални инструменти
Ελληνικά	Οδηγίες χρήσης: Οδοντιατρικά εργαλεία
Magyar	Használati utasítás: Fogászati eszközök
Română	Instrucțiuni de utilizare: instrumente dentare
Nederlands	Gebruiksaanwijzing: Tandheelkundige instrumenten
Svenska	Bruksanvisning: Dentalinstrument
Norsk	Bruksanvisning: Dentale instrumenter

**BioHorizons Implant Systems Inc.**  
2300 Riverchase Center  
Birmingham AL, 35244 USA  
TOLL FREE 888.246.8338  
TEL 205.967.7880  
FAX 205.870.0304  
www.biohorizons.com

L02036 Rev K JAN 2025

**ENGLISH**

This document supersedes all prior revisions. Original language is English.

This document applies to BioHorizons Dental Instruments and Accessories listed below. Every product packaging label contains a description of the enclosed device.

Devices in Scope	
Surgical kits and surgical trays	Screw drivers
Prosthetics kits and prosthetic trays	Implant try-ins
Surgical drills provided non-sterile	Surgical guides
Surgical drills provided sterile (single-use)	Master cylinders and pilot sleeves
Surgical taps provided non-sterile	Connectors and extenders
Surgical taps provided sterile (single-use)	Tapered depth gauges
Tissue punches and tissue groomers	Parallel and fixation pins
Bone profiling burs	Ratchets, tools and handles
Tapered ridge expanders	Carriers
Implant level drivers	Tray grommets
Abutment level drivers	Snap rings
Reusable scan bodies	

**DESCRIPTION**

BioHorizons instruments are used for dental implant procedures such as site development, implant placement, and implant restorations within the specific indications of each implant system. The label on each instrument contains important product information including whether the instrument is supplied sterile or non-sterile and if the instrument is for single-use only. Instruments and kits supplied non-sterile must be cleaned and sterilized prior to first use and each use thereafter. Refer to the following table for material(s) of devices with direct patient contact:

Devices with Direct Patient Contact	Material (main elements)
Surgical drills provided non-sterile	Stainless Steel (Iron, Chromium, Nickel)
Surgical drills provided sterile (single-use)	Stainless Steel (Iron, Chromium, Nickel)
Surgical taps provided non-sterile	Stainless Steel (Iron, Chromium, Nickel)
Surgical taps provided sterile (single-use)	Stainless Steel (Iron, Chromium, Nickel)
Tissue punches	Stainless Steel (Iron, Chromium, Nickel)
Tissue groomers	Stainless Steel (Iron, Chromium, Nickel) and 120/140 mesh synthetic diamond
Bone profiling burs	Stainless Steel (Iron, Chromium, Nickel)
Tapered ridge expanders	Stainless Steel (Iron, Chromium, Nickel)
Implant level drivers	Stainless Steel (Iron, Chromium, Nickel)
Implant try-ins	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Reusable scan bodies	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)

**INTENDED USE**

BioHorizons instruments and surgical kits are intended for use in site development, placement, and restoration of BioHorizons implants and associated components.

### **INDICATIONS FOR USE**

BioHorizons instruments and surgical kits are indicated for use in site development, placement, and restoration of BioHorizons implants and associated components. The BioHorizons Single Use Drills are indicated for use with the BioHorizons Tapered Internal and Laser-Lok 3.0 implants only. No attempt should be made to re-use or re-sterilize the drills. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

### **CONTRAINDICATIONS**

BioHorizons instruments should not be used with patients that have allergies to the specific materials used, including stainless steel and titanium alloy.

### **PATIENT POPULATION**

BioHorizons dental implant systems are intended for use in skeletally mature, non-pediatric edentulous or partially edentulous patients as long as the defined contraindications are nonapplicable.

### **INTENDED USERS**

BioHorizons implant systems are only intended for use by licensed health care professionals, more specifically, BioHorizons implant systems are intended to be used by trained dentists, surgeons and dental technicians in a standard dental surgical setting, which may range from General Dentists' practice offices to Maxillofacial Surgical Operating Rooms, as well as laboratories for dental processes. Use of these products requires specialized knowledge and experience in implant dentistry. BioHorizons implant systems are marked and labeled as a medical device (MD) and are Rx only.

### **DIRECTIONS FOR USE**

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the instructions pertaining to BioHorizons products.

Internal Abutment-level Driver, Hex-chuck Handpiece: Use only with compatible W&H Hexagon Chucking System handpieces. Drills and Taps: Should be replaced when wear is noticed, such as a decrease in cutting efficiency or when signs of discoloration appear. Cutting instruments should be replaced after approximately 12 to 20 osteotomy cycles, depending on bone density. BioHorizons recommends use of a drill usage chart to track drill use and to ensure drills are replaced as directed. Single Use Drills: Insert the drill into a handpiece and twist protective tubing off prior to use. CGS Drills and Guides: Care must be taken to ensure the CGS Drills follow the angle of the guide handle inserts or sleeves to prevent drill binding. Refer to the drill report for proper drill selection and usage. Site Dilators: Push the dilator into the bone while rotating it. In firmer bone, it may be necessary to use a mallet to tap the dilator to proper depth. Laser-Lok Protective Sleeves: Secure Titanium Base Abutments onto the sleeve side labeled Hybrid or Custom Titanium Abutments onto the side labeled Full Ti to protect Laser-Lok region during bonding and polishing.

The BioHorizons Simple Solutions Replica Lab Tools are hand instruments to aid in the fabrication of burnout wax copings for Simple Solution restorative protocols only. They may be used in the dental suite or lab and mimic the geometry and size of the corresponding Simple Solution Abutment and/or Replica without retentive snap feature, allowing the wax coping(s) to be removed without deformation after hardening. The replica lab tools are manufactured of titanium alloy.

### **WARNINGS AND PRECAUTIONS**

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Refer to OSHA standard 29CFR1910.1030 prior to cleaning and sterilization. Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at [www.biohorizons.com](http://www.biohorizons.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

Prosthetic Instruments - Only use a torque wrench for final installation of prosthetic components requiring a specified installation torque. Finger-tightening of prosthetic components can result in insufficient torque and eventual loosening of the component. Over-tightening prosthetic components could break the component or spin the mating implant.

### **COMPLICATIONS AND ADVERSE EFFECTS**

The risks and complications with instruments, prosthetics, and implants include, but are not limited to: (1) allergic reaction(s) to instrument, implant and/or abutment material; (2) breakage of implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) abutment screw and/or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; (11) bone loss possibly resulting in revision or removal; and (12) respiratory and/or intestinal tract damage as a result of component aspiration.

### **HANDLING AND STERILIZATION**

Always handle the device with powder-free gloves and avoid contact with hard objects that may damage the surface. If the device is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove device from the package only after the correct size has been determined. For non-sterile devices, remove and discard any shipping material before initial processing. Non-sterile single use devices and devices intended for reuse must be cleaned and sterilized prior to use. BioHorizons devices have not been validated for automated cleaning. The following cleaning protocol must be used:

#### **Instruments (Excluding Drills):**

- 1) Assembled instruments (e.g., ODSecure Cap Insert/Extract Tool (ODSCT)) should be disassembled before each cleaning and sterilization cycle to avoid debris encapsulation, material discoloration and/or inappropriate drying of components.
- 2) Remove all the instruments from the instrumentation trays and disassemble the trays before each cleaning.
- 3) Prepare a detergent solution in a container using a broad-spectrum cleaning agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations. Refer to the legal manufacturer's instructions for use/preparation of detergent solution.
- 4) Fully immerse the devices in the prepared detergent solution and brush them to remove visible soil.
- 5) Remove the devices from the detergent solution and rinse them thoroughly under running utility tap water.
- 6) Place devices in an ultrasonic unit filled with fresh detergent solution and sonicate for ten (10) minutes.
- 7) Remove the devices from the detergent solution and thoroughly rinse them under running utility tap water.
- 8) Wipe the devices with 70% isopropyl alcohol (IPA).
- 9) Blot the devices dry with clean lint free cloths and allow them to air dry completely.

#### **Drills:**

- 1) Remove all the drills from the instrumentation trays.

- 2) Prepare a detergent solution in a container using a broad-spectrum cleaning agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations. Refer to the legal manufacturer's instructions for use/preparation of detergent solution.
- 3) Drills shall be brushed to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 4) Thoroughly rinse drills thoroughly under running utility tap water.
- 5) Drills shall be placed in a sterile tube filled with fresh detergent solution and sonicated for ten (10) minutes.
- 6) Thoroughly rinse drills under running utility tap water.
- 7) Spray the drills with 70% IPA.
- 8) Blot the drills dry with clean lint free cloths.

For AS123 Hand Unit (300-100): blow pressurized air through one of the four drain holes to purge. For CGS Quick Connect Handle (CGS-QCH): blow pressurized air through one end to purge. For both instruments (300-100 and CGS-QCH): use a medical instrument lubricant such as STERIS' Hinge Free® after each use. Refer to the legal manufacturer's instructions for use/preparation of lubricant.

Instruments that are unable to be cleaned, are discolored, do not properly interface with mating components, and/or do not articulate as designed should be disposed of.

If applicable, return the instruments to the appropriate locations in the instrumentation tray. The BioHorizons Prosthetic Tray, Small (PROS1500) should be cleaned and sterilized empty. For sterilization, place the cleaned device in an FDA cleared sterilization bag or wrap, and run through one of the following qualified sterilization cycles:

Sterilization Cycles			
<b>Reference:</b>	AAMI TIR12:2020	AAMI TIR12:2020	UK HTM 01-01 PartC:2016
<b>Type:</b>	Gravity Steam	Prevacuum Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	20 minutes	20 minutes

The following qualified sterilization cycle may also be used for individual instruments, Tapered Short Guided Kit (TSGKIT), Tapered Short Guided Tray (TSGT), and empty BioHorizons Prosthetic Tray, Small (PROS1500):

Sterilization Cycle	
<b>Reference:</b>	AAMI TIR12:2020
<b>Type:</b>	Gravity Steam
<b>Exposure Time and Temperature:</b>	30 minutes at 121°C (250°F)
<b>Minimum Dry Time:</b>	30 minutes

The following qualified sterilization cycle may be used for the INSERT END of ODSecure Cap Insert/Extract Tool (ODSCT):

Sterilization Cycle	
<b>Type:</b>	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	3 minutes at 135°C (275°F)
<b>Minimum Dry Time:</b>	16 minutes
<b>Pressure:</b>	2.06 bar

The following qualified sterilization cycle may be used for the EXTRACT END and BLUE HANDLE of ODSecure Cap Insert/Extract Tool (ODSCT) and ODSecure Impression Coping Tool (ODSICT):

Sterilization Cycle	
Type:	Prevacuum Steam
Exposure Time and Temperature:	20 minutes at 121°C (250°F)
Minimum Dry Time:	15 minutes
Pressure:	1.1 bar

It is recommended to include a 30-minute cool-down period before removing the devices from the sterilization bag or wrap.

**Attention!** Improper cleaning may lead to inadequate sterilization. Failure to completely dry instruments during autoclaving may leave moisture and cause discoloration and oxidation. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the instruments. Periodic testing, cleaning, and calibration of the autoclave equipment is recommended to ensure the unit remains in proper working order.

Devices to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

### SYMBOLS AND DESCRIPTIONS

The symbols table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic instructions for use
	Manufacturer
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number
	Lot/ batch number
	Unique Device Identifier
	Do not re-use
	Do not re-sterilize

Symbol	Symbol Description
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician
	European Union Authorized Representative
	Do not use if package is damaged. Discard device and package.
	Medical Device
	Non-sterile
	Single sterile barrier system with protective packaging outside
	Single sterile barrier system
	Home
 MR Conditional	Magnetic resonance warning: Device is MR conditional
	United Kingdom Responsible Person.

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