<b>English</b>	Instructions for use: Conical Dental Instruments
<b>Español</b>	Instrucciones de uso: instrumentos dentales cónicos
<b>Русский</b>	Инструкция по применению. Стоматологические инструменты для конических имплантатов
<b>Deutsch</b>	Gebrauchsanweisung: Konische Zahninstrumente
<b>Français</b>	Instructions d'utilisation: Implants dentaires coniques
<b>Italiano</b>	Istruzioni per l'uso: strumenti dentali conici
<b>Português</b>	Instruções de utilização: instrumentos dentários cónicos
<b>Türkçe</b>	Kullanım talimatları: Konik Dental Aletler
<b>العربية</b>	تعليمات الاستخدام: أدوات طب الأسنان المخروطية
<b>Polski</b>	Instrukcja użycia: Stożkowe instrumenty stomatologiczne
<b>Українська</b>	Інструкції з використання: конічні стоматологічні інструменти
<b>Български</b>	Инструкции за употреба: Конични дентални инструменти
<b>Česky</b>	Návod k použití: Kónické zubní nástroje
<b>Ελληνικά</b>	Οδηγίες χρήσης: Κωνικά οδοντιατρικά εργαλεία
<b>Magyar</b>	Használati utasítás: Kúpos fogászati műszerek
<b>Română</b>	Instrucțiuni de utilizare: Instrumente dentare conice
<b>Nederlands</b>	Gebruiksaanwijzing: Conische tandheelkundige instrumenten
<b>Svenska</b>	Bruksanvisning: Koniska dentalinstrument
<b>Norsk</b>	Bruksanvisning: Koniske tanninstrumenter

**BioHorizons Implant Systems Inc.**  
2300 Riverchase Center  
Birmingham AL, 35244 USA  
TOLL FREE 888.246.8338  
TEL 205.967.7880  
FAX 205.870.0304  
[www.biohorizons.com](http://www.biohorizons.com)

L02086 Rev B FEB 2025

**ENGLISH**

This document supersedes all prior revisions. Original language is English.

This document applies to the BioHorizons Conical Dental Instruments listed below. Every product packaging label contains a description of the enclosed device.

<b>Devices in Scope</b>	
Screw drivers	Analogs and replicas
Bone profiling burs	Laboratory Components
Implant level drivers	Ratchets, tools and handles
Multi-Unit Abutment Try-Ins	Snap rings
Reusable scan bodies	Impression Components

**DESCRIPTION**

BioHorizons Conical instruments are used for dental implant procedures such as site development, implant placement, and implant restorations within the specific indications of each Conical dental implant system. The label on each instrument contains important product information including whether the instrument is supplied sterile or non-sterile and if the instrument is for single-use only. Instruments supplied non-sterile must be cleaned and sterilized prior to first intraoral use and each use thereafter. Refer to the following table for material(s) of devices with direct patient contact:

<b>Devices with Direct Patient Contact</b>	<b>Material (main elements)</b>
Bone profiling burs	Stainless Steel (Iron, Chromium, Nickel)
Implant level drivers	Stainless Steel (Iron, Chromium, Nickel)
Multi-Unit Abutment Try-Ins	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
Reusable scan bodies	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)

**INTENDED USE**

BioHorizons instruments and surgical kits are intended for use in site development, placement, and restoration of BioHorizons implants and associated components.

**INDICATIONS FOR USE**

BioHorizons Conical instruments are indicated for use in site development, placement, and restoration of BioHorizons Conical dental implants and associated components.

**CONTRAINDICATIONS**

BioHorizons Conical instruments should not be used with patients that have allergies to the specific materials used, including stainless steel and titanium alloy.

## **PATIENT POPULATION**

BioHorizons dental implant systems are intended for use in skeletally mature, non-pediatric edentulous or partially edentulous patients as long as the defined contraindications are nonapplicable.

## **INTENDED USERS**

BioHorizons implant systems are only intended for use by licensed health care professionals, more specifically, BioHorizons implant systems are intended to be used by trained dentists, surgeons and dental technicians in a standard dental surgical setting, which may range from General Dentists' practice offices to Maxillofacial Surgical Operating Rooms, as well as laboratories for dental processes. Use of these products requires specialized knowledge and experience in implant dentistry. BioHorizons implant systems are marked and labeled as a medical device (MD) and are Rx only.

## **DIRECTIONS FOR USE**

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the instructions pertaining to BioHorizons products.

Bone Profilers: Should be replaced when wear is noticed, such as a decrease in cutting efficiency or when signs of discoloration appear. Bone Profilers should be replaced after approximately 12 to 20 osteotomy cycles, depending on bone density.

## **WARNINGS AND PRECAUTIONS**

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Refer to OSHA standard 29CFR1910.1030 prior to cleaning and sterilization. Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at [www.biohorizons.com](http://www.biohorizons.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

## **COMPLICATIONS AND ADVERSE EFFECTS**

The risks and complications with instruments include, but are not limited to: (1) allergic reaction(s) to instrument material; (2) breakage of the implant required to be explanted using clinician judgement; (3) infection requiring revision of the dental implant; (4) nerve damage that could cause permanent weakness, numbness, or pain; (5) histologic responses possibly involving macrophages and/or fibroblasts; (6) formation of fat emboli; (7) maxillary sinus perforation; (8) labial or lingual plate perforation; (9) bone loss possibly resulting in revision or removal; and (10) respiratory and/or intestinal tract damage as a result of component aspiration/ingestion.

## **HANDLING AND STERILIZATION**

Always handle instruments with powder-free gloves and avoid contact with hard objects that may damage the surface.

For non-sterile devices, remove and discard any shipping material before initial processing of instruments.

BioHorizons instruments have not been validated for automated cleaning.

Non-sterile single-use instruments and reusable instruments that are patient contacting must be cleaned and sterilized prior to use. The following cleaning protocol must be used:

- 1) Assembled instruments (e.g., Guided Implant Driver/Screw) should be disassembled before each cleaning and sterilization cycle to avoid debris encapsulation, material discoloration and/or inappropriate drying of components.
- 2) Remove all the instruments from the instrumentation trays and disassemble the trays before each cleaning.
- 3) Prepare a detergent solution in a container using a broad-spectrum cleaning agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations. Refer to the legal manufacturer's instructions for use/preparation of detergent solution.
- 4) Fully immerse the devices in the prepared detergent solution and brush them to remove visible soil.
- 5) Remove the devices from the detergent solution and rinse them thoroughly under running utility tap water.
- 6) Place devices in an ultrasonic unit filled with fresh detergent solution and sonicate for ten (10) minutes.
- 7) Remove the instruments from the detergent solution and thoroughly rinse them under running utility tap water.
- 8) Wipe the devices with 70% isopropyl alcohol (IPA).
- 9) Blot the instruments dry with clean lint free cloths and allow them to air dry completely.

Instruments that are unable to be cleaned, are discolored, do not properly interface with mating components, and/or do not articulate as designed should be disposed of.

If applicable, return the instruments to the appropriate locations in the instrumentation tray. For sterilization, place the cleaned instrument in an FDA cleared sterilization bag or wrap, and run through one of the following qualified sterilization cycles:

	Sterilization Cycles			
<b>Reference:</b>	AAMI TIR12:2020	AAMI TIR12:2020	AAMI TIR12:2020	UK HTM 01-01 PartC:2016
<b>Type:</b>	Gravity Steam	Prevacuum Steam	Gravity Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	30 minutes at 121°C (250°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	20 minutes	30 minutes	20 minutes

It is recommended to include a 30-minute cool-down period before removing the devices from the sterilization bag or wrap.

**Attention!** Improper cleaning may lead to inadequate sterilization. Failure to completely dry instruments during autoclaving may leave moisture and cause discoloration and oxidation. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the instruments. Periodic testing, cleaning, and calibration of the autoclave equipment is recommended to ensure the unit remains in proper working order.

Devices to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

## SYMBOLS AND DESCRIPTIONS

The symbols table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic instructions for use
	Manufacturer

Symbol	Symbol Description
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number
	Lot/ batch number
	Unique Device Identifier
	Do not re-use
	Do not re-sterilize
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
<b>Rx Only</b>	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician
	European Union Authorized Representative
	Do not use if package is damaged. Discard device and package.
	Medical Device
	Non-sterile
	Single sterile barrier system with protective packaging outside
	Single sterile barrier system
	Home

Symbol	Symbol Description
 MR Conditional	Magnetic resonance warning: Device is MR conditional
	United Kingdom Responsible Person.

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