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L02037CAen Rev B MAR 2025

This document supersedes all prior revisions. Original language is English.

This document applies to BioHorizons Dental Abutments listed below. Every product packaging label contains a description of the enclosed device.

<b>Devices in Scope</b>
BioHorizons Custom Ti abutments
BioHorizons Custom milled abutments
BioHorizons Titanium base abutments
BioHorizons Hybrid base abutments

## DESCRIPTION

BioHorizons CAD/CAM titanium blanks (Custom Ti and Custom Milled) and CAD/CAM titanium bases (Titanium Base and HybridBase) are intended for the restoration of BioHorizons Tapered Internal and Internal dental implants and Zimmer® Dental ScrewVent® and Tapered ScrewVent® dental implants within the specific indications of each implant system. The label on each prosthetic package contains important product information including an indication that the prosthetic is supplied non-sterile. Refer to the following table for material(s) of implantable devices:

<b>Implantable Devices</b>	<b>Material (main elements)</b>
BioHorizons Custom Ti abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons Custom milled abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons Titanium base abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons Hybrid base abutments	Ti-6AL-4V ELI (Titanium, Aluminum, Vanadium)

## INTENDED USE

BioHorizons abutments are intended for the restoration of BioHorizons dental implants in the mandible or maxilla, for single tooth replacement or for fixed bridgework and dental retention, within the specific indications of each implant system.

## INDICATIONS FOR USE

BioHorizons CAD/CAM Abutments are dental abutments placed onto a dental implant to provide support for dental prosthetic restorations. The abutments include: 1) Titanium abutment blanks with a pre-machined implant connection where the upper portion may be custom-milled in accordance with a patient-specific design using CAD/CAM techniques; and 2) Titanium bases with a pre-machined implant connection upon which a CAD/CAM designed superstructure may be fitted to complete a two-piece dental abutment. The abutments include an abutment screw for fixation to the underlying implant. The abutments may be used for single-unit (single-tooth) or multiple-unit (bridges and bars) restorations and are compatible for use with BioHorizons Internal and Tapered Internal implant systems and Zimmer® Dental Screw-Vent® and Tapered Screw-Vent® implants with 3.5mm, 4.5mm and 5.7mm internal hex-connection mating platform diameters. All digitally designed abutments and/or copings for use with BioHorizons CAD/CAM Abutments are intended to be sent to a BioHorizons-validated milling center for manufacture. BioHorizons abutments designed using CAD/CAM techniques must fulfill the BioHorizons allowable range of design parameters.

## CONTRAINdicATIONS

BioHorizons prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity.

## PATIENT POPULATION

BioHorizons dental implant systems are intended for use in skeletally mature, non-pediatric edentulous or partially

edentulous patients as long as the defined contraindications are nonapplicable.

## INTENDED USERS

BioHorizons implant systems are only intended for use by licensed health care professionals, more specifically, BioHorizons implant systems are intended to be used by trained dentists, surgeons and dental technicians in a standard dental surgical setting, which may range from General Dentists' practice offices to Maxillofacial Surgical Operating Rooms, as well as laboratories for dental processes. Use of these products requires specialized knowledge and experience in implant dentistry. BioHorizons implant systems are marked and labeled as a medical device (MD) and are Rx only.

## DIRECTIONS FOR USE

A custom abutment design file is created by the customer using a 3Shape abutment library. The abutment design file is converted to a stereolithography (.stl) file using 3Shape software. The .stl file is converted to a numerical control (.nc) file using an appropriate software application (e.g. SUM3D). After the custom abutment design is uploaded, the milling machine software executes the necessary commands to mill the final abutment. The digitally designed file is to be sent to a registered and listed BioHorizons contract manufacturer for milling. Visit [www.biohorizons.com](http://www.biohorizons.com) for a list of validated milling centers. BioHorizons CAD/CAM Abutments are compatible with commercially available dental CAD/CAM systems, such as 3Shape.

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions for Use that accompany BioHorizons products. The Tissue Level Hybrid Abutment Base is used for single-unit restorations where the underlying implant is placed supra-crestally and a single-stage protocol is desired. Abutment screws must be torqued to 30Ncm to prevent screw loosening. BioHorizons Custom Milled Abutments, Hybrid Abutment Bases and 3.0mm Custom Ti Abutments contain a LIGHT BLUE abutment screw specifically designed for use in those abutments. Use of the LIGHT BLUE abutment screws with incompatible abutments may result in abutment and/or abutment screw failure. Restorations utilizing angled screw channels are only compatible with YELLOW Precision Angled Screws. Angled screw channels can have a maximum angulation of 15° off the vertical axis of the abutment.

Design limitation parameters for the Custom Ti and Custom Milled Abutments are included in the following table.

Feature	Custom Ti Platforms				Custom Milled Platforms	
	Ø3.0mm		Ø3.5mm, Ø4.5mm, Ø5.7mm		Ø3.0mm, Ø3.5mm, Ø4.5mm, Ø5.7mm	
	Min	Max	Min	Max	Min	Max
Wall thickness	0.3mm	3.8mm	0.3mm	3.6mm	0.3mm	3.7mm
Gingival margin diameter	3.0mm	9.7mm	3.0mm	9.7mm	3.3mm, 3.8mm, 4.6mm, 5.8mm respectively	9.5mm
Gingival margin height	0mm	12.1mm	0mm	12.1mm	0.3mm, 0.2mm, 0.1mm, 0.1mm respectively	12.3mm, 12.5mm, 12.6mm, 12.5mm respectively
Length of the abutment post	*	12.1mm	*	12.1mm	*	12.3mm, 12.5mm, 12.6mm, 12.5mm respectively
Taper of the abutment	0°	90°	0°	90°	0°	90°
Angulation	0°	30°	0°	30°	0°	30°

\* Minimum allowable abutment post dimensions are 3.0mm for screw retained abutment and 4.0mm for cement retained abutment.

Titanium base abutments - BioHorizons recommends the use of Metoxit Z-CAD® dental ceramic (or equivalent) to fabricate the custom zirconia superstructure and RelyX Unicem (or equivalent) to affix the zirconia superstructure to the titanium base.

Design limitation parameters for the custom zirconia superstructure portion of the Titanium Base and Hybrid Base Abutments are included in the following table.

Feature	Titanium Base Platforms		Hybrid Base Platforms	
			Ø3.0mm <sup>†</sup> , Ø3.5mm, Ø4.5mm, Ø5.7mm	
	Min	Max	Min	Max
Wall thickness	0.4mm	**	0.4mm	**
Gingival margin diameter	5.1mm	**	3.8mm, 3.8mm, 4.6mm, 5.8mm respectively & 5.1mm, 6.1mm, 7.4mm respectively for Tissue Level	**
Length of the abutment post	3.9mm	**	4.0mm	**
Taper	Variable based on anatomy	**	Variable based on anatomy	**
Angulation of the final titanium base	0°	20°	0°	20°

<sup>†</sup> Not available for Tissue Level.

\*\*Feature maximum sizes are based on patient anatomy and clinician discretion.

Secure Laser-Lok abutments onto the Laser-Lok Protective Sleeve to protect Laser-Lok during handling.

## WARNINGS AND PRECAUTIONS

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of BioHorizons prosthetic components.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics.

If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the modifier of the implant/abutment interface a medical device company subject to FDA rules and regulations.

Titanium base abutments - Altering the size or shape of the titanium base abutment may affect proper fit with the restoration. Leaving bonding agent or particulate in the Laser-Lok surface may affect the connective tissue attachment.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

Devices should not be used in patients with a known allergy or sensitivity to the device material(s).

Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at [www.biohorizons.com](http://www.biohorizons.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

For implantable devices, the summary of safety and clinical performance (SSCP) according to Article 32 of Regulation (EU) 2017/745 can be found in the European database for medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed>.

Implantable Device	Basic UDI-DI Number
BioHorizons Custom Ti abutments	08472360IIBWETABUT007J6
BioHorizons Custom milled abutments	08472360IIBWETABUT007J6
BioHorizons Titanium base abutments	08472360IIBWETABUT007J6
BioHorizons Hybrid base abutments	08472360IIBWETABUT007J6

## MAGNETIC RESONANCE (MR) SAFETY INFORMATION

MR-Conditional

Non-clinical testing has demonstrated BioHorizons Implant Systems are MR conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 4,000 G/cm (40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined above, the implant system is expected to produce a maximum temperature rise of 3.6 °C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant system extends approximately 30 mm from this implant when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

### COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) abutment screw and/or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

### HANDLING AND STERILIZATION

Always handle the device with powder-free gloves and avoid contact with hard objects that may damage the surface. After abutment preparation and before final placement in the restorative site, CAD/CAM abutments must be cleaned and sterilized. BioHorizons devices have not been validated for automated cleaning. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a container using a broad-spectrum cleaning agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations. Refer to the legal manufacturer's instructions for use/preparation for detergent solution.
- 2) Brush the device to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse device under running utility tap water.
- 4) Place device in the sterile tube filled with the prepared detergent solution and sonicate for two (2) minutes.
- 5) Thoroughly rinse device under running utility tap water.
- 6) Spray device with 70% isopropyl alcohol (IPA).
- 7) Blot device dry with clean lint free cloth.

For sterilization, place final-prepped and cleaned abutment in an FDA cleared sterilization bag or wrap, and run through one of the following qualified sterilization cycles:

Sterilization Cycles				
Reference:	AAMI TIR12:2020	AAMITIR12:2020	AAMI TIR12:2020	UK HTM 01-01 Part C:2016
Type:	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
Exposure Time and Temperature:	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
Minimum Dry Time:	30 minutes	30 minutes	20 minutes	20 minutes

It is recommended to include a 30-minute cool-down period before removing the device from the sterilization bag or wrap.

Devices to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

## SYMBOLS AND DESCRIPTIONS

The symbols table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution.
	Electronic instructions for use.
	Manufacturer.
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label. The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number.
	Lot/ batch number.
	Unique Device Identifier.
	Do not re-use.
	Do not re-sterilize.
	Use-by-date.
	Sterile by gamma irradiation.

Symbol	Symbol Description
	Date of manufacture.
<b>Rx Only</b>	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician.
<b>EC REP</b>	European Union Authorized Representative
	Do not use if package is damaged. Discard device and package.
<b>MD</b>	Medical Device.
<b>Non-sterile</b>	Non-sterile.
	Single sterile barrier system with protective packaging outside.
	Single sterile barrier system.
	Home Position.
 MR Conditional	Magnetic resonance warning: Device is MR conditional.
<b>UK RP</b>	United Kingdom Responsible Person.

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc.

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Zimmer®, Dental ScrewVent® et Tapered ScrewVent® sont des marques déposées de Zimmer Dental, Inc.