English	Instructions for use: Custom abutments
Español	Instrucciones de uso: Pilares a medida
Русский	Инструкция по применению: Индивидуальные абатменты
Deutsch	Gebrauchsanweisung: Abutments nach Maß
Français	Instructions d'utilisation: Piliers prothétiques personnalisés
Italiano	Istruzioni per l'uso: Monconi personalizzati
Português	Instruções de utilização: Pilares personalizados
Türkçe	Kullanım talimatları: Özel Abutmentler
简体中文	使用说明：自定义基牙
日本語	使用説明書：カスタム アバットメント
한국어	적응증: 맞춤형 지대치
عربي	إرشادات الاستخدام: دعائم مخصصة
Polski	Instrukcja obsługi: Zaczepy specjalne
Čeština	Návod k použití: Abutmenty na míru
Українська	Інструкція по застосуванню: Індивідуальні абатменти
Български	Инструкции за употреба: Персонализирани опори
Ελληνικά	Οδηγίες χρήσης: Εξατομικευμένα κολοβώματα
Magyar	Használati utasítás: Egyéni kiosztások
Română	Instrucțiuni de utilizare: Bonturi personalizate
Nederlands	Gebruiksaanwijzing: op maat gemaakte abutments
Svenska	Bruksanvisning: Anpassade distanser
Norsk	Bruksanvisning: Spesialtilpassede distanser

BioHorizons Implant Systems Inc.

2300 Riverchase Center
Birmingham AL, 35244 USA
TOLLFREE 888.246.8338
TEL 205.967.7880
FAX 205.870.0304
www.biohorizons.com

This document supersedes all prior revisions. Original language is English. The instructions for use are available free of charge in printed form upon request from BioHorizons or your local distributor.

This document applies to BioHorizons Dental Abutments listed below. Every product packaging label contains a description of the enclosed device.

Devices in Scope
BioHorizons Custom Ti abutments
BioHorizons Custom milled abutments
BioHorizons Titanium base abutments
BioHorizons Hybrid base abutments
BioHorizons Chairside Ti-Bases

1.0 DESCRIPTION

BioHorizons CAD/CAM titanium blanks (Custom Ti and Custom Milled) and CAD/CAM titanium bases (Titanium Base and Hybrid Base) are intended for the restoration of BioHorizons Tapered Internal and Internal dental implants and Zimmer® Dental ScrewVent® and Tapered ScrewVent® dental implants within the specific indications of each implant system. The label on each prosthetic package contains important product information including an indication that the prosthetic is supplied non-sterile. Refer to the following table for material(s) of implantable devices:

Implantable Devices	Material (main elements)
BioHorizons Custom Ti abutments	Ti-6Al-4V ELI (Titanium, Aluminum, Vanadium)
BioHorizons Custom milled abutments	
BioHorizons Titanium base abutments	
BioHorizons Hybrid base abutments	
BioHorizons Chairside Ti-Base titanium components	

The label on each prosthetic package contains important product information including an indication that the prosthetic is supplied non-sterile. Prosthetic connection color coding is defined as shown in the table below:

Internal Hex Prosthetic Color Coding:	Gray	Yellow	Green	Blue
Internal Hex Prosthetic Size:	3.0mm	3.5mm	4.5mm	5.7mm

2.0 INTENDED USE

BioHorizons abutments are intended for the restoration of BioHorizons dental implants in the mandible or maxilla, for single tooth replacement or for fixed bridgework and dental retention, within the specific indications of each implant system.

3.0 CLINICAL BENEFITS

BioHorizons dental devices allow for the replacement of missing teeth and / or the restoration of masticatory function and the corresponding improvement in quality of life.

4.0 INDICATIONS FOR USE

BioHorizons CAD/CAM Abutments are dental abutments placed onto a dental implant to provide support for dental prosthetic restorations. The abutments include: 1) Titanium abutment blanks with a pre-machined implant connection where the upper portion may be custom-milled in accordance with a patient-specific design using CAD/CAM techniques; and 2) Titanium bases with a pre-machined implant connection upon which a CAD/CAM designed superstructure may be fitted to complete a two-piece dental abutment. The abutments include an abutment screw for fixation to the underlying implant. The abutments may be used for single-unit (single-tooth) or multiple-unit (bridges and bars) restorations and are compatible for use with BioHorizons Internal and Tapered Internal implant systems and Zimmer® Dental Screw-Vent® and Tapered Screw-Vent® implants with 3.5mm, 4.5mm and 5.7mm internal hex-connection mating platform diameters. All digitally designed abutments

and/or copings for use with BioHorizons CAD/CAM Abutments are intended to be sent to a BioHorizons-validated milling center for manufacture. BioHorizons abutments designed using CAD/CAM techniques must fulfill the BioHorizons allowable range of design parameters.

BioHorizons Chairside Ti-Base prosthetic components connected to the endosseous dental implants are intended for use as an aid in prosthetic rehabilitations of the maxillary or mandibular arch to provide support for prosthetic restorations.

BioHorizons 3.0 Chairside Ti-Bases are limited for use as an artificial root structure for single tooth replacement of mandibular central and lateral incisors and maxillary lateral incisors when used with 3.4mm diameter implants.

All digitally designed abutments and/or copings for use with Chairside Ti-Bases are to be designed using Sirona CEREC® Software and manufactured using a Sirona CEREC® MC XL milling unit.

5.0 CONTRAINDICATIONS

BioHorizons prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity.

BioHorizons 3.0 Chairside Ti-Bases are contraindicated for use with 3.0mm diameter implants. Tapered Tissue Level implants are not intended for use with the Chairside Ti-Bases.

6.0 PATIENT POPULATION

BioHorizons dental implant systems are intended for use in skeletally mature, non-pediatric edentulous or partially edentulous patients as long as the defined contraindications are nonapplicable.

7.0 INTENDED USERS

BioHorizons implant systems are only intended for use by licensed health care professionals, more specifically, BioHorizons implant systems are intended to be used by trained dentists, surgeons and dental technicians in a standard dental surgical setting, which may range from General Dentists' practice offices to Maxillofacial Surgical Operating Rooms, as well as laboratories for dental processes. Use of these products requires specialized knowledge and experience in implant dentistry. BioHorizons implant systems are marked and labeled as a medical device (MD) and are Rx only.

8.0 DIRECTIONS FOR USE

- 8.1** A custom abutment design file is created by the customer using a 3Shape abutment library. The abutment design file is converted to a stereolithography (.stl) file using 3Shape software. The .stl file is converted to a numerical control (.nc) file using an appropriate software application (e.g. SUM3D). After the custom abutment design is uploaded, the milling machine software executes the necessary commands to mill the final abutment. The digitally designed file is to be sent to a registered and listed BioHorizons contract manufacturer for milling. Visit www.biohorizons.com for a list of validated milling centers. BioHorizons CAD/CAM Abutments are compatible with commercially available dental CAD/CAM systems, such as 3Shape.
- 8.2** Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions for Use that accompany BioHorizons products.
- 8.3** The Tissue Level Hybrid Abutment Base is used for single-unit restorations where the underlying implant is placed supra-crestally and a single-stage protocol is desired.
- 8.4** Abutment screws must be torqued to 30 Ncm to prevent screw loosening.
- 8.5** BioHorizons Custom Milled Abutments, Hybrid Abutment Bases and 3.0mm Custom Ti Abutments contain a LIGHT BLUE abutment screw specifically designed for use in those abutments. Use of the LIGHT BLUE abutment screws with incompatible abutments may result in abutment and/or abutment screw failure.
- 8.6** Restorations utilizing angled screw channels are only compatible with YELLOW Precision Angled Screws. Angled screw channels can have a maximum angulation of 15° off the vertical axis of the abutment.
- 8.7** Custom Ti and Custom Milled Abutments

8.7.1 Design Limitations:

Design limitation parameters for the Custom Ti and Custom Milled Abutments are included in the following table.

Feature	Custom Ti Platforms				Custom Milled Platforms	
	Ø3.0mm		Ø3.5mm, Ø4.5mm, Ø5.7mm		Ø3.0mm, Ø3.5mm, Ø4.5mm, Ø5.7mm	
	Min	Max	Min	Max	Min	Max
Wall thickness	0.3mm	3.8mm	0.3mm	3.6mm	0.3mm	3.7mm
Gingival margin diameter	3.0mm	9.7mm	3.0mm	9.7mm	3.3mm, 3.8mm, 4.6mm, 5.8mm respectively	9.5mm
Gingival margin height	0mm	12.1mm	0mm	12.1mm	0.3mm, 0.2mm, 0.1mm, 0.1mm respectively	12.3mm, 12.5mm, 12.6mm, 12.5mm respectively
Length of the abutment post	*	12.1mm	*	12.1mm	*	12.3mm, 12.5mm, 12.6mm, 12.5mm respectively
Taper of the abutment	0°	90°	0°	90°	0°	90°
Angulation	0°	30°	0°	30°	0°	30°

* Minimum allowable abutment post dimensions are 3.0mm for screw retained abutment and 4.0mm for cement retained abutment.

8.8 Titanium Base Abutments

8.8.1 Design:

Titanium base abutments - BioHorizons recommends the use of Metoxit Z-CAD® dental ceramic (or equivalent) to fabricate the custom zirconia superstructure and RelyX Unicem (or equivalent) to affix the zirconia superstructure to the titanium base.

8.8.2 Design Limitations:

Design limitation parameters for the custom zirconia superstructure portion of the Titanium Base and Hybrid Base Abutments are included in the following table.

Feature	Titanium Base Platforms		Hybrid Base Platforms	
			Ø3.0mm†, Ø3.5mm, Ø4.5mm, Ø5.7mm	
	Min	Max	Min	Max
Wall thickness	0.4mm	**	0.4mm	**
Gingival margin diameter	5.1mm	**	3.8mm, 3.8mm, 4.6mm, 5.8mm respectively & 5.1mm, 6.1mm, 7.4mm respectively for Tissue Level	**
Length of the abutment post	3.9mm	**	4.0mm	**
Taper	Variable based on anatomy	**	Variable based on anatomy	**
Angulation of the final titanium base	0°	20°	0°	20°

† Not available for Tissue Level.

**Feature maximum sizes are based on patient anatomy and clinician discretion.

Secure Laser-Lok abutments onto the Laser-Lok Protective Sleeve to protect Laser-Lok during handling.

8.9 Chairside Ti-Bases

8.9.1 Design:

The scanning, design, and manufacturing processes for a zirconia superstructure are dictated by the Sirona CEREC workflow.

8.9.2 Design Limitations:

Design limitations, including a 20° maximum post angle, are according to the Sirona CAD/CAM System parameters. The critical design parameters of a superstructure intended to mate with the titanium bases originate from Sirona CEREC workflow per K234018 and are as follows:

- Abutment post height ≥4.0 mm (refer to **Figure 1**)
- 20° maximum post angulation

- 0.5mm minimum wall thickness


Figure 1 – Zirconia superstructure design limitations.

8.9.3 Select the corresponding BioHorizons Chairside Ti-Base in the Sirona CEREC Software.

Chairside Ti Base		BioHorizons Ti-Base in the CEREC software
Part No.	Description	
TP3CTB1	3.0mm Chairside Ti-Base, 1mm Emergence	BioHorizons, 3.0
TP3CTB2	3.0mm Chairside Ti-Base, 2mm Emergence	
TP3CTB3	3.0mm Chairside Ti-Base, 3mm Emergence	
PYCTB1	3.5mm Chairside Ti-Base, 1mm Emergence	BioHorizons, 3.5
PYCTB2	3.5mm Chairside Ti-Base, 2mm Emergence	
PYCTB3	3.5mm Chairside Ti-Base, 3mm Emergence	
PGCTB1	4.5mm Chairside Ti-Base, 1mm Emergence	BioHorizons, 4.5
PGCTB2	4.5mm Chairside Ti-Base, 2mm Emergence	
PGCTB3	4.5mm Chairside Ti-Base, 3mm Emergence	
PBCTB1	5.7mm Chairside Ti-Base, 1mm Emergence	BioHorizons, 5.7
PBCTB2	5.7mm Chairside Ti-Base, 2mm Emergence	
PBCTB3	5.7mm Chairside Ti-Base, 3mm Emergence	

8.9.4 BioHorizons Chairside Ti-Base zirconia superstructure:

It is recommended the zirconia superstructure be produced from CEREC Cercon 4D CAD/CAM Blocks (K234018) and sintered according to manufacturer's instruction. The milled superstructure shall be prepared for cementation according to manufacturer's instruction. Cement the milled superstructure to the Chairside Ti-Base titanium component using 3M™ RelyX™ Unicem 2 Automix Self-Adhesive Resin Cement (K022476). The proposed Ti-Base components are assembled through extraoral cement bonding (3M™ RelyX™ Unicem 2 Automix Self-Adhesive Resin Cement (K022476)) with the patient specific CEREC Cercon 4D™ Abutment Block (K234018), to form the complete, two-piece CAD/CAM Titanium Base system abutments.

9.0 WARNINGS AND PRECAUTIONS

- 9.1** Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of BioHorizons prosthetic components.
- 9.2** Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include, but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics.
- 9.3** If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the modifier of the implant/abutment interface a medical device company subject to FDA rules and regulations.
- 9.4** Titanium base abutments - Altering the size or shape of the titanium base abutment may affect proper fit with the restoration. Leaving bonding agent or particulate in the Laser-Lok surface may affect the connective tissue attachment.
- 9.5** Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.
- 9.6** Devices should not be used in patients with a known allergy or sensitivity to the device material(s).
- 9.7** Use of Chairside Ti-Bases with any dental cement, ceramic or other top-half component (and/or hybrid abutment-crown component) materials, scanners, milling units and components, CAD/CAM software, templates, and tools other than those specifically identified in this labeling, may result in incorrect fit and/or damage of the dental restoration.
- 9.8** Small diameter implants are intended for the anterior region of the mouth and are not intended for the posterior region of the mouth due to possible failure of the implant.
- 9.9** Chairside Ti-Bases with 3.4mm implants are indicated for maxillary lateral and mandibular central/lateral incisors only.
- 9.10** Chairside Ti-Base angled abutments are not recommended for the posterior region.
- 9.11** Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at www.biohorizons.com. Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

10.0 COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) abutment screw and/or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

11.0 HANDLING AND STERILIZATION

Always handle the device with powder-free gloves and avoid contact with hard objects that may damage the surface. After abutment preparation and before final placement in the restorative site, CAD/CAM abutments must be cleaned and sterilized. BioHorizons devices have not been validated for automated cleaning. The implant to abutment connection is a critical feature for success of the dental restoration. To avoid potential damage to the abutment connection, a manual cleaning process only was validated which would give greater flexibility of delicate components, and allows for potential inspection after cleaning.

The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a container using a broad-spectrum cleaning agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations. Refer to the legal manufacturer's instructions for use for detergent solution preparation.
- 2) Brush the device to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse device under running utility tap water.
- 4) Place device in the sterile tube filled with the prepared detergent solution and sonicate for two (2) minutes.
- 5) Thoroughly rinse device under running utility tap water.
- 6) Spray device with 70% isopropyl alcohol (IPA).
- 7) Blot device dry with clean lint free cloth.

For sterilization, place final-prepped and cleaned abutment in an FDA cleared sterilization bag or wrap, and run through one of the following qualified sterilization cycles:

Sterilization Cycles				
Reference:	AAMI TIR12	AAMITIR12	AAMI TIR12	UK HTM 01-01 Part C
Type:	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
Exposure Time and Temperature:	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
Minimum Dry Time:	30 minutes	30 minutes	20 minutes	20 minutes

For sterilization of Chairside Ti-Bases bonded with a zirconia superstructure, place product in an FDA cleared sterilization bag or wrap and run through one of the following qualified sterilization cycles:

Sterilization Cycles		
Reference:	AAMI TIR12	UK HTM 01-01 Part C
Type:	Prevacuum Steam	Prevacuum Steam
Exposure Time and Temperature:	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
Minimum Dry Time:	20 minutes	20 minutes

It is recommended to include a 30-minute cool-down period before removing the device from the sterilization bag or wrap.

Devices to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant local regulations.

12.0 MAGNETIC RESONANCE (MR) SAFETY INFORMATION

MR Conditional

BioHorizons implant systems have been demonstrated through non-clinical testing to be magnetic resonance (MR) conditional. A patient with this device can be scanned safely in an MR system under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Static magnetic field of 1.5-Tesla or 3-Tesla, only;
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40-T/m);
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode.

Under the scan conditions defined, the implant system is expected to produce a maximum temperature rise of 3.6°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implant system extends approximately 30-mm from this system when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

Parameter	Condition
Nominal Values of Static Magnetic Field (T)	1.5-Tesla and 3.0-Tesla
Maximum Spatial Field Gradient (T/m and gauss/cm)	40-T/m (4,000-gauss/cm)
Type of RF Excitation	Circularly Polarized (CP) (i.e., Quadrature-Transmission)
Transmit RF Coil Information	There are no transmit RF coil restrictions. Accordingly, the following may be used: body transmit RF coil and all other RF coil combinations (i.e., body RF coil combined with any receive-only RF coil, transmit/receive head RF coil, transmit/receive knee RF coil, etc.)
Operating Mode of MR System	Normal Operating Mode
Maximum Whole Body Averaged SAR	2-W/kg (Normal Operating Mode)
Limits on Scan Duration	Whole body averages SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks).
MR Image Artifact	The presence of this implant system produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant system is located in the area of interest.

To allow medical professionals to identify the specific medical devices a patient has, the MRI safety status of the medical devices, and the conditions for safe use in the MR environment for MR Conditional devices, BioHorizons recommends that clinicians provide the patient with the dental implant(s) and dental abutment(s) / bar(s) device specific peel-offs affixed to the patient card. Patient cards are available free of charge upon request from BioHorizons or for direct printing at <https://ifu.biohorizons.com>

13.0 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

For implantable devices, the summary of safety and clinical performance (SSCP) according to Article 32 of Regulation (EU) 2017/745 can be found in the European database for medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed>.

Implantable Device	Basic UDI-DI Number
BioHorizons Custom Ti abutments	08472360IIBWETABUT007J6
BioHorizons Custom milled abutments	
BioHorizons Titanium base abutments	
BioHorizons Hybrid base abutments	
BioHorizons Chairside Ti-Bases	

14.0 SYMBOLS AND DESCRIPTIONS

The symbols table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic instructions for use
	Manufacturer
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label.

Symbol	Symbol Description
	The four-digit number accompanying the CE mark on applicable devices corresponds to the assigned EU Notified body.
	Reference/ article number
	Lot/ batch number
	Unique Device Identifier
	Do not re-use
	Do not re-sterilize
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
Rx Only	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician
	Authorized representative in the identified country or jurisdiction. The [XX] text of the symbol indicates the two-letter or the three-letter country or jurisdiction code, e.g. EU for European Union.
	Do not use if package is damaged. Discard device and package.
	Medical Device
	Non-sterile
	Single sterile barrier system with protective packaging outside
	Single sterile barrier system
	Home
	Magnetic resonance warning: Device is MR conditional

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc.

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Zimmer®, Dental ScrewVent® og Tapered ScrewVent® er registrerte varemerker tilhørende Zimmer Dental, Inc.