




Tapered Pro implants

Predictable, immediate results

Immediate implant treatment requires predictability. Tapered Pro implants have been developed based on over 10 years of tapered implant success. The unique design elements provide a predictable solution for immediate treatment.



Design features include:

- Tapered body and aggressive threads provide primary stability
- End cutting, self-tapping thread design for controlled implant placement in challenging sites
- Unique Laser-Lok microchannels create connective tissue attachment and retain crestal bone, allowing better control of esthetic outcomes
- Reduced collar diameter preserves vital bone

customer care: **888.246.8338**
shop online at store.biohorizons.com

BIOHORIZONS[®]

Tapered Pro implants & instruments

BioHorizons Tapered Pro implants can be placed with the existing freehand and guided surgery kits. The surgical kits include the instrumentation required for Tapered Pro, Tapered Plus, Tapered Internal and Tapered Tissue Level implants.




TSK4000 – Tapered HD Surgical Kit

CGS4000 – Guided Surgery Kit*

Specifications

Body diameter	3.8mm	4.2mm	4.6mm	5.2mm
Prosthetic connection	3.0mm	3.5mm	3.5mm	4.5mm
Laser-Lok zone	1.8mm	1.8mm	1.8mm	1.8mm
Apical diameter	2.8mm	3.0mm	3.2mm	3.3mm
9.0mm length	BTA3809	BTA4209	BTA4609	BTA5209
10.5mm length	BTA3810	BTA4210	BTA4610	BTA5210
12.0mm length	BTA3812	BTA4212	BTA4612	BTA5212
15.0mm length	BTA3815	BTA4215	BTA4615	BTA5215
18.0mm length	BTA3818	BTA4218	BTA4618	-



*Mount-free Tapered Pro, Tapered PTG, Tapered Internal, Tapered Plus, Tapered 3.0, Tapered Tissue Level implants. 5.2mm and 5.8mm instruments sold separately. BioHorizons® and Laser-Lok® are registered trademarks of BioHorizons. Not all products are available in all countries. ©BioHorizons. All Rights Reserved.

L02056

L02056

REV B APR 2022

REV B APR 2022