









 <p>Tapered Surgical Drill, Grey, 2.5mm TD-25</p>	  <p>Tapered Surgical Drill, Purple, 2.8mm TD-28</p> <p>Tapered Crestal Bone Drill, Purple, 3.3mm CBD-33</p>	  <p>Tapered Surgical Drill, Yellow, 3.2mm TD-32</p> <p>Tapered Crestal Bone Drill, Yellow, 3.8mm CBD-38</p>	  <p>Tapered Surgical Drill, Green, 3.7mm TD-37</p> <p>Tapered Crestal Bone Drill, Green, 4.2mm CBD-42</p>	 <p>Tapered Surgical Drill, Red, 4.1mm TD-41</p>	  <p>Tapered Surgical Drill, Blue, 4.5mm TD-45</p> <p>Tapered Crestal Bone Drill, Blue, 5.0mm CBD-50</p>
--------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



**INTRA-LOCK**  
DENTAL IMPLANTS  
SURGICAL KIT

Ø IMPLANT: 0 IMPLANT, 3.3mm, 3.8mm, 4.2mm, 5.0mm

START: Pilot, 2.5mm, 2.8mm, 3.2mm, 3.7mm, 4.1mm, 4.5mm

Spares: 3.3mm, 3.8mm, 4.2mm, 5.0mm

Ancillary: Straight Pin, Drill Extender, Ratchet Extender

Hex Driver, Hand Wrench, Handpiece, Ratchet

Implant Staging: 0, 7.5, 9, 10.5, 12, 15

Ratchet R-4MM

Intra-Lock Surgical Kit IL-SK

Spare slots can be used to hold ancillary instrumentation like pilot drills and bone profilers. The implant staging area is designed to temporarily secure implant vials during surgery only. Do not attempt to sterilize implant vials in the surgical tray. Images shown not to scale.



## FUSION SURGICAL KIT REFERENCE CHART



L02079



REV A NOV 2021