

Remote Anchorage Surgical Kit

Drill Usage Chart

Part description	Part number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
MARKING DRILL Ø2.0MM, ZYGOMA	DM20Z																				
MARKING DRILL Ø2.0MM, PTERYGOID	DM20P																				
ROUND BUR Ø3.8MM	DBB 38-BH																				
CYLINDRICAL BUR Ø3.8MM	DCB 38-BH																				
CYLINDRICAL BUR Ø4.2MM	DCB 42-BH																				
ZYGOMA PILOT DRILL Ø1.9MM SHORT	DZP19S																				
ZYGOMA DRILL Ø3.8MM SOFT BONE SHORT	DZS38S																				
ZYGOMA DRILL Ø3.8MM DENSE BONE SHORT	DZD38S																				
ZYGOMA DRILL Ø4.2MM SOFT BONE SHORT	DZS42S																				
ZYGOMA DRILL Ø4.2MM DENSE BONE SHORT	DZD42S																				
ZYGOMA PILOT DRILL Ø1.9MM MEDIUM	DZP19M																				
ZYGOMA DRILL Ø3.8MM SOFT BONE MEDIUM	DZS38M																				
ZYGOMA DRILL Ø3.8MM DENSE BONE MEDIUM	DZD38M																				
ZYGOMA DRILL Ø4.2MM SOFT BONE MEDIUM	DZS42M																				
ZYGOMA DRILL Ø4.2MM DENSE BONE MEDIUM	DZD42M																				
ZYGOMA PILOT DRILL Ø1.9MM LONG	DZP19L																				
ZYGOMA DRILL Ø3.8MM SOFT BONE LONG	DZS38L																				
ZYGOMA DRILL Ø3.8MM DENSE BONE LONG	DZD38L																				
ZYGOMA DRILL Ø4.2MM SOFT BONE LONG	DZS42L																				
ZYGOMA DRILL Ø4.2MM DENSE BONE LONG	DZD42L																				

Important note:

All drills should be replaced after they have been used 12 to 20 times, depending on bone density.

customer care: 888.246.8338 • shop online: [biohorizons.com](https://www.biohorizons.com)