

Laser-Lok 3.0 Drill Usage Chart

Description	Part Number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	
2.0mm Starter Drill - Extended Shank	122-403																					
2.5mm Depth Drill - Extended Shank	122-425																					
2.5x10.5mm Depth Drill - Extended Shank	122-425105																					
2.5x12mm Depth Drill - Extended Shank	122-42512																					
2.5x15mm Depth Drill - Extended Shank	122-42515																					
Bone Tap	122-900																					
Crestal Bone Drill	TP3CBD																					

Important Note: All drills should be replaced after they have been used 12 to 20 times, depending upon bone density.

Description	Part Number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	
2.0mm Starter Drill - Extended Shank	122-403																					
2.5mm Depth Drill - Extended Shank	122-425																					
2.5x10.5mm Depth Drill - Extended Shank	122-425105																					
2.5x12mm Depth Drill - Extended Shank	122-42512																					
2.5x15mm Depth Drill - Extended Shank	122-42515																					
Bone Tap	122-900																					
Crestal Bone Drill	TP3CBD																					

Important Note: All drills should be replaced after they have been used 12 to 20 times, depending upon bone density.

Shop on-line at www.biohorizons.com