

make the switch


The Tapered Plus implant system offers all the great benefits of BioHorizons highly successful Tapered Internal system PLUS it features a Laser-Lok treated beveled-collar for bone and soft tissue attachment and platform switching designed for increased soft tissue volume.

platform switching

designed to increase soft tissue volume around the implant connection

Laser-Lok® zone

creates a connective tissue seal and maintains crestal bone



optimized threadform

buttress thread engineered for superior stability over microthreaded implants



restorative ease

conical internal hex connection is color-coded for quick identification and component matching



tapered internal plus



3.0mm body / 3.0 platform

TLX3010 10.5mm length
TLX3012 12mm length
TLX3015 15mm length



3.8mm body / 3.0 platform

TLXP3809 9mm length
TLXP3810 10.5mm length
TLXP3812 12mm length
TLXP3815 15mm length



4.6mm body / 3.5 platform

TLXP4607 7.5mm length
TLXP4609 9mm length
TLXP4610 10.5mm length
TLXP4612 12mm length
TLXP4615 15mm length



5.8mm body / 4.5 platform

TLXP5807 7.5mm length
TLXP5809 9mm length
TLXP5810 10.5mm length
TLXP5812 12mm length
TLXP5815 15mm length

The new BioHorizons Surgical Kit - now available

- use to place Tapered Plus, Tapered Internal, Tapered 3.0 and Laser-Lok 3.0
- 40% smaller and 40% lighter than other kits
- versatile, removable, hinged lid
- easy to disassemble and assemble during cleaning
- implant staging area for implant vials during surgery



Tapered Internal Plus **now available**,
place your order today.

Customer Care: 1.888.246.8338
visit us online at www.biohorizons.com

BIOHORIZONS[®]
SCIENCE • INNOVATION • SERVICE

Tapered Internal 3.0 is not platform switched due to its small diameter.

Tapered Plus is packaged mount-free for quick placement and maximum site visibility. Comes with a cover cap. Titanium Alloy (Ti-6AL-4V) for maximum strength.

Product descriptions have been modified for presentation purposes. Not all products shown are available in all countries. © 2012 BioHorizons. All Rights Reserved.

ML0107 Rev A OCT 2012