

External Drill Usage Chart

Description	Part Number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
2.0mm Starter Drill	122-103																				
2.5mm Depth Drill	122-225																				
2.5mm x 9.0mm Depth Drill	122-12509																				
2.5mm x 10.5mm Depth Drill	122-125105																				
2.5mm x 12mm Depth Drill	122-12512																				
2.5mm x 15mm Depth Drill	122-12515																				
2.0mm Starter Drill, Extended Shank	122-403																				
2.5mm x 7.0mm Depth Drill, Extended Shank	122-42507																				
2.5mm x 9.0mm Depth Drill, Extended Shank	122-42509																				
2.5mm x 10.5mm Depth Drill, Extended Shank	122-425105																				
2.5mm x 12mm Depth Drill, Extended Shank	122-42512																				
2.5mm x 15mm Depth Drill, Extended Shank	122-42515																				
2.5mm Depth Drill, Extended Shank	122-425																				
3.0mm Width Increasing Drill	122-230																				
3.0mm Width Increasing Drill, Extended Shank	122-430																				
3.5mm Crestal Bone Drill	123-000																				
3.5mm D3 Bone Tap	123-300																				
3.4mm Width Increasing Drill	122-232																				
3.4mm Width Increasing Drill, Extended Shank	122-432																				
4.0mm Crestal Bone Drill	124-000																				
4.0mm D3 Bone Tap	124-300																				
3.9mm Width Increasing Drill	122-237																				
3.9mm Width Increasing Drill, Extended Shank	122-437																				
4.4mm Width Increasing Drill	122-242																				
4.4mm Width Increasing Drill, Extended Shank	122-442																				
5.0mm Crestal Bone Drill	125-000																				
5.0mm D3 Bone Tap	125-300																				
4.9mm Width Increasing Drill	122-247																				
4.9mm Width Increasing Drill, Extended Shank	122-447																				
5.4mm Width Increasing Drill	122-252																				
5.4mm Width Increasing Drill, Extended Shank	122-452																				
6.0mm Crestal Bone Drill	126-000																				
6.0mm D3 Bone Tap	126-300																				

Important Note:
All drills should be replaced after they have been used 12 to 20 times, depending upon bone density.

Shop on-line at www.biohorizons.com