

W&H surgical motor kits



Base Kit \$5,277
 WH-90000216
with ISQ Module* \$8,643
 WH-90000217

- Implantmed SI-1015 console
- Handheld motor EM-19
- Foot pedal
- (3) Disposable irrigation tubes
- WI-75 surgical 20:1 contra-angle
- Bur testing gauge
- Service oil and spray cap

*W&H Osstell ISQ Module



WiFi Foot Pedal Upgrade Kit \$5,767
 WH-90000218
with ISQ Module* \$8,933
 WH-90000219

- Implantmed SI-1015 console
- Handheld motor EM-19
- Wireless foot pedal
- (3) Disposable irrigation tubes
- WI-75 surgical 20:1 contra-angle
- Bur testing gauge
- Service oil and spray cap

*W&H Osstell ISQ Module



LED+ Motor Upgrade Kit \$6,142
 WH-90000220
with ISQ Module* \$9,308
 WH-90000221

- Implantmed SI-1015 Console
- LED motor EM-19 LC
- Foot pedal
- (3) Disposable irrigation tubes
- WS-75 L surgical 20:1 contra-angle
- Bur testing gauge
- Service oil and spray cap

*W&H Osstell ISQ Module



Fully Loaded Kit \$6,332
 WH-90000222
with ISQ Module* \$9,368
 WH-90000223

- Implantmed SI-1015 console
- LED motor EM-19 LC
- Wireless foot pedal
- (3) Disposable irrigation tubes
- WS-75 L surgical 20:1 contra-angle
- Bur testing gauge
- Service oil and spray cap

*W&H Osstell ISQ Module



new straight & contra-angle handpieces with LED+ & W&H connection

- Optimal illumination at all stages of the work
- Work with daylight quality
- Natural reproduction of red tones
- Can be dismantled for thorough cleaning
- Work with W&H motor EM-19 LC



WS-75 L \$1,602
 WH-30033000
 • speed reducing 20:1
 • press-button chuck system, for surgical burs and cutters with contra-angled shank Ø 2.35 mm E/KM internal cooling system



WS-56 L \$1,469
 WH-30040000
 • 1:1 ratio
 • press-button chuck system, for surgical burs and cutters with contra-angled shank Ø 2.35 mm



WS-91 L \$1,469
 WH-30049000
 • speed increasing 1:2.7
 • 45° head, triple spray, press-button chuck system, for surgical burs and cutters with FG shank Ø 1.6 mm



WS-92 L \$1,469
 WH-30043000
 • speed increasing 1:2.7
 • triple spray, press-button chuck system, for surgical burs and cutters with FG shank Ø 1.6 mm



S-11 L \$1,581
 WH-30058000
 • 1:1 ratio
 • straight, lever chuck system, for surgical burs and cutters Ø 2.35 mm



	WI-75 E/KM (WS-75*)	WS-75 L G	WS-56	WS-56 L G	WS-91	WS-91 L G	WS-92	WS-92 L G	S-9 S-11	S-9 L G S-11 L G	S-10	S-12	S-15
Transmission ratio	20:1	20:1	1:1	1:1	1:2,7				1:1	1:1	1:1	1:2	1:1
Motor speed max (RPM)	50,000	40,000	50,000	40,000	50,000	40,000	50,000	40,000	50,000	40,000	50,000	40,000	30,000
LED light (independent generator)	N/A	Mini LED+	N/A	Mini LED+	N/A	Mini LED+	N/A	Mini LED+	N/A	Mini LED+	N/A	N/A	N/A
Light quality	N/A	Daylight	N/A	Daylight	N/A	Daylight	N/A	Daylight	N/A	Daylight	N/A	N/A	N/A
Coupling system	ISO 3964 (DIN 13.940)				ISO 3964 (DIN 13.940)				ISO 3964 (DIN 13.940)				
External spray	Single				Triple				External				
Internal cooling system	Kirschner/Meyer		N/A		N/A				N/A				
Chucking system	Push-button chuck				Push-button chuck				Lever chuck				
Rotary instruments	Surgical burs, drills, cutters with contra-angled shank Ø2.35mm				Surgical burs, drills, cutters with FG shank Ø1.6mm				Surgical burs, drills, cutters, Ø2.35mm				
Maximum length									45mm	45mm	70mm	70mm	45mm
Dismantable	No *WS-75 - yes	Yes	Yes	Yes	Yes				Yes				
Indication, e.g.	Implantation, segmental osteotomy		Maxillary and mandibular osteotomy, gementectomy, sequestrectomy		Hemisection, wisdom teeth extraction				Apical resection, osteotomy, apical trepanning, bone modelling		Osteotomy, fenestration, osteosynthesis in distal areas		
BioHorizons product code	WH-10207550 (WH-30065000)	WH-30032000	WH-30038000	WH-30039000	WH-30047000	WH-30048000	WH-30041000	WH-30042000	WH-30053000 WH-30056000	WH-30054000 WH-30057000	WH-30059000	WH-30061000	WH-30060000
Price	\$1,298 (\$1,612)	\$1,912	\$1,374	\$1,762	\$1,464	\$1,762	\$1,464	\$1,762	\$1,584 \$1,226	\$2,018 \$1,882	\$1,703	\$1,703	\$1,435

For more information, contact your local Territory Manager or call Customer Care: **1.888.246.8338**

