



BIOHORIZONS®

Mem-Lok® Comparison

Mem-Lok® RCM

Mem-Lok® Pliable

Mem-Lok® Pericardium

Composition/Source	Type I collagen derived from achilles tendon	Purified porcine peritoneum tissue	Human derived tissue
Handling	Excellent memory, ideal for space maintenance	Strong, dense single collagen layer, conformable	Conformable, adapts easily to graft site contours
Primary Applications	<ul style="list-style-type: none"> - Socket grafts - Sinus grafts - Ridge augmentation 	<ul style="list-style-type: none"> - Socket grafts - Ridge augmentation - Block grafting 	<ul style="list-style-type: none"> - Socket grafts - Ridge augmentation - Block grafting
Storage	15-30°C	15-30°C	15-30°C
Shelf Life	3 years	2 years	5 years
Resorption Time	26-38 weeks	12-16 weeks	16-24 weeks
Can be sutured/tacked?	Yes	Yes	Yes

ordering information

RCM-ML1520	15mm x 20mm	\$119	PBLE-ML1520	15mm x 20mm	\$119	PERI-ML1520INT	15mm x 20mm	\$145
RCM-ML2030	20mm x 30mm	\$139	PBLE-ML2030	20mm x 30mm	\$139	PERI-ML2030INT	20mm x 30mm	\$215
RCM-ML3040	30mm x 40mm	\$194	PBLE-ML3040	30mm x 40mm	\$194	PERI-ML3040INT	30mm x 40mm	\$265

BioHorizons® and Mem-Lok® are registered trademarks of BioHorizons. Not all products shown or described are available in all countries. © BioHorizons. All Rights Reserved.

For more information, contact your local Territory Manager or call Customer Care: **1.888.246.8338**
Order 24/7 at store.biohorizons.com

