



BIOHORIZONS®

Mem-Lok® Comparison

	Mem-Lok® RCM	Mem-Lok® Pliable	Mem-Lok® Pericardium
Processor	N/A	N/A	DCI Donor Services Tissue Bank
Manufacturer	Collagen Matrix, Inc.	Collagen Matrix, Inc.	N/A
Handling	Excellent memory, ideal for space maintenance	Strong, dense single collagen layer, conformable	Conformable, adapts easily to graft site contours
Primary Applications	<ul style="list-style-type: none"> - Socket grafts - Sinus grafts - Ridge augmentation 	<ul style="list-style-type: none"> - Socket grafts - Ridge augmentation - Block grafting 	<ul style="list-style-type: none"> - Socket grafts - Ridge augmentation - Block grafting
Storage	15-30°C	15-30°C	15-30°C
Shelf Life	3 years	18 months	5 years
Resorption Time	26-38 weeks	12-16 weeks	16-24 weeks
Can be sutured/tacked?	Yes	Yes	Yes
Distribution	Global	Global*	US only

ordering information

RCM-ML1520	15mm x 20mm	\$ 119
RCM-ML2030	20mm x 30mm	\$ 139
RCM-ML3040	30mm x 40mm	\$ 194

PBLE-ML1520	15mm x 20mm	\$ 119
PBLE-ML2030	20mm x 30mm	\$ 139
PBLE-ML3040	30mm x 40mm	\$ 194

PERI-ML1520	15mm x 20mm	\$ 145
PERI-ML2030	20mm x 30mm	\$ 215
PERI-ML3040	30mm x 40mm	\$ 265

BioHorizons® and Mem-Lok® are registered trademarks of BioHorizons. Mem-Lok® RCM and Mem-Lok® Pliable are manufactured by Collagen Matrix, Inc. Mem-Lok® Pericardium is processed by DCI Donor Services Tissue Bank. Not all products shown or described are available in all countries. © BioHorizons. All Rights Reserved.

*Call to check for availability

For more information, contact your local Territory Manager or call Customer Care:

888.246.8338