

Mem-Lok® comparison



	Mem-Lok® Amnio	Mem-Lok® Pliable	Mem-Lok® Pericardium	Mem-Lok® RCM
Composition/Source	Human derived, amnion-chorion placenta tissue	Purified porcine peritoneum tissue	Human derived, heart-lining tissue	Type I collagen derived from achilles tendon
Handling	Superior flexibility, lightweight, and adaptable handling experience	Drapes easily, dense single collagen layer,	Strong and tear resistant, adapts easily to graft site contours	Excellent memory and longevity, ideal for space maintenance
Primary Applications	<ul style="list-style-type: none"> • Socket preservation • Sinus perforations • Ridge augmentation • Periodontal regeneration • Gingival recession • Endodontic applications 	<ul style="list-style-type: none"> • Extraction sites • Periodontal defects • Graft containment around immediate implants • Ridge augmentation 	<ul style="list-style-type: none"> • Ridge augmentation • Extraction sites 	<ul style="list-style-type: none"> • Extraction sockets • Sinus augmentation • Ridge preservation • Bone augmentation around implants • Bony defects • Peri-implant bone defect around implants
Storage	15-30°C	15-30°C	15-30°C	15-30°C
Shelf Life	5 years	3 years	5 years	3 years
Resorption Time	8-12 weeks	12-16 weeks	16-24 weeks	26-38 weeks
Can be sutured/tacked?	Yes	Yes	Yes	Yes

Ordering information

Mem-Lok® Amnio

Item Number	Size	Price
AMN-ML1212	12mm x 12mm	\$99
AMN-ML1520	15mm x 20mm	\$159
AMN-ML2030	20mm x 30mm	\$289

Mem-Lok® Pliable

Item Number	Size	Price
PBLE-ML1520	15mm x 20mm	\$129
PBLE-ML2030	20mm x 30mm	\$154
PBLE-ML3040	30mm x 40mm	\$219

Mem-Lok® Pericardium

Item Number	Size	Price
PERI-ML1520INT	15mm x 20mm	\$145
PERI-ML2030INT	20mm x 30mm	\$215
PERI-ML3040INT	30mm x 40mm	\$265

Mem-Lok® Resorbable Collagen Membrane

Item Number	Size	Price
RCM-ML1520	15mm x 20mm	\$129
RCM-ML2030	20mm x 30mm	\$154
RCM-ML3040	30mm x 40mm	\$219

BioHorizons® and Mem-Lok® are registered trademarks of BioHorizons. Not all products shown or described are available in all countries. © BioHorizons. All Rights Reserved. SPMP15053 REV F AUG 2022

For more information, contact your local Territory Manager or call Customer Care: **1.888.246.8338**
Order 24/7 at store.biohorizons.com

