



# BIOHORIZONS®

## Mem-Lok® Comparison

### Mem-Lok® RCM

### Mem-Lok® Pliable

### Mem-Lok® Pericardium

<b>Composition/Source</b>	Type I collagen derived from achilles tendon	Purified porcine peritoneum tissue	Human derived tissue
<b>Handling</b>	Excellent memory, ideal for space maintenance	Strong, dense single collagen layer, conformable	Conformable, adapts easily to graft site contours
<b>Primary Applications</b>	<ul style="list-style-type: none"> <li>- Socket grafts</li> <li>- Sinus grafts</li> <li>- Ridge augmentation</li> </ul>	<ul style="list-style-type: none"> <li>- Socket grafts</li> <li>- Ridge augmentation</li> <li>- Block grafting</li> </ul>	<ul style="list-style-type: none"> <li>- Socket grafts</li> <li>- Ridge augmentation</li> <li>- Block grafting</li> </ul>
<b>Storage</b>	15-30°C	15-30°C	15-30°C
<b>Shelf Life</b>	3 years	3 years	5 years
<b>Resorption Time</b>	26-38 weeks	12-16 weeks	16-24 weeks
<b>Can be sutured/tacked?</b>	Yes	Yes	Yes

### ordering information

<b>RCM-ML1520</b>	15mm x 20mm	\$119	<b>PBLE-ML1520</b>	15mm x 20mm	\$119	<b>PERI-ML1520INT</b>	15mm x 20mm	\$145
<b>RCM-ML2030</b>	20mm x 30mm	\$139	<b>PBLE-ML2030</b>	20mm x 30mm	\$139	<b>PERI-ML2030INT</b>	20mm x 30mm	\$215
<b>RCM-ML3040</b>	30mm x 40mm	\$194	<b>PBLE-ML3040</b>	30mm x 40mm	\$194	<b>PERI-ML3040INT</b>	30mm x 40mm	\$265

BioHorizons® and Mem-Lok® are registered trademarks of BioHorizons. Not all products shown or described are available in all countries. © BioHorizons. All Rights Reserved.

For more information, contact your local Territory Manager or call Customer Care: **1.888.246.8338**  
Order 24/7 at [store.biohorizons.com](http://store.biohorizons.com)

