

Elos Accurate® Scan Body for Multi-unit Abutments

The new Elos Accurate Scan Body for the BioHorizons Multi-unit system now provides an improved digital workflow for both dentists and dental technicians. The single Scan Body design can be used on the BioHorizons Multi-unit analog for scanning and creation of a digital model. This allows for incorporation into CAD/CAM workflows for the final restoration.

The Scan Body is manufactured from biocompatible titanium alloy (Ti-6Al-4V) and polyetheretherketone (PEEK).




Ordering Information



IO 11B-A

Elos Accurate Scan Body for Multi-unit Abutments



C13485

Elos Accurate Scan Body Driver

Library Information

The Elos Accurate Scan Body is included in the digital library provided by Vulcan Custom Dental. The digital files can be used to plan and design restorations in 3Shape and Exocad.

<http://www.vulcandental.com/Home/Downloads>

BioHorizons® is a registered trademark of BioHorizons. Elos Accurate® is a registered trademark of Elos MedTech. Not all products are available in all countries. ©BioHorizons. All Rights Reserved.

ELOS Accurate® Scan Body IFU: <https://elosmedtech.com/IFU/>

customer care: **888.246.8338**
shop online at store.biohorizons.com

L02064 REV C AUG 2022

