




International Biologics Applications Matrix

bone grafting options




	 MinerOss® Cortical & Cancellous Chips	 Grafton® DBM Matrix Plug	 Grafton® DBM Matrix Putty	 MinerOss® XP Cancellous	 MinerOss® X Cancellous	 MinerOss® X Cortical	 MinerOss® X Collagen
periodontal defects	✓	✓	✓	✓	✓		✓
4-wall extraction sockets	✓	✓	✓	✓	✓	✓	✓
3-wall extraction sockets	✓	✓	✓	✓	✓	✓	✓
sinus augmentation	✓	✓	✓	✓	✓	✓	✓
ridge augmentation	✓			✓	✓	✓	
grafting for implant placement	✓			✓	✓	✓	
block grafting				✓	✓		

For more information, contact BioHorizons Customer Care





International Biologics Applications Matrix

soft tissue augmentation

	 NovoMatrix™ Reconstructive Tissue Matrix	 AlloDerm® Regenerative Tissue Matrix (RTM)	 AlloDerm GBR® Regenerative Tissue Matrix (RTM)
root coverage	✓	✓	✓
gingival augmentation	✓	✓	
soft tissue augmentation around implants	✓	✓	
soft tissue ridge augmentation	✓	✓	
graft protection & containment			✓
flap extender to achieve primary closure			✓

dental membranes

	 Mem-Lok® Resorbable Collagen Membrane (RCM)	 Mem-Lok® Pliable	 CYTOPLAST™ Translucent Resorbable Membrane	 BioPlug	 BioStrip
graft protection & containment	✓	✓	✓	✓	
sinus lift	✓				
ridge preservation	✓	✓	✓		✓
bony defects	✓	✓	✓		

MinerOss® and Mem-Lok® are registered trademarks of BioHorizons. AlloDerm® and AlloDerm GBR® are registered trademarks of LifeCell Corporation. Grafton® DBM is a registered trademark of Medtronic, Inc. Cytoplast™ is a registered trademark of Osteogenics Biomedical, Inc. Not all products shown or described are available in all countries.

For more information, contact BioHorizons Customer Care

