


Tapered Ridge Expanders

The BioHorizons Tapered Ridge Expanders provide osseodensification of the site by compacting bone while expanding the ridge. Tapered Ridge Expanders feature a 4mm square and can be used with a ratchet or hand wrench.

- TRE30** 3.0mm Tapered Ridge Expander
- TRE34** 3.4mm Tapered Ridge Expander
- TRE38** 3.8mm Tapered Ridge Expander
- TRE42** 4.2mm Tapered Ridge Expander
- TRE46** 4.6mm Tapered Ridge Expander



surgical protocol for a 3.8mm implant



BioHorizons® is a registered trademark of BioHorizons. Not all products are available in all countries. ©BioHorizons. All Rights Reserved.

For more information, contact your local Territory Manager
or visit store.biohorizons.com

BIOHORIZONS®