



Pro Freehand Surgical Kit (FSK5000)

Includes the instrumentation required to place: Tapered Pro Conical, Tapered Short Conical, Tapered Pro, Tapered Plus, Tapered Internal, Tapered 3.0, Tapered PTG, Tapered Short & Tapered Tissue Level.











Ratchet
130-000



Straight Parallel Pins (2 per kit) 144-100



Drill Extender
122-100



Hand Wrench
300-400



.050" (1.25mm) Hex Driver
135-351



4mm Square Drive Extender
300-205

Spare slots can be used to hold preferred surgeon instrumentation such as depth indicators and extended shank drills.
The implant staging area is designed to temporarily secure implant vials during surgery only. Do not attempt to sterilize implant vials in the surgical tray.

©BioHorizons Implant Systems Inc.
All Rights Reserved. Images shown not to scale.

customer care: 888.246.8338 • shop online at store.biohorizons.com

L02097 REV B JUL 2026