# SmartShape® Healers

Efficiency meets immediacy –
fewer appointments, faster results



biohorizons camlog group

BIOHORIZONS®

BIOHORIZONS®


# SmartShape® Healers

Traditional implant workflows often require multiple visits, repeated component changes, and significant chairtime for both clinicians and patients. Each time a healing abutment is removed for scanning or temporization, peri-implant tissue is disrupted — increasing the risk of bone remodeling and compromising esthetic outcomes.

This repeated "disconnect–reconnect" cycle can extend treatment to five or more visits and introduce restorative inaccuracies. A more efficient approach is needed — one that preserves the biology from day one while reducing procedural complexity. SmartShape is designed to deliver exactly that.





## Traditional workflow • 6 - 8 months




Extract & graft | Place implant | Fabricate custom healer | Scan | Final

**5 visits. Multiple components. Complex workflow.**



## Benefits for surgical specialists

Designed for efficiency with fewer appointments by integrating a healing abutment, scan body, and anatomical profile into a single component.

Simplifies the restorative process for your referring doctors.

# Faster treatment, fewer visits - up to 60% chairtime



## SmartShape® Healers

### The All-In-One Solution

**Integrated Multi-Functional Design:**
Combines the healing abutment, scan body, and anatomical profile into a single component, eliminating the need for multiple restorative procedures.


**Biological Stability & Tissue Preservation:**
Clinical studies suggest that minimizing the removal and reconnection of implant components is associated with improved crestal bone stability and more predictable restorative outcomes.[2,3]

**Simplified Clinical Workflow:**
Fewer components streamline the procedure, reducing chairtime and the overall number of clinical steps.

**Superior Esthetic & Patient Outcomes:**
Maintains a consistent, natural emergence profile while improving the overall patient experience through fewer adjustments and appointments.

## SmartShape workflow • 2 months

**Visit 1**



Extract | Place implant | Place SmartShape | Scan

**Visit 2**



Final

## 2 visits. 0 removals. 60% faster.[1]




## Benefits for restorative dentists


All-in-one solution reduces the number of appointments, streamlining predictable esthetic outcomes.

Consistent emergence profiles from healing to the final restoration reduce seating issues and manual crown adjustments.


# reduction for specialists, referrals and patients[1]




**Occlusal markings**
for identification and orientation




**Matte surface**
for easy scanning




**Anatomic preservation**
maintain soft tissue profile from healer to final







**Occlusal-sealing screw design**
no block-out required



**All-metal design**
for seating verification



**Seamless restorations**
profile-matched custom final supported by Vulcan Custom Dental




## Benefits for labs

Standardized emergence profiles save design time.


Open libraries allow anatomical profile matching to reduce chairside adjustments.

# SmartShape® part numbering & identification

- Anterior 1 (A1)
- Anterior 2 (A2)
- Mid 1 (M1)
- Mid 2 (M2)
- Posterior 1 (P1)
- Posterior 2 (P2)





## Regular height



| | | A1 | A2 | M1 | M2 | P1 | P2 |
|---|---|---|---|---|---|---|---|
| Hexed | 3.0mm | TP3SSHA1 | TP3SSHA2 | TP3SSHM1 | TP3SSHM2 | - | - |
| | 3.5mm | PYSSHA1 | PYSSHA2 | PYSSHM1 | PYSSHM2 | PYSSHP1 | PYSSHP2 |
| | 4.5mm | - | - | - | PGSSHM2 | PGSSHP1 | PGSSHP2 |
| | 5.7mm | - | - | - | - | PBSSHP1 | PBSSHP2 |
| Conical | Narrow | CNSSHA1 | CNSSHA2 | CNSSHM1 | CNSSHM2 | - | - |
| | Regular | CRSSHA1 | CRSSHA2 | CRSSHM1 | CRSSHM2 | CRSSHP1 | CRSSHP2 |

## Tall height




| | | A1 | A2 | M1 | M2 | P1 | P2 |
|---|---|---|---|---|---|---|---|
| Hexed | 3.0mm | TP3SSHA1T | TP3SSHA2T | TP3SSHM1T | TP3SSHM2T | - | - |
| | 3.5mm | PYSSHA1T | PYSSHA2T | PYSSHM1T | PYSSHM2T | PYSSHP1T | PYSSHP2T |
| | 4.5mm | - | - | - | PGSSHM2T | PGSSHP1T | PGSSHP2T |
| | 5.7mm | - | - | - | - | PBSSHP1T | PBSSHP2T |
| Conical | Narrow | CNSSHA1T | CNSSHA2T | CNSSHM1T | CNSSHM2T | - | - |
| | Regular | CRSSHA1T | CRSSHA2T | CRSSHM1T | CRSSHM2T | CRSSHP1T | CRSSHP2T |


***References:*** *1. Designed to reduce total office visits by up to 60%*
*2. Ayyildiz BG, et al, Int J Oral Maxillofac Implants. Published online September 4, 2025.*
*3. Calatrava, J. et al. Clin Implant Dent Relat Res. 2024;26(5):998-1011.*



Not all products are available in all countries.

# SmartShape® abutment heights


The SmartShape Healers are available in two gingival heights, Regular and Tall. Detailed heights of the abutments are listed in the table below.




Anterior




Mid




Posterior

## Conical gingival height





| | | Emergence Height | Anatomical Profile Height | | Max Tissue Height | |
|---|---|---|---|---|---|---|
| | | | Buccal-Lingual | Mesial-Distal | Buccal-Lingual | Mesial-Distal |
| Regular | Anterior | 2mm | 1.25mm | 2.75mm | 3.25mm | 4.75mm |
| | Mid | 2mm | 1.50mm | 2.25mm | 3.50mm | 4.25mm |
| | Posterior | 2.25mm | 0.75mm | | 3mm | |
| Tall | Anterior | 4mm | 1.25mm | 2.75mm | 5.25mm | 6.75mm |
| | Mid | 4mm | 1.50mm | 2.25mm | 5.50mm | 6.25mm |
| | Posterior | 3.75mm | 0.75mm | | 4.5mm | |

## Hex gingival height





| | | Emergence Height | Anatomical Profile Height | | Max Tissue Height | |
|---|---|---|---|---|---|---|
| | | | Buccal-Lingual | Mesial-Distal | Buccal-Lingual | Mesial-Distal |
| Regular | Anterior | 2mm | 1.25mm | 2.75mm | 3.25mm | 4.75mm |
| | Mid | 2mm | 1.45mm | 2.20mm | 3.45mm | 4.20mm |
| | Posterior | 2mm | 1mm | | 3mm | |
| Tall | Anterior | 4mm | 1.25mm | 2.75mm | 5.25mm | 6.75mm |
| | Mid | 4mm | 1.45mm | 2.20mm | 5.45mm | 6.20mm |
| | Posterior | 3.5mm | 1mm | | 4.5mm | |